
Elan^{IC}

Clinician's Manual

ELANIC22L1S ELANIC30R8S
ELANIC22L1SD ELANIC30R8SD

EN	Clinician's Manual	2
FR	Manuel de l'orthoprothésiste	39
DE	Gebrauchsanweisung für Techniker	76
IT	Manuale per il tecnico ortopedico	113
ES	Manual del técnico	150
NL	Handleiding voor artsen	187
PL	Podręcznik użytkownika dla lekarza	224
PT	Manual do profissional de saúde	261

Contents	2
1 Description and Intended Purpose	3
2 Safety Information	6
3 Construction	7
4 Function.....	8
5 Maintenance	9
6 Limitations on Use	9
7 Battery Charging	10
7.1 Setting Up	10
7.2 Charging.....	10
7.3 Warnings	11
7.4 Battery Level Indication	12
7.5 Device Status Indication.....	12
8 Bench Alignment	13
9 Install the Blatchford PC Software Interface.....	14
10 Using the PC Software Interface to set up the Device	15
10.1 Set-up Procedure	15
10.2 Launch the PC Interface Program	17
10.3 <i>Bluetooth</i> ® Connection [Start].....	18
10.4 Status Bar	19
11 Static Alignment	20
11.1 Static Alignment Procedure	20
12 Biomimetic Alignment	21
12.1 Biomimetic Adjustment.....	22
13 Dynamic Alignment	23
13.1 Dynamic Alignment Procedure.....	23
13.2 Resistance Setting Troubleshooting	24
14 Walking Calibration	25
14.1 Walking Calibration Procedure	25
15 Options	26
16 Fitting Advice.....	27
17 Maintenance Instructions	28
17.1 Footshell Removal	28
17.2 Spring Replacement.....	29
18 Technical Data	31
19 Ordering Information	32
19.1 Replacement Parts.....	32
Appendix 1 Quick Guide to LED Indications	33
Manufacturer's Declarations and Legal Information	35

1 Description and Intended Purpose

These instructions are for use by the practitioner.

The term *device* is used throughout these instructions for use to refer to ElanIC.

Ensure that the user has understood all instructions for use, drawing particular attention to the sections regarding maintenance and use of the inductive battery charger connection.

Application

The device is to be used exclusively as part of a lower limb prosthesis.

Intended for single user.

The device is a biomimetic foot and ankle system with microprocessor controlled dorsiflexion and plantarflexion. It is designed to adapt dynamically (via motorized valves) to provide assistance when walking on a variety of surfaces, when changing gait speeds and increased levels of safety and security when walking and standing.

Features

(See also Section 4 *Function*)

- Dynamic damping resistance aims to improve gait characteristics enabling safer walking
- Microprocessor control of resistance in real time
- The device Software adapts to individual user gait characteristics
- Standing Support mode
- Aims to reduce stress at the socket interface
- PC interface setup with **Bluetooth®** link and direct communications
- Power is supplied via an integral rechargeable battery pack
- Recharging through an easy-to-use magnetically self-aligning inductive charger
- Waterproof according to the parameters described in this manual

This device provides limited self-alignment of the prosthesis on varied terrain and following footwear changes. It is intended to improve postural sway, balance stability and inter-limb loading symmetry while easing high pressures at the socket interface.

Activity Level

This device is recommended for users that have the potential to achieve Activity Level 3 who may benefit from enhanced stability and an increase in confidence on uneven surfaces.

There are exceptions and in our recommendation we want to allow for unique, individual circumstances. There may also be a number of users of Activity Levels 2 and 4* who would benefit from the enhanced stability offered by the device.

Users at Activity Levels 2 and 4* will require softer or stiffer springs as appropriate for the individual user rather than as shown in the spring selection guide.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

*for maximum user weight - 100 kg (220 lbs) - always use one higher spring rate category than shown in the Spring Set Selection table.

Contraindications

This device may not be suitable for Activity Level 1 individuals or for competitive sports events, as these types of users will be better served by a specially designed prosthesis optimized for their needs.

The device and its charger must not be used by users of pacemakers or implantable cardioverter-defibrillators (ICDs).

Clinical Benefits

- Increased ground clearance reduces risk of trips and falls
- Improved balance through self-alignment and standing support
- Improved control and safety on slope negotiation
- Improved kinetic gait symmetry
- Reduced loading on the residual limb
- Increased walking speed

Spring Set Selection

Activity Level 3

44-52	53-59	60-68	69-77	78-88	89-100	101-116	117-125	kg	User Weight
(100-115)	(116-130)	(131-150)	(151-170)	(171-195)	(196-220)	(221-255)	(256-275)	(lbs)	
1	2	3	4	5	6	7	8		Foot Spring Set

Note:

If in doubt choosing between two categories, choose the higher rate spring set.

Foot spring set recommendations shown are for transtibial users.

For transfemoral users we suggest selecting a foot spring set one category lower, refer to Section 16 *Fitting Advice* to ensure satisfactory function and range of movement.

Checking the package contents:

- 1. ElanIC
- 2. Inductive Battery Charger Kit (Part No. 409087IC)
- 3. Clinician's Manual
- 4. User Guide
- 5. External USB *Bluetooth*® Adapter
- 6. USB stick with interface software

Also available for device set-up:

- (supplied separately)
- 1. USB stick and interface software
 - 2. External USB *Bluetooth*® Adapter
 - 3. Alignment wedge
 - 4. Inductive Battery Charger Kit (Part No. 409087IC)
 - 5. Blatchford Programming Tablet
 - 6. Mains adaptors for different countries

2 Safety Information

-  This warning symbol highlights important safety information which must be followed carefully.
-   The device and its charger must not be used by users of pacemakers or implantable cardioverter-defibrillators (ICDs). Contact your service provider for advice if the device is to be used with any other electronic medical device.
-  Any changes in performance of the ankle motion e.g. restricted movement, non-smooth motion should be immediately reported to your service provider.
-  Always use a hand rail when descending stairs and at any other time if available.
-  Any excessive changes in heel height after programming will adversely affect ankle function and should be immediately reported to your service provider for reprogramming and calibration.
-  After continuous use or charging the ankle casing may become hot to the touch.
-  Avoid strong magnetic fields, sources of electrical interference.
-  Avoid exposure to extreme heat and/or cold.
-  The device is sealed for water resistance. However the water resistance of the device is limited by environmental factors and exposure time as detailed in Section 6 of this manual.
-  The inductive battery charger kit is not waterproof. Do not charge in atmospheres containing liquids and/or powders.
-  The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling is acceptable.
-  Assembly, programming and repair of the device must only be carried out by suitably qualified practitioners, certified by Blatchford.
-  The user must not adjust or tamper with the setup of the device.
-  The user should be advised to contact their clinician if their condition changes.
-  The device detects standing and walking on inclines, declines and walking at various speeds on normal ground. Non-identified movements such as treadmills and travelators may cause unexpected ankle behavior.
-  Never plug the charger into the ElanIC while still wearing the limb.
-  Only charge the batteries in the ElanIC using the charger supplied with the product. Do not use this charger for charging any other device.
-  Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.
-  To minimise the risk of slipping and tripping, appropriate footwear that fits securely onto the footshell must be used at all times.
-  Be aware of finger trap hazard at all times.
-  Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
-  Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer; otherwise, degradation of the performance of this equipment could result.


See further warnings relating to the charging of the device in Section 7.3

3 Construction

Principal Parts

- Hydraulic Body Assembly including pyramid (aluminum/stainless steel/titanium incorporating batteries (NiMH), electric motors and PCB)
- Carrier Assembly (aluminum/stainless steel)
- Heel & Toe Springs (e-Carbon)
- Spring Attachment Screws (titanium/stainless steel)
- Glide Sock (UHMW PE (Ultra-high Molecular Weight Polyethylene))
- Foot Shell (PU) (Polyurethane)
- Charging Surface & Charger Pod (thermoplastic)

Component Identification

4 Function

With microprocessor control the device adapts dynamically (via motorized valves) the ankle plantarflexion and dorsiflexion resistance settings to match the user's walking mode and terrain.

Activity	Device Mode	Hydraulic Resistance
Walking (level ground)	Additional ground clearance in swing phase.	Slightly dorsiflexed Basic resistance settings
Walking fast	<i>Fast Walk</i> <i>[Assist]</i> Added momentum in the 'step to step' transition of the gait cycle (<i>transtibial only</i>).	Increased plantarflexion resistance, decreased dorsiflexion resistance
Walking uphill (moderate/steep incline)	<i>Ramp Assist</i> Enables easier walking up inclines.	Progressively changes; increasing plantarflexion and decreasing dorsiflexion resistance
Walking downhill	<i>Ramp Brake</i> Helps forward momentum and provides greater stability and security	Progressively changes; decreasing plantarflexion resistance and increasing dorsiflexion resistance
Standing	Basic Resistance settings or <i>Standing Support*</i> assistance	User settings from set-up/ high plantarflexion and dorsiflexion resistance
Low battery power	Sleep	Reverts to Basic Resistance settings from setup (No dynamic control)

Note:

Standing Support mode provides trip and fall risk reduction and prevention features.

*This feature is only available when enabled in the software setup and the battery has sufficient charge.

To enable Standing Support mode see Section 15 *Options> Standing Support Mode*.

5 Maintenance

Maintenance must be carried out by competent personnel.

A service is required every 20 months as a condition of warranty.

It is recommended that the following maintenance is carried out annually:

- Remove the foot shell and glide sock, check for damage, wear or corrosion and replace if necessary.
- Check all screws for tightness or corrosion, clean and reassemble as necessary.
- Visually check the heel and toe springs for signs of delamination or wear and replace if necessary. Some surface damage may occur after a period of use; this does not affect the function or strength of the foot.

(See Section 17 for the foot shell removal and spring replacement procedure.)

The wearer should be advised:

Any changes in performance of this device must be reported to the practitioner.

Changes in performance may include:

- Increase in ankle stiffness
- Reduced ankle support (free movement)
- Any unusual noise
- Recurring rapid flashing of Charging Indicator/Status Indicator LEDs.

The user should be advised that a regular visual check of the foot is recommended, signs of wear that may affect function should be reported to their service provider (e.g. significant wear or excessive discoloration from long term exposure to UV).

Cleaning

Use a damp cloth and mild soap to clean outside surfaces, DO NOT use aggressive cleansers.

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

This device is IP67 rated: waterproof to a maximum depth of 1 meter. Occasional submersion should be limited to a max. of 30 minutes. Dry after using in water.

Avoid exposing the device to abrasive or corrosive environments, such as those containing sand, acids, salt, chlorinated, or chemically treated water e.g. swimming pools. Failure to comply will void the warranty. After any accidental exposure to such environments rinse the device in fresh water and dry thoroughly.

Exclusively for use between -15°C and 50°C
(5°F and 122°F).

Charge only between 0°C to 35°C (32°F to 95°F).

Suitable for submersion

7 Battery Charging



Please advise user of the charging information in this section.
The batteries are not field replaceable and will be replaced as necessary when serviced.

Read the user instructions supplied with the charger before using for the first time.
We recommend the device is charged daily; it should take 5.5 hours to charge fully, depending on residual charge. An exhausted or flat battery may take longer.

7.1 Setting Up

Lay the limb on a flat surface and connect the inductive charger pod to the charging surface on the back of the device.
Use the Charging Indicator LED on the right hand side of the device to confirm that charging has initiated correctly:

 **Do not allow anything (e.g. stocking, clothing) to intrude in between the charging surface and the charger pod while charging is in progress. This will result in decreased charging performance and buildup of heat.**

7.2 Charging

One initial flash- x1 — Connection Established

Followed By

	LED	Charging state
Single Flash Repeated	 —  —  — 	Charging - Battery Empty
Double Flash Repeated	  —   —  	Charging - Charge Level Medium
Triple Flash Repeated	   —   	Charging - Charge Level High
Steady ON/Steady OFF	 / —	Charging Completed
Rapid Flashing	       	High Temperature Shut Down*


*Charging continues when normal charging temperature is reached. (See Section 6.)

Charging Error

Charging error is always indicated by a non-responding Charging Indicator LED on the right hand side of the device (no initial flash or charging state indication).

Further information is provided by the charger pod LED (blue):

LED	Error	Corrective Action
Slowly Pulsating	No Error/Normal charging	Not needed.
Steady Blue Light	Charger operational but not mated correctly to the charging surface of the device.	Try reattaching and check if anything prevents perfect mating of charging surfaces. Remove any intruding material.
No Light	The charger is not powered.	Check that charger is connected to the mains.
Rapid Pulsating	Possible overheating of the charger.	Charging resumes when normal charging temperature is reached. If recurring contact your service provider.

Reset

The device will take several seconds to reset when the charger is disconnected and a single blue flash from the Status Indicator LED on the left hand side of the device will confirm a correct startup.

7.3 Warnings


-   The device and its charger must not be used by users of pacemakers or implantable cardioverter-defibrillators (ICDs). Contact your service provider for advice if the device is to be used with any other electronic medical device.
-  Only charge the batteries in the device using the charger supplied with the product. Do not use this charger for charging any other device.
-  The battery charger can only be used with an electrical outlet that is rated between 100 V and 240 V AC, 50/60 Hz.
-  Never connect the charger to the device while still wearing the limb.
-  During charging no power is delivered to the ankle control system.
-  The ankle and the charger may become warm during charging but should never be too hot to touch. There should be no abnormal smell. If this is the case immediately unplug from electrical supply and contact your service provider.
-  Do not refit the limb until the charger has been disconnected.
-  Always recharge the device before use if it has not been used for an extended period of time.
-  The inductive battery charger kit is not waterproof. Do not charge in atmospheres containing liquids and/or powders.
-  Do not exceed the charging temperature limitations set out in this manual.
-  When charging an empty battery the charger may terminate prematurely (after a few minutes). If this happens, reconnect the charger. If this problem recurs contact your service provider.
-  Do not leave the charger plugged to electrical supply after use.

7.4 Battery Level Indication

To check the battery charge level, tilt the foot down vertically and hold stationary for 4 seconds. The Status Indicator LED on the left hand side of the device will flash blue, 1-3 times (see table) and then repeat to indicate the charge status. Tilt the foot again to repeat the battery status check.

Note: This feature can be switched off, see Section 15 > Options.

Note: After charging the device fully the Status Indicator LED may indicate 2 or 3 flashes. This will depend on the charging time, the age and condition of the battery.

LED	Battery Charge
3 flashes	Maximum
2 flashes	High
1 flash	Low

7.5 Device Status Indication

The same LED on the left hand side of the device is also used to indicate device status:


1. Single Flash
When the charger is disconnected it will flash once to indicate correct start up.
2. Steady ON
If after the charger is disconnected it stays on it indicates that an error was detected by your device.
3. Flashing rapidly
If it flashes rapidly after charger is disconnected, it indicates that Bluetooth® is active - flashing should go away after a few minutes.

Note: The device Status Indicator LED is also used for indicating stages of a **Bluetooth®** connection. (See Section 10.)

See Appendix 1 for a full summary of all possible LED indications of the device and the charger pod.

8 Bench Alignment


Keep the build line between pivots as shown, using shift and/or tilt devices as necessary.

Align transfemoral devices according to the instructions for use that are supplied with the knee.

Tilt Setting

Align limb to achieve range of motion shown.

Align with shoe on and foot fully plantar flexed.

*Approximate ratio

9 Install the Blatchford PC Software Interface

Note: When installing a newer version of the interface first uninstall any older versions.

To install the Blatchford PC Software Interface on a computer running Microsoft Windows (see System Requirements below) follow this procedure:

1. Put the PC software USB stick in a spare USB port.
2. In the folder **EN** double click **setup.exe** the **Blatchford Elan Setup Wizard** screen will display click **[Next]** and follow the on-screen message prompts to continue the installation.
3. When **Installation Complete** displays click **[Close]** to exit the install program.

System Requirements

Bluetooth® 2.0 or later

Supported Operating Systems:

- Windows 10
- Windows 8
- Windows 7
- SP1 Windows 7

Supported Architectures: x86, x64

Hardware Requirements

Recommended Minimum: Intel Dual Core Duo

1 GHz or higher with 512 MB RAM or more



Minimum disk space: 100 MB

Prerequisites

Windows Installer 3.1 or later

Internet Explorer 5.01 or later

i Note: This is a generic Elan interface program. The images shown may differ from your device. This has no effect on compatibility.

i Note: The device can also be programmed using an Android App.

See 938447, *Elan Android Programming App Instruction for Use*.

10 Using the PC Software Interface to set up the Device

Getting Started

The device is a biomimetic foot and ankle system with microprocessor control for plantar and dorsiflexion resistance. This device is designed to adapt dynamically to provide assistance when walking and standing on a variety of surfaces and changing gait speeds.


To set up the ankle resistance setting to an individual user preference the device control system is adjusted and calibrated with the PC interface via a **Bluetooth®** link.

The final walking calibration routine is used by the computer program in the device to adapt to gait characteristics after the alignment and resistance settings have been optimized.


This calibration data is then used by the device software to recognize and adapt the device automatically to changing walking and daily living patterns.

10.1 Set-up Procedure

1. Ensure the batteries are charged (see Section 7).
2. Complete the Bench Alignment for the device (see Section 8).
3. Install the PC software interface (see Section 9) on a Microsoft Windows computer or remote device with a **Bluetooth®** capability.
4. Ensure **Bluetooth®** is enabled on the computer/remote device and disconnected from any other device.
5. To activate the built-in **Bluetooth®** capability of the device, fully dorsiflex the foot then tilt it to approximately 45°, as if performing a battery level check.
6. Hold the device in this tilted position for approximately 4 seconds while the left LED flashes the battery charge status indication (1, 2, or 3 short flashes) twice*.
7. Maintain the tilted position. After a further 5 seconds the left LED illuminates steady ON.
8. Return device to standing position.

Fully dorsiflex the foot.

Tilt and then wait to activate **Bluetooth®**.

Return to standing position to connect/maintain **Bluetooth®** connection.


Note: If the battery charge status indication function has been disabled you will not see two sets of battery status flashes.

9. The left LED now flashes continuously for 7 mins. During this period **Bluetooth®** is on and the device is discoverable.
10. Pair the device to your computer/remote device using the **Bluetooth®** manager program on your computer.


Bluetooth® Device Name:

Device Serial Number:

ELAN-56A7
SN 123456


11. If prompted during installation to enter a pair code, enter pair code > **1234**. Or you may be asked when connected to confirm the ID code of the built-in **Bluetooth®** device. If so, enter > Yes. (The Built-in **Bluetooth®** device name can be found on the label on the back of the device.)
12. Carry out the device Alignment, Set-up and Calibration procedure using the PC interface controls with the user walking on a horizontal and level surface. (See procedure on the following pages).

Bluetooth® Range

The **Bluetooth®** range of operation may vary depending on the remote device used, we recommend a range of 5-10 meters distance. If the connection terminates unexpectedly and the Status Indicator LED changes from steady (connected) to flashing blue (standby), reduce the distance between the ankle and the remote device being used until the connection is restored and can be maintained.

Steady Blue > Connected

Flashing Blue > Standby

10.2 Launch the PC Interface Program

1. To launch the PC Interface program:
Click the Elan program icon in the list of programs **or** from your desktop.

2. When the **Login** screen appears

Click on:

B


to continue to the **Start** screen.

Before setting up the device ensure that the battery is charged, see Section 7.

A

Name and user **Password** are only required where access is restricted.

10.3 Bluetooth® Connection [Start]

- 1. Ensure the device is charged up and the **Bluetooth®** is activated on the device (flashing LED).
- 2. Select the paired port for the device from the pull down menu. **A**


If the device has been used previously you can use **[Auto Connect]**. **B**

When connected, the Status Indicator LED will change to steady and the **Bluetooth®** symbol on the status bar will change to active.

i To return to this screen click on:

C **Start**

i If for **15 minutes** no command is sent from the PC Interface Software to the device the connection will time out automatically.

Select auto connect or communication port to start

D **E**

Auto Connect **Disconnect** **System Reset**

Connected

[Disconnect] **D**

To disconnect the device from the PC/Laptop.

[System Reset] **E**

Resets the software in the device if the **Bluetooth®** connection is interrupted.

10.4 Status Bar

The current status of the device is shown by the indicators on the Status bar.

Mode Indicator	Description
Mode	Displays the current device mode
Angle	Shows foot angle dynamically
Tilt	Shows tilt position dynamically
PF and DF Valve	Displays the current dynamic resistance settings for plantarflexion and dorsiflexion
Battery 	Shows the device is powered up
Power Save/ Sleep 	If displayed the device is in Power Save mode*
Standing Support** 	If displayed Standing Support Mode is active, high ankle resistance is selected.
Bluetooth® 	Shows that the device is connected/disconnected via Bluetooth®

*Power Save/Sleep

After Calibration if the device detects no movement for approximately 10 seconds it activates power saving mode disabling the sensors. Note during the calibration and set-up the time interval is increased to 90 seconds.

**Standing Support Mode

This mode is only available if enabled, the default setting is OFF, see Section 15 *Options*. When the device detects a stationary standing posture high ankle resistance is selected.

To adjust or fine tune the setting on a previously calibrated device:

If the alignment has not been changed it is not necessary to follow the whole set up procedure again. Follow the **Walking Calibration** sequence (see Section 14) and the device software will then add the new information to the previous data and average out the result giving improved settings for the user.

If alignment adjustments are required carry out a full factory reset and follow the full set up procedure.


11 Static Alignment

1. Complete the Bench Alignment set-up (see Section 8).
2. Click on

A

Alignment

to select the *Dynamic Alignment Vector* screen.

11.1 Static Alignment Procedure

1. With the user standing comfortably adjust the device's A-P shift to achieve a neutral balance position*.

See following pages.

Tilt B


2. Use the display to adjust the angular position of the device when standing.
3. Adjust using shift and/or tilt devices as required.

(When standing the *Vector* indicator position shown is recommended)

WARNING: Ensure the user is not resting on the DF (Dorsiflexion) limit when standing in an upright position.

Tilt Setting Position


*Ensure the device is uncalibrated prior to static alignment. See Section 14/*Factory Reset*.

12 Biomimetic Alignment

The aim of alignment is to achieve a *balance point* while standing and set the hydraulically damped range of motion. The aim of damping adjustment is to fine tune the ankle-foot roll-over stiffness characteristics until a comfortable gait is achieved.


Due to the increased range of motion provided by the ankle the user may experience the need for more voluntary control and initially find the ankle disconcerting during setup. This should quickly pass upon completion of satisfactory setup.

12.1 Biomimetic Adjustment

Note: Carry out static alignment while ensuring the user has some means of support such as parallel bars. This is standing alignment only.

Check the Bench Alignment and the heel height.


Use shift for static alignment and standing. The device should encourage some degree of self adjustment to achieve a sense of balance for the user during standing.

13 Dynamic Alignment

Click on:

A Resistance Setup

to select the *Basic Resistance Settings* screen.

13.1 Dynamic Alignment Procedure


1. Instruct the user to walk for approximately 3 minutes to become familiar with the resistance settings.
2. Adjust the Plantarflexion and Dorsiflexion hydraulic resistance settings to optimize the device to suit the user's preference for normal walking. (Note the scale is from 4 to 7.)
3. Allow some time for practice and familiarization with the ankle resistance settings.
4. Verify that the hydraulics are moving when walking using the Vector Indicator. Reduce the resistance settings if minimal movement is shown.

Power Save/Sleep


If the device detects no movement for approximately 90 seconds it activates power saving mode, the Power Save symbol appears on the Status Bar and all the indicators will stop moving. Moving the device will automatically bring it back to normal operation.

13.2 Resistance Setting Troubleshooting


14 Walking Calibration

Click on:

to select the *Walking Calibration* screen.

i Calibration may be repeated at any time to fine tune the set-up and is recommended after any change in alignment and/or basic resistance settings.

14.1 Walking Calibration Procedure

1. Move the slider to select the amputation level,

B Transfemoral *TF* or Transtibial *TT*

2. Start calibration.

2.1 To calibrate the device to the user's individual gait characteristics instruct the user to walk on a horizontal, level surface and at a comfortable pace.

2.2 After several steps, click on:

D The step counter counts down the number of remaining steps required.

(Note: 8-10 uninterrupted, regular steps on level ground are required for calibration).

i The device adapts to the individual user's gait characteristics from the Calibration stages. (See Section 4 for more information)

Factory Reset

Should it be necessary to begin a completely new set-up for the device,

Caution:
 Clicking on [**Factory Reset**] will *delete* all the Dynamic Alignment adjustment and any walking calibration data stored in the device.

Click on

If this option is chosen return to Section 13 to begin set-up.

15 Options

Click on:

to select the *Options* screen.

Battery Indicator

Shows battery charge level with a flashing blue LED when foot is plantarflexed.

Standing Support Mode

Default is OFF

Increases resistance to reduce fatigue when standing for long periods.

Standing mode is automatically disabled when battery charge becomes low to conserve remaining power.

Fast Walk Mode (TT users only)

Enhances resistance settings to assist fast walking.

Ramp Assist Mode

Enhances resistance settings to assist ascending a ramp/incline.

Ramp Brake Mode

Enhances resistance settings with a braking effect when descending a ramp/incline.

Mode Default Settings

Mode	TT	TF
Battery	ON	ON
Standing Support	OFF	OFF
Fast Walk	ON	n/a
Ramp Assist	ON	ON
Ramp Brake	ON	ON

Note:

Disabling [OFF] the Fast Walk and Ramp Modes will disable the response to certain conditions.

Fast Walk Mode is not recommended for TF users.

16 Fitting Advice

The correct alignment (A-P position), range of motion (distribution of plantar to dorsiflexion) and adjustment of the hydraulic resistance settings are critical in achieving a smooth roll-over and correct gradient adaptation.

The device is supplied assembled with heel and toe springs of the same category. If after following the instructions below you have problems with the function please contact the sales team in your area for advice.

Any of the following will have a negative effect on function and stability.

- Incorrect spring selection
- Incorrect A-P shift alignment
- Incorrect distribution of plantar and dorsiflexion range

Prior to setup or if further adjustment is required ensure the ankle is not calibrated by carrying out a factory reset and that standing support mode is disabled. See Section 14 for details of factory reset and Section 15 to switch off standing support mode.

	Symptoms	Remedy
1	• Sinking at heel strike • Difficulty in achieving a smooth progression to mid stance • User feels they are walking up hill or forefoot feels excessively long	• Increase plantarflexion resistance • Check A-P shift alignment; ensure foot is not too anteriorly positioned • Check distribution of plantar and dorsiflexion movement; ensure that the plantarflexion range is not excessive • Check spring category is not too soft, if so fit a higher rate spring
2	• Progression from heel strike to mid stance is too rapid • Difficulty in controlling the energy return from the foot at the heel strike (reduced knee stability) • User feels heel is too hard, fore foot is too short	• Reduce plantarflexion resistance • Check A-P shift alignment; ensure foot is not too posteriorly positioned • Check distribution of plantar and dorsiflexion movement; ensure that there is adequate plantarflexion range • Check the spring category is not too high for the weight and activity of the patient, if so fit lower rate spring
3	Heel contact and progression feel OK but: • Forefoot feels too soft • Forefoot feels too short • User feels they are walking down hill, possibly with reduced knee stability • Lack of energy return	• Increase dorsiflexion resistance • Check A-P shift alignment; ensure foot is not too posteriorly positioned • Check distribution of plantar and dorsiflexion movement; ensure that there is not excessive dorsiflexion range • Check the spring category is not too soft for the weight and activity of the patient, if so fit higher rate spring

	Symptoms:	Remedy
4	• Forefoot feels too rigid • Forefoot feels too long • Feels like walking up hill	• Reduce dorsiflexion resistance • Check A-P shift alignment; ensure foot is not too anteriorly positioned • Check distribution of plantar and dorsiflexion movement; ensure that there is sufficient dorsiflexion range • Check the spring category is not too rigid for the weight and activity of the patient, if so fit lower rate spring

17 Maintenance Instructions

- 

Be aware of finger trap hazard at all times.
- 

Do not lever against the charging surface. Use appropriate tooling and protective materials (e.g. towel, cloth) to prevent damage.
- 


Use appropriate health and safety equipment at all times including extraction facilities.

17.1 Footshell Removal

1
Insert shoe horn behind heel spring.

2
Rotate shoe horn as shown to remove shell.

17.2 Spring Replacement

3

Remove heel spring and screws.


4

Remove toe spring screw, replace toe. Upon reassembly, use Loctite 243 (926012) and torque to 35 Nm. Ensure toe spring is central to the carrier.

5

Re-assemble with replacement heel spring. Use Loctite 243 (926012) and torque to 15 Nm.

6

Cover appropriate lines on carrier with permanent black marker to leave spring set number showing.

7

If a foam cosmesis is to be fitted, roughen top surface of foot shell to provide ideal bonding surface.

8


Lubricate toe and heel if required. (Footshell is pre-lubricated).
Fit sock as shown.


- 9** Slide carrier/heel spring assembly into the foot shell.


- 10** Toe spring location in foot shell

- 11** Use a suitable lever to encourage the heel spring into location in the foot shell.

- 12** Ensure heel spring is engaged into slot.

Heel spring location slot

13


18 Technical Data

Operating and Storage Temperature Range:	-15°C to 50°C (5°F to 122°F)
Charging Temperature Range:	0°C to 35°C (32°F to 95°F)
Component Weight [Size 26N]:	1.24 kg (2 lb 12 ozs)
Recommended Activity Level:	3
Maximum User Weight:	125 kg (275 lb)
Proximal Alignment Attachment:	Male Pyramid (Blatchford)
Range of Hydraulic Ankle Motion (excludes additional range of motion provided by heel and toe springs)	6 degrees plantarflexion to 3 degrees dorsiflexion
Build Height:	[sizes 22-26] 172 mm
[See diagram below]	[sizes 27-30] 177 mm

Battery:	
Battery Type	Rechargeable
Battery Life	Variable , but at least 24 hrs from full charge. Charge daily.
Battery Charger:	
Charging time to full charge	5.5 hours
Input Voltage	100V to 240V AC
Frequency	50 to 60 Hz

Fitting Length

19 Ordering Information

Order Example

ELANIC	25	L	N	3	S
	Size	Side (L/R)	Width* (N/W)	Spring Set Category	Sandal Toe

Available from size 22 to size 30:

ELANIC22L1S to ELANIC30R8S
ELANIC22L1SD to ELANIC30R8SD
(add 'D' for a dark tone foot shell)

*Sizes 25-28 only. For all other sizes, omit the Width field.

e.g. ELANIC25LN3S, ELANIC22R4S, ELANIC27RW4SD

19.1 Replacement Parts

Spring Kits				
Rate	Foot Sizes			
	Small (S)	Medium (M)	Large (L)	Extra Large (XL)
	22–24	25–26	27–28	29–30
Set 1	539801S	539810S	539819S	539828S
Set 2	539802S	539811S	539820S	539829S
Set 3	539803S	539812S	539821S	539830S
Set 4	539804S	539813S	539822S	539831S
Set 5	539805S	539814S	539823S	539832S
Set 6	539806S	539815S	539824S	539833S
Set 7	539807S	539816S	539825S	539834S
Set 8	539808S	539817S	539826S	539835S



Foot Shell (for dark add 'D')			Item	Part. No.
Size/Side	Narrow	Wide		
22L	539038S	-	Glide Sock (Sizes 22–26)	531011
22R	539039S	-	Glide Sock (Sizes 27–30)	532811
23L	539040S	-		
23R	539041S	-	Inductive Battery Charger Kit	409087IC
24L	539042S	-	–ElanIC Mains Adaptor - UK	950185
24R	539043S	-	–ElanIC Mains Adaptor - US	950184
25L	539044SN	539044SW	–ElanIC Mains Adaptor - Europe	950183
25R	539045SN	539045SW		
26L	539046SN	539046SW	ElanIC Mains Adaptor - Australia	950186
26R	539047SN	539047SW		
27L	539048SN	539048SW	External USB <i>Bluetooth</i> ® Adapter	950129
27R	539049SN	539049SW	USB Stick with Interface Software	406413E
28L	539050SN	539050SW		
28R	539051SN	539051SW	Blatchford Programming Tablet	019179
29L	-	539052S		
29R	-	539053S		
30L	-	539054S	Alignment Wedge	940093
30R	-	539055S		

Appendix 1 Quick Guide to LED Indications

Device Status Indicator LED (Left Hand Side)


Note: Take into consideration that the battery level indication feature may be disabled.

Note: The built-in **Bluetooth®** device can be activated only when the ankle is in a fully dorsiflexed position.

Tilted Position

LED	Indicating	Action
Single Flash x2	Battery level low	Charge to use
Double Flash x2	Battery level high	Charge for maximum performance
Triple Flash x2	Battery level high	Ready for use
Steady ON after battery level indication (foot dorsiflexed)	Bluetooth® is ready for turning ON	Turn upright to turn Bluetooth® ON
Steady ON without battery level indication (foot dorsiflexed)	Bluetooth® is ready for turning ON (Battery level indication Disabled)	Turn upright to turn Bluetooth® ON


Upright Position

LED	Indicating	Action
Single flash (After disconnecting charger)	Correct start up after charging	None
Steady ON (After disconnecting charger)	Error detected	Contact service provider
Continuous flashing after disconnecting charger	Bluetooth® is ready for connection	Wait few minutes for deactivation
Continuous flashing after Bluetooth® turned On	Bluetooth® is ready for connection	Connect PC software or Wait for deactivation
Steady ON after Bluetooth® connection initiated	Bluetooth® connection to PC software successful	Use PC software
Continuous flashing after Bluetooth® connection established	Bluetooth® connection is broken Bluetooth® is ready for connection	Try reestablish connection Try reducing distance between devices

Charging Indicator LED (Right Hand Side)

Note: Always interpret the indications given by the Charging Indicator LED and the Charger Pod LED in connection with each other.


Note: The charger will turn itself off after 6 hours of operation. This will manifest itself by the Charger Pod LED switching from pulsating to steady ON.

LED	Indicating	Action
Start Up		
OFF	Charger connection failed	Try reconnect
Single flash	Correct start up after charger connected	See further indications below
Single flash + LED stays OFF	Correct start up after charger connected	Battery full
Further Indications After Correct Start Up		
Single flash repeated	Charging empty battery	None
Double flash repeated	Charging medium level battery	None
Triple flash repeated	Charging high level battery	None
Rapid flashing	High temperature shut down	Charging resumes after cooling
End of Charging		
Steady ON	Charging completed	Ready for use

Charger Pod LED

LED	Indicating	Action
Slow pulsating	Normal Operation of charger	None
Steady ON	Operational but no charging	Check mating to charging surface
OFF	Not powered	Check mains connection
Rapid pulsating	Possible overheating	Charging may resume after cooling down

Manufacturer’s Declarations and Legal Information

Electromagnetic Environment

The device is designed for operation in the following electromagnetic environments:

- Operation in a professional healthcare facility (e.g. hospital)
- Operation in areas of home healthcare (e.g. use at home, use outdoors)

Electromagnetic Emissions

Standard	Test	Test Level Compliance
IEC 60601-1-2: 2014	Radiated Emissions	Group 1, Class B
	Conducted Power Line Emissions	Group 1, Class B
	Harmonic Current Emissions	Exempt
	Voltage and Fluctuations & Flicker	Exempt
FCC CFR47:2018 part 15 (b)	Radiated Emissions - Class B	Group 1, Class B
	Conducted Emissions - Class B	Group 1, Class B

Electromagnetic Interference Immunity

Standard	Test	Test Level Compliance
IEC 60601-1-2: 2014	ESD (Electrostatic discharge)	±2 kV, ±4 kV, ±8 kV, ±15 kV Air ±2 kV, ±4 kV, ±6 kV, ±15 kV Contact
	Radiated Immunity	80-2700 MHz, 10V/m
	Electrical Fast Transient/Burst	+500V, –500V, +1 kV, –1 kV, ±2kV
	Surge	+500V, –500V, +1 kV, –1 kV
	Power Frequency Magnetic Field	3V rms 6V rms (ISM And Amateur Radio bands)
	Voltage Dips and Short Interruptions	100% for 10 ms / 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 100% for 20 ms / 1 cycle at 0° 30% for 500 ms / 25 cycles at 0° 100% for 5000 ms

Bluetooth® Module

Manufacturer’s Part Number	RN4678-V
FCC ID	A8TBM78ABCDEFGH
Frequency	2.4 GHz
Effective Radiated Power	Class 2: +1.5 dBm Typical

FCC Compliance Statement

This device complies with part 15 of the FCC Rules.

Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

Note... Contains pre-approved **Bluetooth®** module: FCC ID: A8TBM78ABCDEF GH.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation. This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party: Blatchford Inc.

1031 Byers Road, Miamisburg, OH 45342
800-548-3534 (toll free) | 937-291-3636
customerservice@blatchfordus.com

RSS Compliance Statement

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

1. this device may not cause interference, and
2. this device must accept any interference, including interference that may cause undesired operation of this device.

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Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website Blatchford Inc. 1031 Byers Road, Miamisburg, OH 45342.

800-548-3534 (toll free) | 937-291-3636 customerservice@blatchfordus.com

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchfordmobility.com

Medical device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Warranty

The device is warranted for 36 months—foot shell 12 months—glide sock 3 months.

See Blatchford website for the current full warranty statement.

Scheduled service every 20 months is required, for maintenance of the warranty.

For extended warranty contact service provider.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licences and exemptions.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

This symbol indicates that the product contains electrical/electronic components and/or batteries that should not be disposed of in general waste or be incinerated at the end of the product's life.

At the end of the product's life, all electrical/electronic components and/or batteries should be recycled or disposed of in accordance with the current regulations for the handling of WEEE (Waste Electrical and Electronic Equipment), or equivalent local regulations. The remainder of the product should also be recycled where possible in accordance with local waste recycling regulations.

To help prevent potential harm to the environment or to human health from uncontrolled waste disposal, Blatchford offers a take-back service. Please contact Customer Services for details.

Retaining the Packaging Label

You are advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgments

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG Inc., and any use of such marks by Blatchford is under license. Other trademarks and trade names are those of their respective owners.

Elan, ElanIC and Blatchford are registered trademarks of Blatchford Products Limited.

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© Blatchford Products Limited 2025. All rights reserved.
Blatchford Europe GmbH
Blatchford Products Ltd.