# Silcare Breathe Active Locking Liner

## Instructions for Use

SBATTLP22-SBATTLP40

|    |                            |     |
|----|----------------------------|-----|
| EN | Instructions for Use       | 2   |
| FR | Instructions d'utilisation | 16  |
| DE | Gebrauchsanweisung         | 30  |
| IT | Istruzioni per l'uso       | 44  |
| ES | Instrucciones de uso       | 58  |
| NL | Gebruiksaanwijzing         | 72  |
| PL | Instrukcje użytkowania     | 86  |
| PT | Instruções de utilização   | 100 |
| CS | Návod k použití            | 114 |

Contents .....2

1 Description and Intended Purpose .....3

2 Safety Information .....5

3 Construction .....6

4 Function.....6

5 Maintenance .....6

    5.1 Cleaning the Device .....7

    5.2 Cleaning the Valve .....7

    5.3 Cleaning the Residual Limb.....7

6 Limitations on Use .....8

7 Choosing the Correct Size.....8

8 Casting/Scanning the Socket with the Casting Dummy.....9

9 Trimming the Device..... 10

10 Donning the Device ..... 10

11 Fitting the Valve ..... 11

12 Fitting the Spacer..... 11

13 Fitting Advice..... 12

    13.1 Moisture in the Distal End of the Device ..... 12

    13.2 Loss of Vacuum ..... 12

    13.3 Separation of the Valve from the Distal Cap..... 12

14 Technical Data ..... 13

15 Ordering Information ..... 14

---

# 1 Description and Intended Purpose

These instructions are for the practitioner.

Please read these instructions carefully before fitting the device.

The term *device* is used throughout these instructions for use to refer to Silcare Breathe Active Locking Liner.

Ensure that the user understands all instructions for use, drawing particular attention to the maintenance and safety information.

## Application

This device is an interface component for use only as a component of a lower limb prosthesis. Intended for single user.

It is intended for moderate to high activity users with well covered residual limbs.

This device is a biocompatible, perforated locking liner that allows air and moisture to escape through its perforations during stance. During the gait cycle, a silicone valve provides vacuum suspension by controlling the airflow through the perforations on the distal cap.

To optimize comfort ensure that the user is instructed in the correct way to handle and don/doff the device. See Section 10 *Donning the Device*. Also ensure that they are made aware of how to maintain the device and to keep it clean for maintaining hygiene as set out in these instructions. See 5 *Maintenance* section.

## Activity Level

The device is recommended for Activity Level 3 and 4 users, although may be suitable for all four activity levels.

The silicone used in the Silcare Walk liners is softer and more flexible than that used in Silcare Active liners. The softer silicone offers more comfort and flexibility in all directions but may feel less connected or secure at higher activity levels. Silcare Walk liners are therefore generally recommended for the less active users (Activity Level 1 and 2).

The silicone used in Silcare Active liners is firmer and slightly less flexible, offering a more secure interface at higher activity levels. Silcare Active liners are therefore generally recommended for the more active users (Activity Level 3 and 4).

Note that some high activity users with poor residual limb tissue coverage and/or high sensitivity may prefer the comfort of the softer silicone Silcare Walk liner to the firmer Silcare Active liner. Conversely, some low activity users with good residual limb tissue coverage and/or low sensitivity may prefer the security of the firmer Silcare Active liner to the softer Silcare Walk liner.

### Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

### Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

### **Activity Level 3**

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

### **Activity Level 4**

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

## **Contraindications**

- Deep scars on the distal end of residual limbs can cause failure of the vacuum suspension. These deep scars reduce skin contact with the silicone inside the device and make a less effective seal
- Conical residual limbs can cause failure of the vacuum suspension. Conical residual limbs reduce skin contact with the silicone inside the device and make a less effective seal
- Short residual limbs can cause failure of the silicone side of the device. After donning, if the unperforated fabric is on or above the patella tendon, the silicone is subject to more stress during knee flexion
- Users with poor hand or cognitive function might find donning and cleaning difficult
- Poor hygiene

## **Clinical Benefits**

- Provides cushioning for the residual limb in the socket
- Distributes in-socket pressure more evenly, compared to other materials and alternative cushioning solutions
- Improvements in residual limb health problems and wound healing, compared to non-perforated liners
- Improved heat dissipation, compared to other temperature regulation solutions
- Removes sweat from skin interface
- Patients reported a preference for their perforated liners, compared to non-perforated liners
- Reduces the need to remove prosthesis throughout the day to dry residual limb, compared to non-perforated liners
- The locking liner provides a means of suspension

## 2 Safety Information

This warning symbol highlights important safety information.

The user should be advised to contact their practitioner if their condition changes.

Any deterioration in residuum condition or any change in sensation should be reported to the practitioner. If the residual limb has tissue damage, the use of the device should be stopped.

Make sure that any damaged skin or open wound is properly and suitably dressed to prevent direct contact with the device.

Wearers with sensitive skin, diabetics and vascular cases should be extra vigilant and may need to apply lubricant to sensitive areas. We recommend a routine visual check and if required the user should consult their medical practitioner.

For other medical conditions the user should follow the advice and recommendation of a physician or a medical practitioner regarding skin care.

Enlarged perforations can trap the skin and cause blisters. If the perforations enlarge, stop using the device.

If distal oedematous swellings corresponding to distal perforations in the liner occur, use of the liner should be discontinued and the swellings reported to the practitioner.

**Do not** use alcohol sprays, household cleaners, or abrasives. These cleaning materials could damage the device and irritate skin.

**Do not** pull or stretch the fabric. Fingernails, sharp jewelry, and the locking pin can tear the fabric. If the fabric is torn, stop using it and contact a Blatchford sales representative.

Sockets with sharp proximal edges can tear the device.

Take care when handling the device to avoid possible contamination with materials such as fibreglass which will stick to the device and cause skin irritation.

When donning a sock, clothing and the prosthetic limb, be aware that the device can build up a static charge.

To avoid danger of suffocation, keep the device away from babies and children.

Keep the device away from direct heat sources.

**Only** use the device in combination with corrosion resistant components.


**Do not** overtighten the locking pin.

---

## 3 Construction

### Principal Parts

- Fabric (polyamide and lycra)
- Valve (silicone)
- Washer (nylon)
- Umbrella (nylon)
- Spacer\*\* (silicone)
- Retaining Screw\* (nylon)
- Casting Dummy (silicone)

\* For transportation only. **Do not** use for fitting.

\*\* For use with some locks. (See *Fitting the Spacer* section.)

---

## 4 Function

The device provides control, secure connection and cushioning between the limb and the socket. The perforations of the device allow moisture to escape during stance.

During stance, the valve opens to allow moisture and air to escape through the perforations in the distal cap. During swing, the valve closes and stops air from re-entering the device, thus creating a vacuum, which both improves proprioception and reduces pistoning.

---

## 5 Maintenance

Advise users to report the following to their practitioner:

- Tears in the fabric or the silicone
- Tissue damage on the residual limb
- Distal oedematous swellings
- Enlarged perforations
- Changes in either body weight or activity level
- Deterioration/changes to the residual limb
- Changes in the performance of the device, for example:
  - Moisture in the distal end of the device
  - Loss of vacuum

Note... Perspiration discolors some socket materials.

## 5.1 Cleaning the Device

Wash the inside of the device daily to avoid any build up of bacteria.

**Take care handling the device when it is inside out to avoid picking up dust, grit and other contamination which could cause irritation to the skin.**

**Dry thoroughly before use.**

**Do not tumble dry.**

### Washing by Hand

1. Invert the device so that its silicone side points externally.
2. Clean the silicone with a solution of water and unperfumed, pH balanced soap.
3. Invert the device so that its silicone side points internally. Clean the valve and carefully clean the distal cap.
4. Fill the device with warm water while both holding the proximal end closed and squeezing the distal end to flush the water through the perforations.
5. Rinse the device with clean water to remove all residues.
6. Either pat dry the device with a lint-free cloth, or leave it to air dry. Take care when drying and handling the device.

Note... Always dry the device with the silicone side of the device pointing internally. Otherwise, the fabric can stretch and become distorted.

### Machine Wash

Suitable for machine washing at 30 °C.

**Do not invert the device for machine washing.**

## 5.2 Cleaning the Valve

Make sure the distal end perforations are not clogged and that there is no contamination/debris is trapped under the valve.

1. During hand wash squeeze water through the distal perforations.
2. Carefully, slightly lift valve and wipe clean underneath with a clean piece of cloth.
3. Check correct operation/function of valve.

## 5.3 Cleaning the Residual Limb

Inspect residual limb before and after prosthetic limb use or at least daily.

**Any deterioration in residuum condition should be reported to the practitioner. If the residual limb has tissue damage, the use of the device should be stopped.**

1. Clean the skin daily with an unperfumed, pH balanced soap.
2. Rinse the skin with clean water to remove all residues.
3. Apply lotion to dry skin, as recommended by a practitioner.

**Make sure that any damaged skin or open wound is properly and suitably dressed to prevent direct contact with the device.**

# 6 Limitations on Use

## Intended Life

A local risk assessment should be carried out based upon activity and usage.

## Environment

Avoid exposing the device to corrosive elements such as acids, industrial detergents, bleach or chlorine. The use of creams or lotions with this device should be used with caution as these can cause the device to soften and overstretch or distort.

Keep away from sharp objects (such as jewelry, fingernails).

Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).


Suitable for shower use

# 7 Choosing the Correct Size

It is important that the correct size is specified to ensure proper fit and comfort in use.

1. Measure the circumference of the residual limb 4 cm up from its distal end.
2. If the circumference of the residual limb matches one of the sizes in the table, choose the next size down.
3. If the circumference of the residual limb is between one of the sizes in the table, choose the smaller size.

|                             |           |         |           |         |         |         |         |         |         |         |
|-----------------------------|-----------|---------|-----------|---------|---------|---------|---------|---------|---------|---------|
| Circumference Measured (cm) | 22.2-23.5 | 23.6-25 | 25.1-26.5 | 26.6-28 | 28.1-30 | 30.1-32 | 32.1-34 | 34.1-36 | 36.1-40 | 40.1-42 |
| Available Sizes             | 22        | 23.5    | 25        | 26.5    | 28      | 30      | 32      | 34      | 36      | 40      |

## 8 Casting/Scanning the Socket with the Casting Dummy

When casting the socket, use the casting dummy to make sure that the distal end of the socket has space for the valve to open. Otherwise, moisture will collect inside the device.

### Before you begin

Allow the user to wear the device for 10 minutes.

(See *Donning the Device* section.)

938413

### If a scanning method is used


Fit/trim the casting dummy to the size of the distal end of the liner, see 938413.

To reduce the amount of distal pressure within the final socket, the model should be lengthened by 10-12 mm (depending upon redundant tissue etc.).

### If a casting method is used

Fit/trim the casting dummy to the size of the distal end of the liner, see 938413.

We recommend using a casting method that accentuates weight bearing areas such as the medial condylar flare in conjunction with applying tension to the distal pin; otherwise excessive distal contact may result in the final socket, giving rise to oedematous swellings corresponding to the distal perforations as shown below.

Additionally, when fabricating a diagnostic socket, we recommend a 5 mm pelite spacer is used to further extend/lengthen the socket. When you are satisfied with the fit of the socket (by adding any necessary socks) the spacer can, if necessary, be removed if excess distal contact gives rise to distal oedematous swellings.

 **The softness of the silicone can mask the amount of distal contact that is occurring, this will only be evident when oedematous swellings are observed.**

 **The use of a diagnostic socket is strongly recommended to assess the amount of distal pressure/contact that is occurring.**

## 9 Trimming the Device

The device may be trimmed to suit as required but never bellow the socket trim line.

We recommend using a rotating disc cutter to achieve a bevelled edge.

Leave enough material on the device so that its proximal edge can cover the edge of the socket.

Avoid trimming through the perforations.

Trim the proximal edge of the device to suit the shape of the user's socket.

If necessary, trim the proximal edge of the device in a wave pattern to reduce shear stress.

**Where possible try to avoid abrupt changes of contour and sharp edges which could cause tears in the silicone and fabric.**

## 10 Donning the Device

### Before you begin

Fit the valve. (See *Fitting the Valve* section.)

**Care must be taken during donning/doffing not to damage the device with fingernails, sharp jewelry or the locking pin.**




**Do NOT pull or stretch the device.**

1. Invert the device so that its silicone side points externally.
2. Align the locking pin with the long axis of the residual limb.
3. Roll the device on the residual limb while releasing any trapped air.

**Note:** The unperforated fabric should be below the patella tendon.

After donning the device, allow the user to wear it for 10 minutes. If the user feels numbness, tingling or any unusual sensation within 10 minutes, do the following:





4. Do off the device, and wait until normal sensation returns.
5. Don the device again.
6. If the user feels numbness, tingling or any unusual sensation again, stop using the device.

## 11 Fitting the Valve


**⚠ Do not overtighten the locking pin.**

We recommend using a locking pin that has a shoulder of between 13.5 mm and 19 mm in diameter.

## 12 Fitting the Spacer

Fit the spacer between the lock and the valve only if the lock prevents the valve from opening.

## 13 Fitting Advice

### 13.1 Moisture in the Distal End of the Device

If moisture collects inside the distal end of the device, refer to the following table:

| Cause                                                                          | Solution                                                                                                                            |
|--------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| The distal perforations are clogged.                                           | Clean the device. (See <i>Maintenance</i> section.)                                                                                 |
| The socket does not have enough space in its distal end for the valve to open. | Put the spacer between the pin lock and the valve. (See <i>Fitting the Spacer</i> section.)<br>If unsuccessful cast another socket. |
| The internal shape of the lock-body is stopping the valve from opening.        | Put the spacer between the pin lock and the valve. (See <i>Fitting the Spacer</i> section.)                                         |

### 13.2 Loss of Vacuum

If the device loses vacuum, refer to the following table:

| Cause                                                                   | Solution                                                              |
|-------------------------------------------------------------------------|-----------------------------------------------------------------------|
| The valve is damaged.                                                   | Stop using the device, and contact a Blatchford sales representative. |
| Debris beneath the valve preventing correct seal.                       | Clean the valve. (See <i>Cleaning the Valve</i> section.)             |
| The locking pin is not aligned with the long axis of the residual limb. | Re-align the locking pin with the long axis of the residual limb.     |

### 13.3 Separation of the Valve from the Distal Cap

If the valve falls off the distal cap, refer to the following table:

| Cause                              | Solution                                                                                            |
|------------------------------------|-----------------------------------------------------------------------------------------------------|
| The valve is not fitted correctly. | Make sure that the washer is between the pin and the valve. (See <i>Fitting the Valve</i> section.) |

Note... Do not over-tighten the locking pin/lock rod.

### 13.4 Oedematous Swellings Corresponding to Distal Perforations


| Cause                     | Solution                                                                 |
|---------------------------|--------------------------------------------------------------------------|
| Excessive distal contact. | Reduce end contact by addition of extra socks or lengthen/remake socket. |

# 14 Technical Data

|                                                               |                                              |
|---------------------------------------------------------------|----------------------------------------------|
| Principal Materials                                           | polyamide, lycra, silicone, nylon            |
| Shore Hardness                                                | 40 shore 00                                  |
| Component Weight (size 28)                                    | 695 g (1 lb 8 oz)                            |
| Activity Level                                                | 3-4                                          |
| Size Range<br>(See <i>Choosing the Correct Size</i> section.) | 22-40 cm                                     |
| Length<br>(See diagram below)                                 | 435 mm                                       |
| Internal Length<br>(See diagram below)                        | 420 mm                                       |
| Matrix Length                                                 | Approx. 10 cm                                |
| Operating and Storage Temperature Range                       | -15 °C to 50 °C (5 °F to 122 °F)             |
| Distal End Attachment                                         | M10<br>Requires Locking Pin With a Shoulder* |
| Locking Pin Shoulder Diameter                                 | 13.5-19 mm                                   |

\*Locking Pin Not Supplied

## Dimensions

# 15 Ordering Information

## Order Example

|                    |          |             |          |          |                  |
|--------------------|----------|-------------|----------|----------|------------------|
| <b>SB</b>          | <b>A</b> | <b>TT</b>   | <b>L</b> | <b>P</b> | <b>25</b>        |
| Silcare<br>Breathe | Active   | Transtibial | Locking  | Parallel | Size*<br>(22-40) |

Available from size 22 to size 40\*:  
SBATTLP22 to SBATTLP40

\*SBATTLP23 is for size 23.5. SBATTLP26 is for size 26.5.

| Part                  | Part Number |
|-----------------------|-------------|
| Small Valve Kit       | 559015      |
| Medium Valve Kit      | 559016      |
| Large Valve Kit       | 559017      |
| Extra Large Valve Kit | 559018      |
| Casting Dummy         | 559019      |

## Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

## CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk)

Medical Device

Single Patient – multiple use

## Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

## Warranty

This device warranted for 6 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses, and exemptions.

See the Blatchford website for the current full warranty statement.

## Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## Environmental Aspects

This product is fabricated from silicone rubber and fabric that cannot be easily recycled: please dispose of it responsibly as general waste, as per local waste handling regulations.

## Retaining the Packaging Label

You are advised to keep the packaging label as a record of the device supplied.

## Trademark Acknowledgments

Silcare Breathe and Blatchford are registered trademarks of Blatchford Products Limited.

## Manufacturer's Registered Address

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© Blatchford Products Limited 2024. All rights reserved.
**Endolite India Ltd.**
Blatchford Europe GmbH