# Orion<sup>3</sup>

## User Guide

ORION3

EN	User Guide	2
FR	Guide de l'utilisateur	17
DE	Bedienungsanleitung	32
IT	Guida per il paziente	47
ES	Guía del usuario	62
NL	Gebruikershandleiding	77
PL	Podręcznik użytkownika	92
PT	Guia do utilizador	107
CS	Návod k obsluze	122

Contents .....	2
1 Description and Intended Purpose .....	3
2 Safety Information .....	4
3 Construction .....	6
4 Key Points.....	7
5 Transport and Storage.....	8
6 Battery Charging .....	9
6.1 Magnetic Charger Connection.....	11
6.2 Battery Status .....	11
7 User Modes.....	12
7.1 Flexion Lock Mode.....	12
7.2 Cycling Mode.....	13
7.3 Training Mode.....	13

## Battery Safety


The device contains Li-ion type rechargeable batteries. Please read and comply with the safety information see Section 2 and Section 6. Failure to do so may cause the battery to become hot, explode, or ignite and cause serious injury.

### IMPORTANT: Battery Charging

To charge the device ensure that the On/Off Switch is in the ON [I] position. When the device is charging the Battery Charger LED will be flashing green.

The battery should only be charged indoors in a dry environment within the temperature range 10 °C to 40 °C (50 °F to 104 °F).

After exposure to water/moisture the Orion3 must be fully drained of any water/moisture ingress and external surfaces wiped dry prior to charging, in particular the Magnetic Charging Terminal.

## Low Temperature Operation

IMPORTANT: If the device is subjected to temperatures lower than -10°C (14°F) it must be returned to be inspected for possible damage to the battery pack.

Following extended periods of inactivity at temperatures below 0°C (32°F) the knee may swing less freely. In such an eventuality, it is recommended to exercise caution while descending stairs and avoid descending stairs leg-over-leg using the knee yield until the user feels the swing return to normal.

## Magnetic Connector Cleaning

Regular cleaning of the magnetic charging surfaces is required to ensure no liquids, contamination or debris prevents correct connection. (See Section 4 for cleaning advice.)

## Package Contents

1. Knee prosthesis
  2. Detachable Cosmetic Cover
  3. User Guide
  4. Cable Bag and accessories
    - i) Battery Charger
    - ii) Battery Charger AC power socket outlet adapter(s)
- 

## 1 Description and Intended Purpose

This document describes the IP (Ingress Protection)-rated version of the Orion3 knee.

These instructions are for the user. Please keep these instructions.

The term *device* is used throughout this manual to refer to Orion3.

Please read all these instructions for use, paying particular attention to the sections regarding safety and maintenance.

### Application

The device is to be used exclusively as part of a lower limb prosthesis.

Intended for single user.

The device is a microprocessor-controlled knee prosthesis that provides swing and stance control and adapts to variable cadence. It has a programmable stance yield setting for controlled descent on ramps and stairs and during active gait with a cadence range from slow walking to fast walking. The knee also provides an additional standing support mode for enhanced security when standing.

The set-up procedure is through an easy to use App interface via a **Bluetooth**® connection to the on-board module.

### Features

- Microprocessor-controlled stance and swing
- Standing support mode that controls knee flexion for standing
- App based set-up via a built-in **Bluetooth**® link
- Up to 3 days of normal use from the internal rechargeable Li-ion battery pack
- 3 User modes
- Activity monitoring of different activities
- Water and dust resilient to IP55 (see Section 4)

In addition the device offers enhanced stability:

- Active control of heel rise during walking
- Dynamic stair descent with easy initiation of knee flexion with customizable yield rate (resistance level increases with knee flexion angle)
- Supportive resistance to flexion as soon as the knee stops flexing
- Increased supportive resistance to flexion to prevent 'stumbles' associated with instability at heel strike
- Supportive yield for sitting down

## 2 Safety Information

 This warning symbol highlights important safety information which must be followed carefully.

 Any changes in the performance or function of the knee e.g. restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

 Always use a hand rail when descending stairs and at any other time if available.

 Any excessive changes in heel height after programming will adversely affect knee function and should be immediately reported to your service provider for reprogramming and calibration.

 After continuous use the knee casing may become hot to the touch.

 Avoid strong magnetic fields, sources of electrical interference, atmospheres containing liquids and/or powders.

 Do not place near any heat source. Do not leave in direct sunshine or inside a car in hot weather.

 The device is not intended for use when immersed in water or as a shower prosthesis.

 Do not walk in water with this device. The water resistance may prevent the limb from extending correctly.

 The device should not be used while showering or bathing due to the potential for seals to be compromised by cleaning substances or detergents and also due to the risk of slipping.

 If the device is exposed to corrosive substances such as salt water, chlorinated water or other solutions (including shampoos, shower gels, soaps and body fluids) it should be rinsed as soon as possible using fresh water, taking care not to exceed the water resistance limitations defined by the IP55 rating (see <sup>1</sup> below).

 When rinsing the device only use low pressure showers or water jets that do not exceed a flow rate of 12.5 litres/minute and pressure of 30 kPa (4.4 psi.)

 If the device has been exposed to significant water ingress there is potential for a slip hazard, depending on the floor surface, while the device is draining.

 Ensure that no solid or foreign objects are inserted into the device, including the drainage holes which should only be cleaned or unblocked as recommended in Section 4 *Key Points*.

 Before connecting the battery charger to the device, ensure that the magnetic mating surfaces of the charging connector and charging terminal are clean, dry and free of any debris or conductive materials that could prevent correct charging operation. (See further warnings in Section 6.1 and cleaning advice in Section 4.)

<sup>1</sup> Water resistance tested from all directions using a water jet from a 6.3 mm (1/4") dia. nozzle at a distance of 3 m (9.8 ft) with fresh water @12.5 litres/minute and 30 kPa (4.4 psi) for a maximum of 3 minutes.


-  Never connect the charger to the device while still wearing it.
-  The charger is not IP rated and is only designed for use in a dry and dust free environment. Never charge the device in or near a wet or damp environment.
-  Be aware of finger trap hazard at all times.
-  The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling is acceptable.
-  The device is not suitable for swimming or paddling/walking in water.
-  Contact your service provider for advice if the device is to be used with a pacemaker or any other electronic medical device.
-  Assembly, programming, maintenance and repair of the device must only be carried out by a suitably qualified clinician that has attended an approved training course.
-  The user must not adjust or tamper with the setup of the device.
-  The user should be advised to contact their practitioner if their condition changes.
-  The device detects standing and walking at various speeds. Non-identified movements may cause unexpected behavior.
-  Only charge the batteries in the device using the charger supplied with the product. Do not use this charger for charging any other device.
-  Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.
-  The charging connector and terminal are sources of strong magnetic fields. Always keep a safe distance (at least 10 cm) between these magnets and all objects that can be damaged or affected by magnetism (such as watches, pacemakers, credit cards, or any magnetically stored media).

# 3 Construction

## Principal Parts

- Shin (CF Composite)
- Chassis (Aluminum Alloy)
- Pyramid (Stainless Steel)
- Hydraulic/Pneumatic Cylinder (Various materials including Aluminum Alloy, Stainless Steel, hydraulic oil, motors and PCBs)
- Knee Pad/Battery Case and Drainage Cover (Polyurethane)
- Batteries (Lithium-ion)

## Component Identification

---

## 4 Key Points

**The device must not be used with large variations in heel height.**

**Do not attempt stairs descent leg-over-leg on this device unless you have been trained for stairs descent.**

Any changes in body weight and/or activity level of the user must also be reported to the practitioner.

High flexion resistance settings may induce some noise during stairs descent, this will not affect function.

Any changes in performance of this device must be reported to the practitioner.

Changes in performance may include:

- Increase or decrease in knee stiffness
- Instability
- Reduced knee support (free movement)
- Any unusual noise
- Any unrecognized beeps or warning lights

The practitioner must also be informed of the following:

- Any changes in body weight and/or activity level
- Significant wear or excessive discoloration from long term exposure to UV

### Cleaning

Use a damp cloth and mild soap to clean outside surfaces, DO NOT use aggressive cleansers. When rinsing the device only use low pressure showers or water jets that do not exceed a flow rate of 12.5 litres per minute and pressure of 30 kPa (4.4 psi). Thoroughly drain and dry before use.

To clean the magnetic mating surfaces of the charging connector and charging terminal use a clean and dry cotton swab or a soft bristle brush. Make sure that no cotton fibers or bristles are left on the surfaces.

**Never clean or expose the device to moisture while charging.**

**Never apply any cleaning liquid, detergents or products to the magnetic mating surfaces.**

## Environment

The device is water resistant to IP55<sup>1</sup> and can withstand temporary exposure to a low-pressure jet of fresh water with no harmful effects, however the device is not to be used for showering, bathing, paddling/walking in water, swimming or prolonged submersion.

Avoid exposing the device to corrosive elements such as acids, detergents (including washing products such as shampoos, shower gels and soaps) and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.

The charging surfaces of the connectors should be kept away from other magnetic or conductive materials.

Exclusively for use between -10 °C and 50 °C (14 °F to 122 °F).

It is recommended that only Blatchford products be used in conjunction with the device.

<sup>1</sup> Water resistance tested from all directions using a water jet from a 6.3 mm (1/4") dia. nozzle at a distance of 3 m (9.8 ft) with fresh water @12.5 litres/minute and 30 kPa (4.4 psi) for a maximum of 3 minutes.

---

## 5 Transport and Storage

1. Ensure that the device is in a straight (fully extended) position when not in use.
2. If the limb is to be stored upright then it must be secured to prevent it falling over or flexing.
3. If the device is shipped either on its own or built into a limb it should be maintained in an extended position and packaged in a suitable box to allow this.
4. If a unit is to be returned, contact Blatchford prior to sending if there is any sign of physical damage to the battery pack.
5. For transport purposes the unit must be switched off, and packed securely so that a) it cannot be inadvertently switched on and become operative while in transit; b) the unit is thoroughly protected from damage in transit.

## 6 Battery Charging


The device contains Li-ion type rechargeable batteries. Please read and comply with the safety information. Failure to do so may cause the battery to become hot, explode, or ignite and cause serious injury.

Please familiarize yourself with the charging information in this section.

The batteries are not field replaceable and will be replaced as necessary when serviced.

If the device is subjected to temperatures lower than  $-10^{\circ}\text{C}$  ( $14^{\circ}\text{F}$ ) it must be returned to be inspected for possible damage to the battery pack.

1. The battery should only be charged in temperatures between  $10^{\circ}\text{C}$  to  $40^{\circ}\text{C}$  ( $50^{\circ}\text{F}$  to  $104^{\circ}\text{F}$ ).

2. Only charge the batteries in the device using the charger supplied with the product. Do not use this charger for charging any other device.

3. The battery charger can be used with an electrical outlet that is rated between 100V and 240V AC, 50/60 Hz.

4. Never connect the charger while still wearing the limb.

5. During charging no power is delivered to the device control system.

6. The knee module may become warm during charging but should never be too hot to touch. There should be no unusual smell, if there is immediately switch **OFF[O]** the battery, unplug the charger from the electrical supply and contact your service provider. .

7. Do not refit the knee until the charger has been disconnected and the device has completed its start-up sequence (see **Step by Step Guide** later in this section).

8. Always recharge the device before use if it has not been used for an extended period of time.


9. The battery should only be charged indoors in a dry and dust free environment.

Fully charge the battery before first use.

We recommend that the device should be fully charged every night. It should take up to 8 hours to charge fully. We also recommend battery charging whenever convenient.

### Step by Step Guide


1. Lay the knee on a flat surface.
2. Ensure the mating surfaces of the magnetic charging connector and terminal are clean (see Section 4).
3. Magnetically join the charging connector to the charging terminal at the posterior of the knee and check alignment is correct (see Section 6.1).
4. Switch the Battery On/Off switch to **ON [I]**.
5. Check the Battery Charger LED is flashing green to confirm the battery is charging.
6. Refer to table below to ensure the knee is charging.
7. Disconnect from the charger when charged sufficiently.
8. Wait for 30 seconds while observing the device LED changing color as indicated:


9. After a successful reset sequence 2 short beeps will be heard. (if this alert is not disabled; see\* below)
10. If the reset sequence is unsuccessful or the LED remains Red switch the battery **OFF [O]** then **ON [I]** and repeat from step 8.

Battery Charger LED	Charging State
Yellow	Standby
Green > fast flashing	Fast charging
Green > slow flashing	Close to full charge
Green	Maintain (trickle charge)
Yellow > fast flashing	Error (Ensure Battery Switch is <b>ON [I]</b> )

Note... If the battery charger LED is yellow and the knee is making a clicking noise check that the Battery switch is in the **ON [I]** position.

Note... Do not move or don the limb while reset is processing otherwise this will prolong the reset time.

\*The user can enable/disable the audible reset confirmation by pressing and holding the [+ ] and [- ] buttons simultaneously for 2 seconds. The change will be confirmed by 2 short beeps.

## 6.1 Magnetic Charger Connection

The charging connector can be connected to the device through the charging terminal located at the posterior of the knee.

The connector and the terminal magnetically lock together and correctly self-align on connection (Figure 1).

Note... The connector and the terminal are not polarized so orientation is not critical (Figure 2).


 Please check for any misalignment between mating surfaces as this may result in no charging. This can happen especially when the charging setup puts strain on the charging cable (Figure 3).

 Before charging ensure charging surfaces are clean, dry and free of any debris.

 Do not allow any ferrous material to attach itself to the charging surfaces.


 Keep the magnetic surfaces away from materials that may be adversely affected by their proximity to strong magnetic fields (e.g. watches, pacemakers, bank cards or magnetically stored media).

 Prevent build up of debris on the charging surfaces as this may prevent charging or damage the charger. Disconnect the charger from the mains socket and remove any debris before reconnecting the charger to the device.

## 6.2 Battery Status

To check the battery condition press and hold the [-] button on the knee, one short beep will be heard and the LED will indicate the battery condition.

Battery Charge Indicator			
LED		Beeps	Battery Condition
Green	3 flashes	--	fully charged
Green	2 flashes	--	good
Green	1 flash	--	OK
Yellow	slow flash	8 x Short beeps	low charge
Red	slow flash	5 x Long beeps	no charge

When fully charged the device battery is intended to provide 3 days of normal use (conditional on the type of use). We recommend that the knee be fully charged every night. Switching off the device when not in use will also prolong the period of use

Note... The device will take 30 seconds to reset after it is switched back **ON [I]**, refer to Section 6 'Battery Charging (continued)' for the reset sequence.

---

## 7 User Modes

The device has 3 user modes. Please note that these modes are only available if they have been enabled by your practitioner. Take extra care while using these user modes.

 **WARNING: Care must be taken to only use those modes as instructed by your practitioner.**

### 1. Flexion Lock Mode

Sets the limit of flexion to a predetermined angle set by the practitioner.

The maximum lock angle is 45 degrees from fully extended.

### 2. Cycling Mode (Free swing)

Releases stance resistance, the device will then swing freely.

 **WARNING: In Cycling Mode there is no stance control active, this mode must be de-activated before resuming walking to avoid falling and possible injury.**

### 3. Training Mode

 **WARNING: In Training Mode, the device has a predefined fixed pneumatic resistance setting which does not adapt to changes in walking speed. Also, Standing Mode is not active in Training mode.**

## 7.1 Flexion Lock Mode

To activate/de-activate Flexion Lock mode use the [-] button on the device, as follows:

### To activate Flexion Lock Mode:

- Fully extend the knee.
- Press the [-] button 3 times (each time pressing until the device responds with a single short beep).
- 3 medium beeps then sound to confirm that Flexion Lock is ON.

### To de-activate Flexion Lock Mode:

- Press and hold the [-] button.
- 2 long beeps then sound to confirm that Flexion Lock is OFF and normal operation has resumed.

 **WARNING: After exiting Flexion Lock Mode, always test/confirm that the device is operating normally again before walking.**

## 7.2 Cycling Mode

 **WARNING: In cycling (free swing) mode there is NO stance resistance.**

To activate/de-activate Cycling Mode use the [-] button on the device, as follows:

### To activate Cycling Mode:

- Sit on bicycle before activating Cycling Mode as once activated there is no support for walking.
- Press the [-] button 4 times (each time pressing until the device responds with a single short beep).
- 4 medium beeps then sound to confirm that Cycling Mode is ON.

### To de-activate Cycling Mode:

- Press and hold the [-] button.
- 2 long beeps then sound to confirm that Cycling Mode is OFF and normal operation has resumed.

 **WARNING: After exiting Cycling Mode, always test/confirm that the limb is operating normally again before walking.**

## 7.3 Training Mode

 **WARNING: In Training Mode, the device has a predefined fixed pneumatic resistance setting which does not adapt to changes in walking speed. Also, Standing Mode is not active in Training mode.**

To activate/de-activate Training Mode use the [-] button on the device, as follows:

### To switch ON Training Mode:

- Press the [-] button 5 times (each time pressing until the device responds with a single short beep).
- 5 medium beeps then sound to confirm that Training Mode is ON.

### To switch OFF Training Mode:

- Press and hold the [-] button.
- 2 long beeps then sound to confirm that Training Mode is OFF and normal operation has resumed.

 **WARNING: After exiting Training Mode, always test/confirm that the limb is operating normally again before walking.**

---

# Manufacturer's Declarations and Legal Information

---

## FCC Compliance Statement

This device complies with part 15 of the FCC Rules.

Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

Note... Contains pre-approved **Bluetooth®** module: FCC ID: PVH0946 IC: 5325A-0946

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation. This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party: Blatchford Inc.

1031 Byers Road, Miamisburg, OH 45342  
800-548-3534 (toll free) | 937-291-3636  
customerservice@blatchfordus.com

## RSS Compliance Statement

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

1. this device may not cause interference, and
2. this device must accept any interference, including interference that may cause undesired operation of this device.

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Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website Blatchford Inc. 1031 Byers Road, Miamisburg, OH 45342.

800-548-3534 (toll free) | 937-291-3636 customerservice@blatchfordus.com

## Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

## CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchfordmobility.com](http://www.blatchfordmobility.com)

Medical Device

Single Patient – multiple use

## Warranty

This device is warranted for 36 months.

Extended warranties of an additional 12, 24, or 36 months are available.

A service between 30 and 36 months is required for all extended warranties.

For the maximum 36-month extended warranty, a second service must be carried out between 54 and 60 months.

We reserve the right to void extended warranties if the device is not serviced at the required intervals.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

The above may vary by market; consult your local representative for details.

See the Blatchford website for the current full warranty statement.

## Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## Environmental Aspects

This symbol indicates that the product contains electrical/electronic components and/or batteries that should not be disposed of in general waste or be incinerated at the end of the product's life.

At the end of the product's life, all electrical/electronic components and/or batteries should be recycled or disposed of in accordance with the current regulations for the handling of WEEE (Waste Electrical and Electronic Equipment), or equivalent local regulations. The remainder of the product should also be recycled where possible in accordance with local waste recycling regulations.

To help prevent potential harm to the environment or to human health from uncontrolled waste disposal, Blatchford offers a take-back service. Please contact Customer Services for details.

## Trademark Acknowledgements

Orion3 and Blatchford are registered trademarks of Blatchford Products Limited.

The **Bluetooth**® word mark and logos are registered trademarks owned by Bluetooth SIG Inc., and any use of such marks by Blatchford is under License. Other trademarks and trade names are those of their respective owners.

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© Blatchford Products Limited 2025. All rights reserved.
Blatchford Europe GmbH
Blatchford Products Ltd.