Linx

Clinician's Manual

LINX22L1S–LINX30R8S
LINX22L1SD–LINX30R8SD

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Contents	2
Package Contents:	3
1 Description and Intended Purpose	4
1.1 Spring Set Selection	6
2 Safety Information	7
3 Construction	8
4 Function.....	9
5 Maintenance	11
6 Limitations on Use	12
7 Battery Charging	13
7.1 Magnetic Charger Connection	15
7.2 Battery Status	15
8 Alignment	16
8.1 Assemble the Complete Limb	16
8.2 Adjust the Pylon Length.....	16
8.3 Cut the Pylon to Length.....	17
8.4 Re-assemble the Limb.....	17
8.5 Bench Alignment	18
9 Device Set-up	19
9.1 Using the Linx Programming App	19
9.2 Button Programming.....	23
10 User Modes.....	28
10.1 Flexion Lock Mode.....	28
10.2 Cycling Mode.....	29
10.3 Training Mode.....	29
11 Technical Data	30
Appendix 1 Installing the Linx Programming App (iOS).....	32
Appendix 2 Installing the Linx Programming App (Android)	33

Warnings

The term *device* is used throughout this manual to refer to Linx.


Battery Safety

The device contains Li-ion type rechargeable batteries. Please read and comply with the safety information (see Section 2 and Section 7). Failure to do so may cause the battery to become hot, explode, or ignite and cause serious injury.

IMPORTANT: Battery Charging

To charge the device ensure that the On/Off Switch is in the **ON [I]** position. When the device is charging the Battery Charger LED will be flashing green.

The battery should only be charged within the temperature range 10 °C to 40 °C (50 °F to 104 °F).

Battery On/Off Switch

Low Temperature Operation

IMPORTANT: If the limb is subjected to temperatures lower than -10 °C (14 °F) it must be returned to be inspected for possible damage to the battery pack.

Following extended periods of inactivity at temperatures below 0 °C (32 °F) the limb may swing less freely. In such an eventuality, it is recommended to exercise caution while descending stairs and avoid descending stairs leg-over-leg using the knee yield until the user feels the swing return to normal.

Package Contents:

1. Knee
2. Ankle, pylon, spring set and foot shell
3. Detachable Cosmetic Cover
4. USB stick with Linx Android App
5. Clinician's Manual
6. User Guide
7. Cable bag and accessories
 - i) Battery charger
 - ii) Battery charger AC power socket outlet adapter(s)
 - iii) Pylon cutting jig
 - iv) Shin extension cable
 - v) Allen tool hex 2 mm

1 Description and Intended Purpose

This document describes the V2 version of the Linx limb.

These instructions are for the practitioner.

Ensure that the user has understood all instructions for use, drawing particular attention to the sections regarding safety and maintenance.

Application

The device is to be used exclusively as part of a lower limb prosthesis.

Intended for a single user.

The *device* is a fully integrated, microprocessor-controlled knee and foot system. It provides simultaneous control over both foot and knee, co-ordinating their responses for a variety of terrains and situations in one complete biomimetic limb. Individual adjustments can be made for slopes and stairs, and there is a standing mode which locks the device to assist when standing.

The device is configured for a specific user via the easy-to-use Linx programming app (available for IOS and Android), connecting via a Bluetooth® wireless technology connection to an on-board module. The app provides comprehensive programming and fine-tuning functions to produce a highly-refined device set-up that accurately responds to the user's needs and operating environment.

If the app is unavailable, the device can also be configured using button programming, via buttons at the rear of the knee. However, the app offers a greater degree of programming control for a more refined set-up than is possible via button programming.

Features

- All the functions of microprocessor control in both knee and hydraulic ankle
- Centralised control of knee and foot to give the user a better coordinated response for standing, level walk, ramp ascent and descent.
- App-based set-up via a built-in **Bluetooth®** link
- Proven to reduce stress at the socket interface
- An easy access single charging and point
- Longer battery life, 3 days average use from the internal rechargeable Li-ion battery pack
- 3 User modes
- Capability to program whole limb via buttons on the rear of the knee
- Activity monitoring of different activities

The enhanced stability system offers:

- **Standing Mode**, control of knee flexion and ankle for standing
- Controlled ramp descent with *Braking* effect of knee and ankle
- Enhanced ramp ascent with optimized effect of knee and ankle
- Active control of heel rise during walking

- Dynamic stair descent with an easy start of knee flexion coupled with increasing support to gently lower the user during flexion
- Supportive resistance to flexion as soon as the knee stops flexing
- Increased supportive resistance to flexion to prevent 'stumbles' associated with instability at heel strike.

Clinical Benefits

- Increased ground clearance and balance reduces risk of trips and falls
- Improved safety and ground compliance over uneven surfaces and slopes
- Improved kinetic gait symmetry
- Reduced energy expenditure
- Increased walking speed
- Improved adaptation to different walking speeds

Activity Level

This device is recommended for users that have the potential to achieve Activity Level 3 who may benefit from enhanced stability and an increase in confidence on uneven surfaces.

There are exceptions and in our recommendation we want to allow for unique, individual circumstances. There may also be a number of users of Activity Levels 2 and 4* who would benefit from the enhanced stability offered by the device, but this decision should be made with sound and thorough justification.

Users at Activity Levels 2 and 4 will require softer or stiffer springs as appropriate for the individual user rather than as shown in the spring selection guide.

*(Maximum user weight 100 kg and always use one higher spring rate category than shown in the Spring Set Selection table, Section 1.1).

For bilateral transfemoral please contact Blatchford technical support for training and assessment when descending stairs as this activity carries higher risk of injury.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Contraindications

This device may not be suitable for Activity Level 1 individuals or for competitive sports events, as these types of users will be better served by a specially designed prosthesis optimized for their needs.

The device must not be used with large variations in heel height.

1.1 Spring Set Selection

Activity Level 3

44-52	53-59	60-68	69-77	78-88	89-100	101-116	117-125	kg	User Weight
100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-275	lbs	
1	2	3	4	5	6	7	8		Foot Spring Set

Note... If in doubt choosing between two categories, choose the higher rate spring set.

2 Safety Information

The caution symbol highlights safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

Any excessive changes in heel height after programming will adversely affect limb function and should be immediately reported to your service provider for reprogramming and calibration.

After continuous use the ankle and knee casing may become hot to the touch.

Avoid strong magnetic fields, sources of electrical interference, atmospheres containing liquids and/or powders.

Do not place near any heat source. Do not leave in direct sunshine or inside a car in hot weather.

The device is not intended for use when immersed in water or as a shower prosthesis. If the limb comes into contact with water wipe dry immediately.

Ensure no water enters the charger plug socket and the cover is on at all times.

Before connecting the battery charger to the device, ensure that the magnetic mating surfaces of the charging connector and charging terminal are clean, dry and free of any debris or conductive materials that could prevent correct charging operation. (See further warnings in Section 7.1 and cleaning advice in Section 5).

Never connect the charger to the device while still wearing the limb.

The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling is acceptable.

Be aware of finger trap hazard at all times

Contact your service provider for advice if the device is to be used with a pacemaker or any other electronic medical device.

Assembly, programming, maintenance and repair of the device must only be carried out by a suitably qualified clinician that has attended an approved training course.

The user must not adjust or tamper with the setup of the device.

The user should be advised to contact their clinician if their condition changes.

The device detects standing and walking on inclines, declines and walking at various speeds. Non-identified movements may cause unexpected behavior.

Only charge the batteries in the device using the charger supplied with the product. Do not use this charger for charging anything else.

Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

Make sure that no liquids enter the device during donning or normal use.

The charging connector and terminal are sources of strong magnetic fields. Always keep a safe distance (at least 10 cm) between these magnets and all objects that can be damaged or affected by magnetism (such as watches, pacemakers, credit cards, or any magnetically stored media).


To minimise the risk of slipping and tripping, appropriate footwear that fits securely onto the footshell must be used at all times.

3 Construction

Principal Parts

- Knee Module including pyramid (aluminium/St. Stl./titanium, carbon)
- Ankle Module (aluminium/St. Stl.)
- Pylon (aluminium)
- Heel & Toe Springs (e-Carbon)
- Spring Attachment Screws (titanium/St.Stl.)
- Glide Sock (UHM PE)
- Foot Shell, cables and flexi-boards (polymers and PU)
- Batteries (lithium ion)

4 Function

The device is a complete external transfemoral limb prosthesis, which utilizes sensors and actuators controlled via microprocessors in an integrated hydro-pneumatic knee and hydraulic foot system.

The Linx app-based interface programs the device to the individual user's gait characteristics during the **Calibration** stages of the set-up sequence. This set-up can then be fine-tuned to meet specific requirements. Alternatively, programming can be achieved via the buttons at the rear of the knee, but with less fine tuning options than are available with the app.

The limb provides control of the stance and swing phase of the user's gait while the stance control is optimized to provide secure support during ambulation on all terrains.

The enhanced stability of the system offers a **Standing Mode** with control of knee flexion and ankle for standing.

There is increased ground contact and better resistance to knee flexion during ramp descent. There is an optimized assistance for ramp ascent provided by the limb.

It has a dynamic action during stair descent with an easy start of knee flexion coupled with increasing support to gently lower the user during flexion. It also provides supportive resistance to flexion as soon as the knee stops flexing and increased supportive resistance to flexion to prevent 'stumbles' associated with instability at heel strike. Heel rise during walking is actively controlled via hydraulic damping.

The device is powered by a rechargeable lithium-ion battery pack, charged via a magnetically secured and self-aligning charger connection on the posterior of the knee.

Limb Operating Modes

Activity	Device Operating Mode/function	Knee module	Foot module
Walking (level ground)	Normal	Alternates stance (high yield resistance) to free swing (low resistance) and pneumatic swing control adjusts for cadence variation. Excessive heel rise is also actively damped hydraulically preventing excessive heel rise associated with sudden cadence variation.	Slightly dorsiflexed during swing Basic resistance settings in plantarflexion and dorsiflexion
Walking up ramp (moderate/ steep incline)	Ramp Ascent Enables easier walking up incline.	Progressively changes hydraulic damping and pneumatic assist to assist up the ramp; adjusts for cadence variation.	Progressively changes* stiffer plantarflexion and lower dorsiflexion resistance to assist up the ramp.
Walking down ramp	Ramp Descent Helps reduce forward momentum and provides greater stability and security.	Moves from a high yield resistance, through an intermediate resistance then to free swing (low resistance) and adjusts for cadence variation.	Progressively changes* decreasing plantarflexion and increasing dorsiflexion resistance to provide a braking action down the ramp.
Standing	Standing Mode Increases stability when standing still.	Very high resistance in knee and ankle (activates after 1-1.5 seconds).	
Any	Low battery power	High resistance (No dynamic control)	Reverts to basic plantarflexion & dorsiflexion resistance settings. (No dynamic control)

*Ramp detection depends on gait, starting, stopping, accelerating and decelerating and any of these can delay response.

5 Maintenance

A service is required every 36 months.

Maintenance must be carried out by competent personnel.

It is recommended that the following maintenance is carried out annually:

- Remove the foot shell and glide sock, check for damage or wear and replace if necessary.
- Check all screws for tightness (see Section 3), clean and reassemble as necessary.
- Visually check the heel and toe springs for signs of delamination or wear and replace if necessary. Some surface damage may occur after a period of use; this does not affect the function or strength of the foot.

The user must be handed the User Guide supplied, and be advised of the following:

Before connecting the charger to the device, always check that the magnetic mating surfaces of the charging connector and charging terminal are clean, dry and free of debris that could prevent correct charging operation; if not, clean, dry and remove any debris as necessary before connection (see Cleaning below).

Any changes in performance of this device must be reported to the clinician.

Changes in performance may include:

- Increase or decrease in knee and/or ankle stiffness
- Instability
- Knee starting to extend too aggressively
- Reduced ankle/knee support (free movement)
- Any unusual noise
- Any unrecognized beeps or warning lights.

The user should be advised that a regular visual check of the foot is recommended, signs of wear that may affect function should be reported to their service provider (e.g. significant wear or excessive discoloration from long term exposure to UV).

Cleaning

Use a damp cloth and mild soap to clean outside surfaces, DO NOT use aggressive cleansers. Ensure water/liquids do not enter the charging socket. Thoroughly dry before use.

To clean the magnetic mating surfaces of the charging connector and charging terminal use a clean and dry cotton swab or a soft bristle brush. Make sure that no cotton fibres or bristles are left on the surfaces.

Never apply water or any other cleaning liquid to the magnetic mating surfaces.

6 Limitations on Use

Intended life:

A local risk assessment should be carried out based upon activity and usage.

Lifting loads:

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment:

Avoid exposing the device to corrosive elements such as water, acids and other liquids.

Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.

Keep magnetic or electrically conductive materials away from the magnetic mating surfaces of the charging connector and charging terminal.

Exclusively for use between -10 °C and 50 °C (14 °F to 122 °F).

Suitable for outdoor use


7 Battery Charging

WARNING: The device contains Li-ion type rechargeable batteries. Please read and comply with the safety information. Failure to do so may cause the battery to become hot, explode, or ignite and cause serious injury.

Please advise the user of the charging information in this section. The batteries are not field replaceable and will be replaced as necessary when serviced. If the device is subjected to temperatures lower than -10°C (14°F) it must be returned to be inspected for possible damage to the battery pack.

1. The battery should only be charged in temperatures between 10 °C to 40 °C (50 °F to 104 °F).

2. Only charge the batteries in the device using the charger supplied with the product. Do not use this charger for charging any other device.

3. The battery charger can be used with an electrical outlet that is rated between 100V and 240V AC, 50/60 Hz.

4. Never connect the charger to the device while still wearing the limb

5. During charging no power is delivered to the device's control system.

6. The knee module may become warm during charging but should never be too hot to touch. There should be no unusual smell, if there is immediately switch OFF [O] the battery, unplug the charger from the electrical supply and contact your service provider.

7. Do not refit the device until the charger has been disconnected and the device has completed the full start-up procedure.

8. Always recharge the device before use if it has not been used for an extended period of time.


7 Battery Charging (continued)

Fully charge the battery before first use.


We recommend that the device should be fully charged every night. It should take 8 hours to charge fully (from fully discharged). We also recommend battery charging whenever convenient.

Step-by-Step Guide:

- 1. Lay the limb on a flat surface.
- 2. Ensure the mating surfaces of the magnetic charging connector and terminal are clean (see Section 5).
- 3. Magnetically join the charging connector to the charging terminal at the posterior of the knee and check alignment is correct (see Section 7.1).
- 4. Switch the Battery On/Off switch to **ON [I]**.
- 5. Check the Battery Charger LED is flashing green to confirm the battery is charging.
- 6. Refer to table below to ensure Limb is charging.
- 7. Disconnect from the charger when charged sufficiently.
- 8. Wait for 30 seconds while observing the device LED changing color as indicated:


Note... Do not move or don the limb while reset is processing otherwise this will prolong the reset time.

- 9. After a successful reset sequence, 2 short beeps will be heard (if this alert is not disabled; see* below).
- 10.If reset sequence is unsuccessful or the LED remains Red, switch the battery **OFF [O]** then **ON [I]** and repeat from step 8.

Battery Charger LED	Charging State
Yellow	Standby
Green > fast flashing	Fast charging
Green > slow flashing	Close to full charge
Green	Maintain (trickle charge)
Yellow > fast flashing	Error (Ensure Battery Switch is ON [I])

Note... If the charger LED is yellow and the Limb is making a clicking noise check that the Battery switch is in the **ON [I]** position.

*The user can enable or disable the audible reset confirmation by pressing and holding the [+] and [-] buttons simultaneously for 2 seconds. The change will be confirmed by 2 short beeps.


7.1 Magnetic Charger Connection

The charging connector can be connected to the device through the charging terminal located at the posterior of the knee.

The connector and the terminal magnetically lock together and correctly self-align on connection (Figure 1).

Note... The connector and the terminal are not polarized so orientation is not critical (Figure 2).


-  Please check for any misalignment between mating surfaces as this may result in no charging. This can happen especially when the charging setup puts strain on the charging cable (Figure 3).
-  Before charging ensure charging surfaces are clean, dry and free of any debris.
-  Do not allow any ferrous material to attach itself to the charging surfaces.
-  Keep the magnetic surfaces away from materials that may be adversely affected by their proximity to strong magnetic fields (e.g. watches, pacemakers, bank cards or magnetically stored media).
-  Prevent build up of debris on the charging surfaces as this may prevent charging or damage the charger. Disconnect the charger from the mains socket and remove any debris before reconnecting the charger to the device.

7.2 Battery Status

To check the battery condition press and **hold** the [-] button on the knee, one short beep will be heard and the LED will indicate the battery condition.

Device

Device LED

Battery Charge Indicator			
LED		Beeps	Battery Condition
Green	3 flashes	--	fully charged
Green	2 flashes	--	good
Green	1 flash	--	ok
Yellow	slow flash	8 x Short beeps	low charge
Red	slow flash	5 x Long beeps	no charge

When fully charged the device's battery is intended to provide up to 3 days of normal use (conditional on the type of use). We recommend that the knee be fully charged every night. Switching off the device when not in use will also prolong the period of use

Note... The device will take 30 seconds to reset after it is switched back **ON** [I] ; refer to Section 7 'Battery Charging (continued)' for the reset sequence.

8 Alignment

The device will need to be adjusted for the height of the user, and the footwear to be used.

Follow these instructions to adjust the length of the pylon in the device.

WARNING

Take care not to damage the internal connecting cable between the knee and ankle modules when assembling or disassembling the device.

8.1 Assemble the Complete Limb

Ensure the connecting cable is not trapped, carefully slide the knee onto the top of the pylon tube with a slight twisting action (there is no need to connect the cable at this stage).


Ensure that the pylon is fully inserted into the knee housing then align the knee and ankle and tighten the clamp bolt to 10 Nm Torque.

Take the measurement required for the user.

8.2 Adjust the Pylon Length

Undo the clamp screw and taking care not to damage the connecting cable remove the Knee module.

To release the pylon undo the 4 pyramid set screws on the ankle module, carefully ease the shin cable and grommet from the pyramid and remove the pylon.

8.3 Cut the Pylon to Length


WARNING

Do NOT attempt to cut the pylon to length without completely removing it from the limb. Failure to do so will result in the connecting cable being damaged.

Using the pylon cutting jig to hold the pylon while cutting it to the required length.

Note that using the cutting jig ensures that the Pylon is not cut shorter than the minimum recommended length.

De-burr and remove any sharp edges from the cut end of the pylon and ensure any debris is cleaned off. De-grease it using a suitable de-greasing agent.

Pylon
Cutting Jig
941257

8.4 Re-assemble the Limb

Connect the shin extension cable to the ankle module cable (for the shortest pylon the extension cable will not be required).



To fit the pylon to the ankle, thread the connecting cable through the small hole in the bottom of the pylon pyramid and gently push the pylon down into the ankle housing ensuring it engages the grommet seal.

Tighten the ankle set screws to secure the pylon.

Lay the ankle and knee flat on the work top and connect the knee module cable to the ankle module cable.

Taking care not to damage the connecting cable carefully slide the knee onto the top of the pylon tube with a slight twisting action.

Align the knee and ankle and tighten the clamp screw.

Clamp Screw



10 Nm

8.5 Bench Alignment

Set the socket in 5° flexion plus any measured flexion contracture.

Tilt Setting: Ankle

Set the ankle in 3° dorsiflexion and 6° plantar flexion. Check that the ankle is not on the tilt limits and that the user is standing in a comfortable upright position before commencing the dynamic alignment.

Ensure that when the device is in full flexion no part of the socket or associated componentry comes into contact with the knee.

*Approximate ratio

Align with shoe on and foot fully plantar flexed.

9 Device Set-up

WARNING

Programming should only be carried out by a suitably qualified practitioner. Incorrect configuration could result in an accident and severe injury.

Prior to setting up a device, install the Linx Programming App on a suitable iOS or Android programming device, e.g. phone or tablet that supports **Bluetooth®** (see Appendix 1 for iOS and Appendix 2 for Android).

The device is set-up, i.e. the limb's onboard program settings are configured and fine-tuned to suit a specific user's requirements, using the offboard Linx Programming App over a **Bluetooth®** connection (see Section 9.1).

However if access to the App is unavailable, then the device may be set-up directly using Button Programming via buttons on the back of the knee (see Section 9.2).

Before setting up the device ensure that the battery is charged and it is switched on.

9.1 Using the Linx Programming App

These instructions describe how to use the Linx Programming App to set-up the device

The App is connected to the device via a **Bluetooth®** module in the limb.

The App allows the Clinician to:

- Set-up the device for an individual user
- Fine-tune the settings
- View the sensor readings dynamically (via the **Monitor** screen)
- Save and restore knee settings.

Setting up a device for the first time

The device will need to be adjusted to suit the height of the user and the footwear to be used, see Section 8.


To launch the Linx Programming App:

1. Click the app icon on the home screen.

2. When the welcome screen appears accept the license agreement and enter your Password (specified during authorization).
3. **Linx Connect** screen appears.
4. Make the **Bluetooth®** connection between the app and the device (see Section 9.1.2).

First time use of the App shows the following screen:

9.1.1 Linx Programming App Screen Layout

Thereafter, the Hub screen will be shown:

The Hub screen provides easy access to the main sections of the device interface: *[Program]*, *[Fine Tune]*, *[User Modes]*, *[Monitor]* and *[Review]*.

For help tap on the  button on the Hub screen or on the title bar on any screen.

Status Bar

Indicates the device operational mode, yield level, assist level (swing control), plantarflexion level, dorsiflexion level and battery level.

9.1.2 Making a Bluetooth® Connection: Connect to the Device

Switch the Battery On/Off switch to **ON** [I] and wait until the reset sequence has been completed successfully.

Press and hold the [+] button on the limb until the LED is purple.

Press the **[Find New Limb]** button.

Select the device from the list of limbs found (pick the entry with the strongest signal if more than one limb is found).

When connected press and hold the [+] button on the limb and wait for the Hub screen to be displayed, or for the **Program** screen to be displayed if this is the first time the App has been used.

Bluetooth® Communication (Press and hold [+] for ½ second)	
LED	Status
Purple	Bluetooth® module enabled
Blue	Connected to App Device

What to do if you cannot connect:

- Check device is switched on.
- Check battery level.
- Ensure **Bluetooth®** is enabled on limb (purple LED).
- Check **Bluetooth®** is enabled on programming device.
- Move closer to limb.
- Wrong limb selected: select another limb from list of limbs.

9.1.3 Setting up the Device

Program

1. Select [**Program**] from the Hub.

2. Check bench alignment:

It is important the device is aligned correctly before calibrating (see Section 8), then press:

Next

3. Adjust yield resistance for sitting in a chair, then press:

Next

4. Press [**Start Walk**] and walk until two beeps are heard. (Adjust dynamic alignment and recalibrate if required).

If realignment is not required, press [**Continue**].

5. Adjust foot plantarflexion and dorsiflexion settings as required, then press [**Continue Limb Calibration**].

6. Press [**Continue Walk**] and walk until 2 beeps are heard. Extra beeps will be heard upon saving of settings.

Programming is now complete, but you can now return to the hub screen and fine tune settings as necessary, and/or save your settings, etc.

Fine Tune

1. Select [**Fine Tune**] from the Hub.

Fine Tune allows refinement of the following settings:

- **Walking Speeds**
- **Stance Settings**
- **Ramp Settings**
- **Stairs Descent**
- **Standing Support.**
- **Foot Settings**

Save and Restore

1. Tap on [**Save and Restore**] button from the Hub:

2. Tap on [**Save from Limb**] button to save settings from connected limb.
3. Select a table entry to review and restore previously saved settings.

9.2 Button Programming

9.2.1 Overview

Button programming, via the [+] and [-] buttons at the rear of the knee, provides an alternative way to program the device in the event of the programming app being unavailable.

There are four limb adjustment/calibration procedures that can be carried out using button programming:

- Yield adjustment, via **Yield Adjust** mode (LED illuminates red)
- Stance release calibration, via **Stance Release Calibration** mode (LED illuminates blue)
- Walk calibration, via **Walk Calibration** mode (LED illuminates yellow)
- Foot calibration, via **Foot Calibration** mode (LED illuminates purple)

Note... There is also a Factory Reset mode available via a separate button-programming mode for resetting all parameters to the factory default settings.

9.2.2 Program Description

The program is in two separate parts:

- **Button-programming mode for calibration/adjustment.**
- **Button-programming mode for factory reset.**


Note... Each part can only be accessed by a specific sequence of button presses.

The button-press sequences for accessing the **Button-programming mode for calibration/adjustment** and carrying out all four calibration/adjustment procedures are shown in Flow Diagram A on page 24 and described in Section 9.2.3 on page 25.


The button-press sequences for accessing the **Button-programming mode for factory reset** and carrying out a factory reset are shown in Flow Diagram B on page 27 and described in Section 9.2.4 on page 27.

Refer also to the Key below.

Key to Flow Diagrams

		Red = Blue = Yellow = Purple = Flashing* = Green = Green flashing = Off = Yield Adjust Mode Stance Release Calibration Mode Walk Calibration Mode Foot Calibration Mode (or Factory Reset) Prepare to Enter... (*color = as modes above) Part of Power On sequence Calibration complete Exited programming mode	
LED and [+] and [-] button location			
Symbols		Button Press Symbols	Audible Alert Symbols
Program function =		Press [-] once =	1 F
Clinician Activity =		Press [+] and [-] simultaneously =	2
Normal flow direction =		Press [+] once =	1
Optional flow direction =		Press and long hold [-] =	1
			2

Flow Diagram A—Button-Programming Mode for Calibration/Adjustment

Once in a button-programming mode, the [+] and [-] buttons are used to navigate through the various stages of the particular program. Short audible beeps are used to confirm button presses.

LED colors/effects and medium or long audible beeps are used to indicate current position within a program, and also act as prompts for the clinician to carry out certain activities.

Navigation Notes

To enter **Button-Programming mode for Calibration/Adjustment**:

- Press both [+] and [-] buttons simultaneously *twice* in quick succession.

To enter **Button-programming mode for factory reset**:

- Press both [+] and [-] buttons simultaneously *three times* in quick succession.

Note... Each simultaneous press of the [+] and [-] buttons together has to be held long enough for the device to register the combined press as a single action (confirmed with two short beeps), before releasing.

To move forward or backward through stages of the program:

- Press of the [+] or [-] button respectively once.

After each single press is registered, the device confirms with one short beep.

Entering a “Prepare to Enter...” mode is confirmed by two medium or two long beeps, and the LED flashing in the color associated with the corresponding adjustment/calibration mode.

Entering an adjustment/calibration mode is indicated by between 3 and 6 medium beeps depending on the mode, and by the LED illumination changing from flashing to steady.

To exit button programming at any point in the program:

- Press and hold [-] button until confirmed by two long beeps, and then release. The LED extinguishes to confirm exit.

If there is no button activity within a certain time period, the program (adjustment/calibration mode, or factory reset mode) will automatically time-out and exit.

9.2.3 Calibration/Adjustment

Normal Flow Sequence

Refer to Flow Diagram A.

Power on the device and, following the solid arrow flow, go through each programming stage in the order given to enter **Button-programming mode for calibration/adjustment**, and complete all four calibration/adjustment procedures and their associated activities.

Note that each calibration/adjustment procedure has an associated LED color, and each procedure has two distinct programming modes:

- **Prepare to Enter...** mode, indicated by the LED flashing in the color associated with that procedure
- **Calibration/adjustment** mode, indicated by the LED illuminating steady in the color associated with that procedure.

Note also that the activity associated with the “Prepare to Enter...” mode must be completed before moving on to the associated calibration/adjustment mode and carrying out its associated calibration or adjustment.

When all four calibration/adjustment procedures have been successfully completed in order, the program enters the **Calibration Complete** stage and the LED flashes green.

To finish and exit Button Programming:

- Press [+] once. (or press and hold [-])
Two long beeps sound and the LED extinguishes to confirm exit.

Skipping a Procedure

You can break from the normal program flow sequence and skip a calibration/adjustment procedure e.g. if you have already successfully completed that procedure beforehand.

To skip a procedures and move to another one:

- Press the [-] button to navigate one mode backwards through the cycle of “Prepare to Enter” modes, following the dashed arrow path in Flow Diagram A.
- Repeat as necessary to navigate to the required “Prepare to Enter” mode.

WARNING

You must have completed all four adjustment/calibration procedures before exiting the program, otherwise the limb will not be fully functional

Note... There is also an automatic time-out feature, such that if no button activity or programming is observed within a certain time period, the device automatically exits button programming mode.

If you want to exit before all four calibration/adjustment procedures are completed, you can do this at any point in the flow by pressing [-] once and holding for a time until two long beeps sound and the LED extinguishes conforming exit, however the limb will not be fully functional.

Yield Adjustment Notes

The Yield Adjustment procedure requires a manual operation to adjust the **Yield Level** (resistance) of the limb for sitting, and when complete you must press [+] to move to the next **Prepare to enter...** stage, whereas for the other three calibration procedures, calibration is an automatic operation initiated when the program is in a calibration mode and the user begins to walk.

When calibration is complete, the program automatically moves onto the next **Prepare to enter...** , or **Calibration Complete** stage, as appropriate.

In **Yield Adjust** mode, (LED illuminated steady red) the limb **Yield Level** can be adjusted manually up or down in stepped increments as follows:

To Decrease Yield Level:

With the user standing, supported (e.g. by holding on to a walking rail):

- Flex the knee through approximately 15°, maintain the angle momentarily then fully extend it to register a single step decrease in **Yield Level**.
(Confirmed by one short beep);
- Repeat as necessary until the required resistance is achieved.

To Increase Yield Level:

With the user standing, supported:

- Flex the knee through approximately 40°, maintain the angle momentarily then and then fully extend it to register a single step increase in **Yield Level**.
(Confirmed by two short beeps);
- Repeat as necessary until the required stiffness is achieved.

9.2.4 Factory Reset

WARNINGS

-  1. **Factory Reset will delete all calibration data and return settings to the factory defaults. The knee will be in a free state with no stance support during Factory Reset.**
-  2. **After a factory reset do not attempt to walk on the limb until it has been suitably reprogrammed.**

Refer to Flow Diagram B.

To enter the **button-programming mode for factory reset**:

- Power on the device, (wait for the power on LED color sequence to complete, ending in LED OFF then,
- Press and hold **[+]** and **[-]** simultaneously *three times* in quick succession (waiting for two short beeps to sound before releasing each press).

The limb now enters **Prepare for Factory Reset** mode indicated by three medium beeps and LED flashing alternating purple and red.


Note... This is the last point in the flow at which you can cancel Factory reset (to cancel and exit, press and hold **[-]**, two long beeps confirms exit.)

To proceed with Factory Reset:

- Ensure the user is sitting down with their knee flexed.
- Press the **[+]** button to move to the **Start Factory Reset** stage (5 short beeps and LED steady purple confirms).
- Flex and extend the leg until two long beeps are heard.

Factory reset is now complete and the program exits Factory Reset mode automatically and the LED extinguishes.

Flow Diagram B—Button-Programming for Factory Reset

10 User Modes

The User Modes are an optional feature of the device and should only be enabled if the user can operate them competently and safely. The user should be advised to take extra care when using these modes.

User access to the User Modes can be **enabled** /**disabled** via the Linx Programming App: by selecting [**User Modes**] from the Hub.

When enabled, the user can then switch the required User Mode On or Off using the [-] button on the device. During set-up the User Modes can be tested in the App.

The device has three User Modes:

1. Flexion Lock Mode

Sets the limit of flexion to a predetermined angle set by the clinician.

The maximum lock angle is 45° from fully extended.

Note... Before the user can **activate** a User Mode it must be **enabled** in the user interface.

10.1 Flexion Lock Mode

To enable/disable Flexion Lock Mode via the Linx Programming App:

- Select [**User Modes**] from the Hub.

To enable/disable Flexion Lock Mode:

- Tap Flexion Lock [**Enable**] to enable.
- Tap Flexion Lock [**Disable**] to disable.

To adjust the lock angle:

- Flex the knee to the required lock angle and tap on the [**Set Knee Angle**] button.

To test Flexion Lock using the app:

- Swipe right to expose the User Modes test buttons then tap on the [**Flexion Lock Test**] button.

2. Cycling (Free swing) Mode

Releases stance, the limb will then swing freely.

WARNING: In Cycling Mode there is no stance active, this mode must be de-activated before resuming walking to avoid falling and possible injury.

3. Training Mode

In this mode, the limb has a predefined fixed pneumatic resistance setting, set by the clinician, which does not adapt to changes in walking speed. Also, Standing Mode is not available or operational in Training Mode.

User Activation of Flexion Lock Mode

To activate/de-activate Flexion Lock use the [-] button on the knee module, as follows:

To switch On Flexion Lock:

- Fully extend the knee.
- Press the [-] button 3 times (each time pressing until the device responds with a single short beep).
- 3 medium beeps then sound to confirm that Flexion Lock is ON.

To switch Off Flexion Lock:

- Press and hold the [-] button.
- 2 long beeps then sound to confirm that Flexion Lock is OFF and normal operation has resumed.

WARNING: After exiting Flexion Lock Mode, always test/confirm that the limb is operating normally again before walking.

10.2 Cycling Mode

WARNING: In Cycling Mode there is **NO yield resistance**.

To enable/disable Cycling Mode via the Linux Programming App:

- Select [*User Modes*] from the Hub.

To enable/disable Cycling Mode:

- Tap Cycling Mode [*Enable*] to enable.
- Tap Cycling Mode [*Disable*] to disable.

To test Cycling Mode using the app:

- Swipe right to expose the User Modes test buttons, ask the User to flex the knee to 45 degrees or greater (typically by sitting) then tap on the [*Cycling Mode Test*] button.

WARNING: After exiting Cycling Mode, always test/confirm that the limb is operating normally again before walking.

10.3 Training Mode

WARNING: In Training Mode, the limb has a predefined fixed pneumatic resistance setting which does not adapt to changes in walking speed. Also, Standing Mode is not active in Training mode.

To enable/disable Training Mode via the Linux Programming App

- Select [*User Modes*] from the Hub.

To enable/disable Training Mode:

- Tap Training Mode [*Enable*] to enable.
- Tap Training Mode [*Disable*] to disable.

To test Training Mode using the app:

- Swipe right to expose the User Modes test buttons, then tap on the [*Training Mode Test*] button.

WARNING: After exiting Training Mode, always test/confirm that the limb is operating normally again before walking.

User activation of Cycling Mode

To activate/de-activate Cycling Mode use the [-] button on the knee module, as follows:

To switch On Cycling Mode:

WARNING: In Cycling Mode there is **NO yield resistance**.

- Press the [-] button 4 times (each time pressing until the device responds with a single short beep).
- 4 medium beeps then sound to confirm that Cycling Mode is ON.

To switch Off Cycling Mode:

- Press and hold the [-] button.
- 2 long beeps then sound to confirm that Cycling Mode is OFF and normal operation has resumed.

User activation of Training Mode:

To activate/de-activate Training Mode use the [-] button on the knee module, as follows:

To switch ON Training Mode:

- Press the [-] button 5 times (each time pressing until the device responds with a single short beep).
- 5 medium beeps then sound to confirm that Training Mode is ON.

To switch OFF Training Mode:


- Press and hold the [-] button.
- 2 long beeps then sound to confirm that Training Mode is OFF and normal operation has resumed.

11 Technical Data

Operating and Storage Temperature Range:	-10 °C to 50 °C (14 °F to 122 °F)
Charging Temperature Range:	10 °C to 40 °C (50 °F to 104 °F)
Component Weight (size 26):	2.6 kg (5 lb 12 oz)
Activity Level:	2, 3, 4
Maximum User Weight:	125 kg (275 lb)
Proximal Alignment Attachment:	T-Slot
Range of Hydraulic Ankle Motion (excludes additional range of motion provided by heel and toe springs)	6 degrees plantarflexion to 3 degrees dorsiflexion
Build Height: (See diagram)	470–565 mm (18.5–22.2")
With Long Pylon Kit 339965 (available separately)	up to 679 mm (26.7")

Battery:	
Battery Type	Rechargeable Li-ion
Battery life	3 days normal use from a full charge. Charge daily
Battery Charger:	
Charging time to full charge	8 hours
Input Voltage	100V to 240V AC
Frequency	50 to 60Hz

Fitting Length

Transport & Storage:

1. Ensure that the device is stored with the knee fully extended when not in use.
2. If the limb is to be stored upright then it must be secured to prevent it falling over or flexing.
3. If the device is shipped either on its own or built into a limb it should be maintained in an extended position and packaged in a suitable box to allow this.
4. If a unit is to be returned, contact Blatchford prior to sending if any obvious physical damage has occurred to the battery pack.
5. For transport purposes the unit must be switched off, and packed securely so that
 - a) it cannot be inadvertently switched on and become operative while in transit;
 - b) the unit is thoroughly protected from damage in transit.

Appendix 1 Installing the Linx Programming App (iOS)

To use the App you will require an **Authorization Code** from your Blatchford representative.
To apply for an **Authorization Code**, first you will need to install the App.

Note... Authorization Codes will only be issued to suitably qualified practitioners who have successfully completed a Blatchford approved Linx training course.

System Requirements:

Any iOS programing device that supports iOS 9.0 or later.

1. Install the App

From the Apple® App Store on your iOS device search for Linx Programming App.

Note... If you are installing onto an iPad ensure that iPhone Only is selected in the App Store search filter.

Download the app from the Apple® App Store, then follow the instructions opposite to obtain an **Authorization Code**.

2. To Obtain an Authorization Code

Launch the App:

- 2.1 Tap on the app icon on the home screen.
- 2.2 When the **LinxV2 Hub** screen appears, click on the **[Authorize]** button.
- 2.3 On the **Authorization Method** page, tap the **[by Code Entry]** button.

- 2.4 Send the 8-digit **Authorization Request Code** with your name and organization details to the appropriate email address for your region from the following:

Region	Support mailbox
US	support@blatchfordus.com
UK	supportUK@blatchford.co.uk
DE	info@blatchford.de
FR	support@blatchford.fr
RoW	support_ROW@blatchford.co.uk

or contact your allocated Blatchford representative.

- 2.5 You will then be contacted by a Blatchford representative who will supply you with a corresponding 12-digit Authorization Code.
- 2.6 Enter the code in the box provided and click on the **[Submit]** button
- 2.7 After accepting the license agreement create a password (this will be needed for all subsequent times you enter the App).

Appendix 2 Installing the Linx Programming App (Android)

To use the App you will require an **Authorization Code** from your Blatchford representative.
To apply for an **Authorization Code**, first you will need to install the App.

Note... Authorization Codes will only be issued to suitably qualified practitioners who have completed a Blatchford approved Linx training course.

Minimum System Requirements

Prior to installation ensure your programming device hardware and software specifications satisfy the following minimum systems requirements:

Hardware


- Processor: ARMv7 processor or higher
- RAM: 1.0 GByte
- Storage: 10 MByte
- Screen Size: For optimum app display and performance a screen size of 7" is recommended.

Supported Operating System

- Version 8.0 (1st Oreo release), or higher

1. Install the App

- 1.1 Remove any older version of this software from your Android programming device before proceeding.
- 1.2 Insert the USB stick part no. 236520V2A into your PC.
- 1.3 Send the **<app-release>.apk** file [Fig. 1] on the USB to your Android programming device at a location of your choice (e.g. "My Files" directory). You can do so by:
 - connecting your programming device to your PC and using your file browser to locate the file on the USB and copying across to the Android programming device,

- e-mailing the file as an attachment from your PC to your Android programming device and then downloading the received attachment, or
- sending the file from your PC over **Bluetooth®** file transfer to the Android programming device.

- 1.4 On the programming device, navigate to the saved **<app-release>.apk** file and tap on Linx app icon [Fig. 1] to begin the installation.
- 1.5 The programming device manager may open a window which asks: "Do you want to install this application?" and gives options Install and Cancel. Tap Install to accept the install.

Note... During the installation process, you may be asked to permit the installation from Unknown Sources. Ensure that this permission is given for this app only

- 1.6 Once installed, the programming device will display an acknowledgement message to say the installation is successful. Press **Open** to open the app (or press **Done** if you don't want to open the app at the time), then follow the instructions in the next section to obtain an **Authorization Code**.

2. To Obtain an Authorization Code

Launch the App:

- 2.1 Tap on the app icon on the home screen.
- 2.2 When the *LinxV2 Hub* screen appears, click on the *[Authorize]* button.
- 2.3 On the *Authorization Method* page, tap the *[By Code Entry]* button.
- 2.4 Send the 8-digit *Authorization Request Code* with your name and organization details to the appropriate email address for your region from the following:

Region	Support Mailbox
US	support@blatchfordus.com
UK	supportUK@blatchford.co.uk
DE	info@blatchford.de
FR	support@blatchford.fr
RoW	support_ROW@blatchford.co.uk

or contact your allocated Blatchford representative.

- 2.5 You will then be contacted by a Blatchford representative who will supply you with a corresponding 12-digit Authorization Code.
- 2.6 Enter the code in the box provided and click on the *[Submit]* button.
- 2.7 After accepting the license agreement create a password (this will be needed for all subsequent times you enter the App).

FCC Compliance Statement

This device complies with part 15 of the FCC Rules.

Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

Note... Contains pre-approved **Bluetooth®** module: FCC ID: PVH0946 IC: 5325A-0946

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/ TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation. This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Responsible party: Blatchford Inc.

1031 Byers Road, Miamisburg, OH 45342

800-548-3534 (toll free) | 937-291-3636

customerservice@blatchfordus.com

RSS Compliance Statement

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

1. this device may not cause interference, and
2. this device must accept any interference, including interference that may cause undesired operation of this device.

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Caution: Exposure to Radio Frequency Radiation.

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website Blatchford Inc. 1031 Byers Road, Miamisburg, OH 45342.

800-548-3534 (toll free) | 937-291-3636 customerservice@blatchfordus.com

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices.

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Warranty

This device is warranted for 36 months, the foot shell for 12 months, and the glide sock for 3 months.

Extended warranties of an additional 12, 24, or 36 months are available.

A service between 30 and 36 months is required for all extended warranties.

For the maximum 36-month extended warranty, a second service must be carried out between 54 and 60 months.

We reserve the right to void extended warranties if the device is not serviced at the required intervals.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

The above may vary by market; consult your local representative for details.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

This symbol indicates that the product contains electrical/electronic components and/or batteries that should not be disposed of in general waste or be incinerated at the end of the product's life.

At the end of the product's life, all electrical/electronic components and/or batteries should be recycled or disposed of in accordance with the current regulations for the handling of WEEE (Waste Electrical and Electronic Equipment), or equivalent local regulations. The remainder of the product should also be recycled where possible in accordance with local waste recycling regulations.

To help prevent potential harm to the environment or to human health from uncontrolled waste disposal, Blatchford offers a take-back service. Please contact Customer Services for details.

Retaining the Packaging Label

You are advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgments

Linx and Blatchford are registered trademarks of Blatchford Products Limited.

Android™ trademark is a registered trademark of Google LLC. All rights reserved.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG Inc., and any use of such marks by Blatchford is under License. Other trademarks and trade names are those of their respective owners.

iPad, iPhone, and Apple are trademarks of Apple Inc., registered in the U.S. and other countries.

iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license.

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© Blatchford Products Limited 2026. All rights reserved.
Endolite India Ltd.
Blatchford Europe GmbH
Blatchford Products Ltd.