



# Adult Footwear Collection

# Adult Footwear Collection

**Our new footwear range has been designed using the latest materials and manufacturing methods to produce the highest quality footwear.**

Mendivil modular and stock footwear are designed to focus on the prevention and treatment of foot pathologies. Continuous evaluation of product lines are carried out by INESCOP resulting in innovative high quality products.

The Stock Therapeutic range from Fior & Gentz and YDA is designed specifically for diabetic feet.

All footwear contained in this catalogue is distributed by Blatchford in accordance with supplier agreements.

To meet the clinical requirement of a patient, Blatchford have the right to modify CE marked footwear and by default this will be labeled as a custom made device.

Blatchford maintains a QMS, which is ISO 9001 and ISO 13485 accredited.

## Contents

	Page
<b>Mendivil Stock Footwear</b>	4
Ladies	6
Gents	10
<b>Cadense Stock Footwear</b>	12
<b>PerfectStride</b>	14
<b>YDA Stock Footwear</b>	18
<b>Fior &amp; Gentz Therapeutic Stock Footwear</b>	20
<b>Camino Modular Collection</b>	34
Measurement Guide	36
Last Additions	37
Materials	39
Sole Units and Punch Patterns	44
Unisex House	46
Ladies	48
Gents	62
<b>Measurement Guides</b>	74
<b>Fitting Problems and Solutions</b>	76

# Adult Footwear Collection **Information**

The range includes Camino Modular and Made to Measure footwear, plus stock footwear from Mendivil, Cadense, PerfectStride, YDA and Fior & Gentz.

Modular footwear is available in a larger style and size range than stock footwear and is made to a trial fit stage as default. Different last sizes are available and additions can be added to create a Modular shoe for optimal treatment of more complex pathology. Three clear options are detailed in the measurement guide on pages 36-37, which simplifies the ordering process. The modular range is shown on pages 48-72.

The stock range allows you to order a completed shoe ready to supply to your patient. The stock ranges are on pages 4-33. We currently have one gents' style available for delivery within 24 hours. All other styles can be ordered with delivery to the clinic within five working days. While stocked footwear enables faster turnaround times, it is limited to specific foot types, as modifications to the last size cannot be accommodated.

All Camino styles in this catalogue can be requested for the production of Made to Measure footwear. However, these are not the limit and other styles can be ordered. The styles that are marked as being ideal for Made to Measure footwear have plain fronts and will allow for some size differences without making it too obvious. For Made to Measure footwear a draft and measurements to British Standards must be provided, along with a cast, diagrams and photos as necessary. Please note that complex Made to Measure footwear may require non-standard sole and heel units, stiffeners and insoles.

A full adaptation service is available for additions such as rocker soles, flares, sockets etc on all types of footwear excluding the Cadense range.

## 1. Ordering and Enquiries:

Ordering footwear from Blatchford is simple and easy. You can contact our Customer Services directly:

Telephone: +44 (0) 114 2637900

Email: [orthoticscs@blatchford.co.uk](mailto:orthoticscs@blatchford.co.uk)

Facsimile: +44 (0) 114 2637901

Post: Blatchford Limited, 11 Atlas Way, Atlas North, Sheffield, S4 7QQ.

## 2. Delivery Promise:

Our goal is to provide a fast and efficient delivery. We aim to despatch orders for Made to Measure footwear in 30 working days and Modular in 25 working days to first fitting. Footwear returned from fitting to completion will be processed in 10 working days for both Made to Measure and Modular. This does not apply to where specific delivery schedules are in place.

## 3. Catalogues/Product Information:

If you require additional product information please contact our customer service advisors who will be happy to help with your enquiry.

## 4. Shortages/Damaged Goods:

If any items are received damaged or there are shortages in the delivery this must be notified in writing to customer services within 5 days of receipt of goods.

## 5. Returns:

In the unlikely event of any footwear being supplied faulty or having been incorrectly supplied, please contact customer services so that the correct goods are despatched immediately. We are happy to accept Stock Footwear that is returned for credit if they are not suitable for the patient.

### **For items to be accepted for credit the following conditions apply:**

- Stock footwear must be unworn and returned in less than 30 days from the date of the delivery
- The product must not have been adapted in any way
- Split sized footwear direct from stock is returnable for credit
- Are returned in the original packaging
- Are in a re-saleable condition

### **Non-returnable:**

- Made to Measure and Modular footwear are non returnable for credit

Blatchford reserves the right to reject items returned for credit if the above conditions are not met. Blatchford is unable to accept return of items that are 'made to order' for credit, unless the goods are faulty or it can be demonstrated they have been damaged in transit. Products and prices are subject to variation at any time and are subject to availability. E & OE.



# Stock Footwear Collection



# Ladies Mendivil Stock Footwear

All models are made with stretch materials which adapt to foot prominences bunions and hammer toes, they have no internal seams which makes them ideal for diabetics.

5 Day Delivery from Stock in Spain



## Palma (Product Code 86220-GR+size)

**Colour:** Grey

**Size Range:** 36 - 42

**Last Width:** X-Wide 68

Ladies models with uppers made of elastic fabric adaptable to prominences of the foot with bunions or hammer toes. Lightweight, non-slip and long-lasting EVA sole units. Indicated for diabetic foot and arthritis.

5 Day Delivery from Stock in Spain



## Vigo (Product Code 86240-BK+size)

**Colour:** Black

**Size Range:** 36 - 42

**Last Width:** X-Wide 68

Ladies sport footwear with seamless uppers made of top quality leathers and stretch fabric adaptable to prominences of the foot with bunions or hammer toes. Lightweight, non-slip and long-lasting EVA sole units. Indicated for diabetic foot and arthritis.

5 Day Delivery from Stock in Spain



## Granada (Product Code 86230-NV+size)

**Colour:** Blue

**Size Range:** 36 - 42

**Last Width:** X-Wide 68

Ladies models with uppers made of elastic fabric adaptable to prominences of the foot with bunions or hammer toes. Lightweight, non-slip and long-lasting EVA sole units. Indicated for diabetic foot and arthritis.

5 Day Delivery from Stock in Spain



## Malaga (Product Code 86250-BR+size)

**Colour:** Brown

**Size Range:** 36 - 42

**Last Width:** X-Wide 68

Ladies sport footwear with seamless uppers made of top quality leathers and stretch fabric adaptable to prominences of the foot with bunions or hammer toes. Lightweight, non-slip and long-lasting EVA sole units. Indicated for diabetic foot and arthritis.

5 Day Delivery from Stock in Spain



## Zamora (Product Code 86260-LG+size)

**Colour:** Lago

**Size Range:** 36 - 42

**Last Width:** X-Wide 68

Ladies sport footwear with seamless uppers made of top quality leathers and stretch fabric adaptable to prominences of the foot with bunions or hammer toes. Lightweight, non-slip and long-lasting EVA sole units. Indicated for diabetic foot and arthritis.

5 Day Delivery from Stock in Spain



## Osuna (Product Code 86270-BE+size)

**Colour:** Beige

**Size Range:** 36 - 42

**Last Width:** X-Wide 68

Ladies sport footwear with seamless uppers made of top quality leathers and stretch fabric adaptable to prominences of the foot with bunions or hammer toes. Lightweight, non-slip and long-lasting EVA sole units. Indicated for diabetic foot and arthritis.

5 Day Delivery from Stock in Spain



## Valencia (Product Code 85430-BK+size)

**Colour:** Black

**Size Range:** 35 - 42

**Last Width:** Wide 42

Ladies footwear with seamless uppers made of natural stretch leathers and/or stretch fabric adaptable to prominences of the foot with bunions or hammer toes. Lightweight, non-slip and long-lasting EVA sole units. Indicated for diabetic foot and arthritis.

# Ladies Mendivil Stock Footwear

All models are made with stretch materials which adapt to foot prominences bunions and hammer toes, they have no internal seams which makes them ideal for diabetics.

5 Day Delivery from Stock in Spain



## Getxo (Product Code 86290-GR+size)

**Colour:** Grey

**Size Range:** 35 - 41

**Last Width:** Wide 47

Ladies footwear with seamless uppers made of natural stretch leathers and/or stretch fabric adaptable to prominences of the foot with bunions or hammer toes. Lightweight, non-slip and long-lasting EVA sole units. Indicated for diabetic foot and arthritis.

5 Day Delivery from Stock in Spain



## Girona (Product Code 1820-NV+size)

**Colour:** Navy

**Size Range:** 35 - 42

**Last Width:** Wide 47

Ladies sandal made of natural leather and stretch fabric over bunion pockets (medial and lateral) adaptable to prominences of the foot with a removal 10 mm shock absorbing insole. PU sole unit.

5 Day Delivery from Stock in Spain



## Melilla (Product Code 86300-BK+size)

**Colour:** Black

**Size Range:** 35 - 41

**Last Width:** Wide 47

Ladies footwear with seamless uppers made of natural stretch leathers and/or stretch fabric adaptable to prominences of the foot with bunions or hammer toes. Lightweight, non-slip and long-lasting EVA sole units. Indicated for diabetic foot and arthritis.

5 Day Delivery from Stock in Spain



## Elda (Product Code 86280-BR+size)

**Colour:** Brown

**Size Range:** 35 - 41

**Last Width:** Wide 47

Ladies footwear with seamless uppers made of natural stretch leathers and/or stretch fabric adaptable to prominences of the foot with bunions or hammer toes. Lightweight, non-slip and long-lasting EVA sole units. Indicated for diabetic foot and arthritis.

# Mendivil Stock Last

## Ladies Measurement Guide


Please note these are not the same as the Modular Range.

Shoe		Wide 47			Wide 42			X-Wide 68		
Size	Length	Joint Width	Joint Circ.	Instep Circ.	Joint Width	Joint Circ.	Instep Circ.	Joint Width	Joint Circ.	Instep Circ.
35	226	82	230	231	84	236	237	-	-	-
36	232	84	234	236	86	240	242	90	243	245
37	239	86	238	241	88	244	247	92	248	251
38	246	88	242	246	90	248	252	94	253	257
39	252	90	246	251	92	252	257	96	258	263
40	259	92	250	256	94	256	262	98	263	269
41	266	94	254	261	96	260	267	100	268	275
42	272	-	-	-	98	264	272	102	273	281

Measurements for all sections are shown in mm.

# Gents Mendivil Stock Footwear

## Featuring Outlast® Technology



How Outlast® Technology Works:

## Mendivil Stock Last Gents Measurement Guide

Please note these are not the same as the Modular Range.

Shoe		Wide 94			X-Wide 95		
Size	Length	Joint Width	Joint Circ.	Instep Circ.	Joint Width	Joint Circ.	Instep Circ.
39	252	94	255	258	96	260	263
40	259	96	260	263	98	265	268
41	266	97	265	268	99	270	273
42	272	98	270	273	100	275	278
43	279	100	275	278	102	280	283
44	285	102	280	283	104	285	288
45	293	103	285	288	105	290	293
46	300	105	290	293	107	295	298

Measurements for all sections are shown in mm.

All models have no internal seams which makes them ideal for diabetics.

## Madrid (Product Code 85310-BK+size) 5 Day Delivery from Stock in Spain

Full Outlast® lining.

Gents footwear with seamless uppers made of natural leather and special "Outlast" lining which ensures thermal comfort of the foot. Lightweight, non-slip and long-lasting EVA sole unit. Specially indicated for diabetic foot.



**Colour:** Black  
**Size Range:** 39 - 46  
**Last Width:** X-Wide 95



## Barcelona (Product Code 85550-BK+size) Next Day Delivery Available from UK Stock

Natural stretchable leather apron.

Popular seller.

Gents footwear with seamless uppers made of natural stretch leather adaptable to prominences of the foot. Lightweight, non-slip and long-lasting EVA sole unit. Specially suitable for diabetic and arthritic feet.

**Colour:** Black  
**Size Range:** 39 - 46  
**Last Width:** Wide 94



## Pamplona (Product Code 1850-BK+size) 5 Day Delivery from Stock in Spain

Gents sandal made of top quality leathers with a removable 10mm shock absorbing insole. PU sole unit.

**Colour:** Black  
**Size Range:** 39 - 46  
**Last Width:** X-Wide 95



# Cadense Stock Footwear

Available direct from UK stock in 5 days

Built for  
safer,  
smoother  
walking.

## Cadense - Adaptive Footwear Engineered for Confidence.

Cadense is transforming the way people walk - with shoes specifically designed for individuals experiencing foot drop, MS, stroke recovery, Parkinson's, and other mobility challenges. Built with patented Variable Friction Technology, Cadense Adaptive Shoes offer a breakthrough approach to safer, smoother walking - empowering over 30,000 customers (and counting) to move confidently through life. Designed in California and New York. Backed by clinical research. Loved by doctors and customers alike.



## Cadense Original

### Features

- **Variable Friction Technology:** Patented sole design helps glide during swing phase and stop during stance
- **Adaptive Nylon Pucks:** Retract under pressure for seamless toe-off
- **Easy On/Off:** Hook and loop closure, wide opening, and lightweight build
- **Removable Insoles:** Accommodates custom orthotics
- **Supports AFO use:** Works with most AFOs that fit inside wide width trainers
- **Available Colours:** Shown here



All-Black **Code:** ABK



Black **Code:** BLK



White **Code:** WHT



Slate **Code:** SLT  
(Gents Sizes Only)



Light Blue **Code:** BLU  
(Ladies Sizes Only)

**Size Range:** Ladies 3.5-10.5/Gents 7.5-14.5

**Widths:** Regular and Wide

**Code Example:** SHOR(W/M) (R/W) (Size) (Colour Code)  
**IF Ladies Regular 7.5 All-Black required, code is:** SHORWR7.5ABK

## Cadense Havena Women's Adaptive Slipper

### Features

- **Removable Backstrap:** Wear with the strap for extra support, or remove for a more relaxed fit. No tools required.
- **Adjustable Front Hook and Loop:** A wide front flap allows for a personalised fit.
- **Width:** One flexible width, with adjustability to accommodate a variety of foot shapes.
- **Lining:** Plush shearling lining keeps feet warm and cozy - comfortable barefoot or with socks.
- **Materials:** Crafted with buttery-soft genuine suede for a premium feel and long-lasting wear.
- **Weight:** Approximately 10 oz per slipper (Women's size 8 used as reference).
- **Sole Height:** Approximately 1.5 inches, designed to support our patented Variable Friction Technology while providing cushioning and comfort.



Chestnut **Code:** CNT



Sand **Code:** SND

**Size Range:** Ladies 3.5-10.5

**Widths:** Regular

**Code Example:** HAV(W) (R/W) (Size) (Colour Code)  
**IF Ladies 7.5 Sand required, code is:** HAVW7.5SND



# Cadense Footwear Measurement Guide

Cadense footwear is available in Standard and Wide options.

The tables below are for guidance only and shows the last measurements for Standard sizes.

It converts US sizes to the nearest UK sizes. Please note: not all UK sizes are available.

The upper material is designed to adapt to minor deformities for a comfortable fit.

Wide option does not mean a wider joint, it is more accommodative to provide extra room in the upper for slight swelling or minor toe deformities.

### Our recommendation:

Order the same size the patient normally wears.

If swelling or deformity has led the patient to purchase a larger size than usual, choose Wide option in their original shoe size.

## Measurements

Cadense Size Chart Ladies Regular Last						
Size UK	Length (mm)	Joint Width	Joint Circ.	Instep Circ.	Toe Depth	Heel Width
3.5	230	82	216	215	24	50
4	234	83	219	218	24	51
4.5	238	84	222	221	24	52
5	242	85	225	224	24	52
5.5	246	86	229	228	25	53
6	250	87	232	231	25	54
6.5	254	88	235	234	25	55
7	258	89	238	237	25	55
7.5	262	90	241	240	25	56
8	266	91	245	244	26	57
8.5	270	92	248	247	26	58
9	274	93	251	250	26	58
9.5	278	94	254	253	26	59
10.5	286	96	260	259	27	61

Cadense Size Chart Gents Regular Last						
Size UK	Length (mm)	Joint Width	Joint Circ.	Instep Circ.	Toe Depth	Heel Width
7.5	268	92	248	259	30	60
8	272	93	251	262	30	60
8.5	276	94	254	265	30	61
9	280	95	257	268	30	62
9.5	284	96	260	271	30	63
10	289	97	264	275	31	63
10.5	293	98	267	278	31	64
11	297	99	270	281	31	65
11.5	301	100	273	284	31	66
12.5	310	102	279	290	32	67
13.5	318	105	286	297	32	69
14.5	327	107	292	303	32	70

## Colour Availability

Regular				Wide			
Black	White	Light Blue	All-Black	Black	White	Light Blue	All-Black
✓	✓	✓	✓				
✓	✓	✓	✓				
✓	✓	✓	✓	✓	✓		✓
✓	✓	✓	✓				
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓				
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓				
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓				
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓				
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓				
✓	✓	✓	✓	✓	✓	✓	✓

Regular				Wide			
Black	White	Slate	All-Black	Black	White	Slate	All-Black
✓	✓	✓	✓	✓		✓	✓
✓	✓	✓	✓				
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓				
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓				
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓				
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓				
✓	✓	✓	✓	✓	✓	✓	✓
✓	✓	✓	✓				
✓	✓	✓	✓	✓	✓	✓	✓

# PerfectStride

EXPERIENCE THE FUTURE OF FOOTWEAR



## PerfectStride

Discover the comfort and performance of our PerfectStride shoes, designed to provide the ultimate solution for your daily activities.

PerfectStride creates a unique walking experience and delivers:

### Features

- A confident stride with the slip-safe rubber outsole, ensuring superior grip and stability wherever you walk.
- Stability and support offered by the EVA midsole for a balanced stride.
- An optimal fit, thanks to its generous volume, accommodating multiple foot shapes and sizes.
- Utmost walking comfort and reduced strain, delivered by the extra padded heel counter design.
- High energy return stride due to its Carbon Spring System, revolutionising your every step.
- Recommended for anyone that suffers from diabetes and a vast range of other foot and ankle pathologies.

### Measurements

PerfectStride					
Size EU	Length (mm)	Joint Width	Joint Circ.	Instep Circ.	Toe Depth
36	242	97	250	238	27
37	249	99	255	243	28
38	256	100	259	247	28
39	262	102	264	252	28
40	269	103	268	256	29
41	275	105	273	261	29
42	282	106	277	265	30
43	289	108	282	270	30
44	295	109	286	274	30
45	302	111	291	279	31
46	309	112	295	283	31
47	315	114	300	288	32

## PerfectStride Water-repellent Knit with TPU Film



White TPU, Knit Textile  
**Code:** MESC008-0100-\*Size



Black TPU, Knit Textile  
**Code:** MESC008-0200-\*Size



Cool Grey TPU, Knit Textile  
**Code:** MESC008-0701-\*Size

**Size Range:** 36-47

**Upper:** Water-repellent Knit with TPU Film

**Lining:** Textile (Airmesh)

**Sole:** EVA/Rubber with Carbon Insole

**Footbed:** Removable, 5mm

## PerfectStride Leather with Engraved Tongue



White Leather  
**Code:** MESC007 -1000-\*size



Black Leather  
**Code:** MESC007 -2000-\*size



Blue Leather  
**Code:** MESC007 -8000-\*size

**Size Range:** 36-47

**Upper:** Leather with Engraved Tongue

**Lining:** Textile (Airmesh)

**Sole:** EVA/Rubber with Carbon Insole

**Footbed:** Removable, 5mm

# PerfectStride Hook/Loop



Black

Code: MESC007-\*Size

Size Range: 36-46 - Unisex

Upper: Leather

Lining: Textile

Sole: EVA/Rubber with Carbon Insole

Footbed: 5mm

## Measurements

PerfectStride					
Size EU	Length (mm)	Joint Width	Joint Circ.	Instep Circ.	Toe Depth
36	242	97	250	238	27
37	249	99	255	243	28
38	256	100	259	247	28
39	262	102	264	252	28
40	269	103	268	256	29
41	275	105	273	261	29
42	282	106	277	265	30
43	289	108	282	270	30
44	295	109	286	274	30
45	302	111	291	279	31
46	309	112	295	283	31

# PerfectStride Free

Specialised Footwear for People with Diabetes and Peripheral Neuropathy

## Features

- Seamless joint area to minimize irritation.
- Optimised rocker angle (~20°) for improved foot movement efficiency non-slip rubber track.
- Includes two adjustment inlays to regulate interior volume.
- Carbon spring reinforcement over partial length of the foot plate.
- C8mm wider across the metatarsal heads, plus 2 x 5mm EVA removable footbeds.



## Measurements

PerfectStride					
Size EU	Length (mm)	Joint Width	Joint Circ.	Instep Circ.	Toe Depth
36	242	98	262	253	37
37	249	99	267	257	38
38	256	101	271	261	38
39	262	102	276	266	39
40	269	104	280	270	39
41	275	105	285	275	40
42	282	107	289	279	40
43	289	108	294	283	40
44	295	110	298	288	41
45	302	111	303	292	41
46	309	113	307	297	41

# XL Series For Healthy Walking

## PerfectStride Free XL Hook/Loop Series

Wide opening shoe.  
Carbon Spring System with gait control.



Black  
Code: MSSC009-\*Size

**Size Range:** 36-46 - Unisex

**Upper:** Artificial Leather

**Lining:** Textile

**Sole:** EVA/Rubber with Carbon Insole

**Footbed:** Removable, 2 x 5mm

## PerfectStride Free XL Hook/Loop Open Toe Series

Fully opening sandal.  
Carbon Spring System with gait control.



Black  
Code: MSSC010-\*Size

**Size Range:** 36-46 - Unisex

**Upper:** Artificial Leather

**Lining:** Textile

**Sole:** EVA/Rubber with Carbon Insole

**Footbed:** Removable, 2 x 5mm

# YDA Trainers

Available direct from UK stock in 5 days

YDA shoes stabilize and assist every step you take. The SPRING SYSTEM® inside the sole encourages a mechanical rolling motion and creates an upward thrust that improves the movement performance. YDA footwear stimulates foot and leg muscles, improves posture and the circulation, protects the foot from injury and reduces the most common foot pain. It is suitable for the management of the diabetic foot, whether or not at risk of developing foot ulcerations (Diabetes + Neuropathy and/or Vasculopathy) and Morton's Neuroma.

## YDA Black

Ladies

Code: PH-YDA-WBS-\*Size

Size Range: 35-41

Gents

Code: PH-YDA-MBS-\*Size

Size Range: 42-47



## YDA Blue

Ladies

Code: PH-YDA-WBL-\*Size

Size Range: 35-41

Gents

Code: PH-YDA-MBL-\*Size

Size Range: 42-47



## YDA Grey

Ladies

Code: PH-YDA-WGY-\*Size

Size Range: 35-41

Gents

Code: PH-YDA-MGY-\*Size

Size Range: 42-47



## YDA Cairn Powder Pink

Ladies

Code: PH-YDA-CROSA\*Size

Size Range: 36-41



## YDA Unisex

Ladies & Gents

Code: PH-YDA-CBlack\*Size

Size Range: 36-46



## YDA Measurement Guides

YDA Trainer Size Chart Ladies/Gents

Size	Length (mm)	Joint Width (mm)	Joint Circ. (mm)	Toe Depth (mm)
35	23.4	8.1	22.8	35
36	24	8.3	23.3	35
37	24.6	8.5	23.8	35
38	25.3	8.7	24.3	35
39	26	8.9	24.8	35
40	26.6	9.1	25.3	35
41	27.3	9.3	25.8	37
42	28	9.5	26.3	37
43	28.6	9.7	26.8	37
44	29.3	9.9	27.3	37
45	29.9	10.1	27.8	37
46	30.6	10.3	28.3	37
47	31.3	10.5	28.8	37

YDA MAC Footwear - Ladies

Size	Length (mm)	Joint Width (mm)	Joint Circ. (mm)	Toe Depth (mm)
36	24.5	8.6	24	37
37	25.2	8.8	24.5	37
38	25.8	8.9	25	37
39	26.5	9.1	25.5	37
40	27.2	9.2	26	37
41	27.8	9.4	26.5	37

# YDA MAC Footwear

Available direct from UK stock in 5 days

Protective footwear with propulsive semi-rigid sole to redistribute loads. MAC3 footwear comply with international guidelines for managing the at risk and high risk diabetic foot (diabetes + neuropathy and/or vasculopathy and/or deformity and/or pre-ulcerative signs).



## YDA MAC 3 Glitter

Ladies

**Code:** PH-MAC-Glitter-\*Size

**Size Range:** 36-41



## YDA MAC 3 Softy

Ladies

**Code:** PH-MAC-SOFTY-\*Size

**Size Range:** 36-41

# Fior & Gentz

## Stock Footwear





## Insole Set Available direct from UK stock

### Features

- Antimicrobial top cover
- Manufactured from EVA, providing different shore hardnesses:  
Grey = 46°    White = 23°    Antimicrobial = 28°

#### Insole Set Standard

Size	Product Code (Left Set)	Product Code (Right Set)
35	FO0053-35/L	FO0053-35/R
36	FO0053-36/L	FO0053-36/R
37	FO0053-37/L	FO0053-37/R
38	FO0053-38/L	FO0053-38/R
39	FO0053-39/L	FO0053-39/R
40	FO0053-40/L	FO0053-40/R
41	FO0053-41/L	FO0053-41/R
42	FO0053-42/L	FO0053-42/R
43	FO0053-43/L	FO0053-43/R
44	FO0053-44/L	FO0053-44/R
45	FO0053-45/L	FO0053-45/R
46	FO0053-46/L	FO0053-46/R
47	FO0053-47/L	FO0053-47/R
48	FO0053-48/L	FO0053-48/R
49	FO0053-49/L	FO0053-49/R
50	FO0053-50/L	FO0053-50/R

#### Insole Set Wide

Size	Product Code (Left Set)	Product Code (Right Set)
35	FO0054-35/L	FO0054-35/R
36	FO0054-36/L	FO0054-36/R
37	FO0054-37/L	FO0054-37/R
38	FO0054-38/L	FO0054-38/R
39	FO0054-39/L	FO0054-39/R
40	FO0054-40/L	FO0054-40/R
41	FO0054-41/L	FO0054-41/R
42	FO0054-42/L	FO0054-42/R
43	FO0054-43/L	FO0054-43/R
44	FO0054-44/L	FO0054-44/R
45	FO0054-45/L	FO0054-45/R
46	FO0054-46/L	FO0054-46/R
47	FO0054-47/L	FO0054-47/R
48	FO0054-48/L	FO0054-48/R
49	FO0054-49/L	FO0054-49/R
50	FO0054-50/L	FO0054-50/R



## EVA Insole Set

- Manufactured from EVA, providing different shore hardnesses:  
Black = 40    Grey = 25    White = 50

#### Insole Set Standard (Fits Flensburg, Heidelberg, Leipzig, Kassel and Nurnberg)

Size	Product Code
35	FO0033-35
36	FO0033-36
37	FO0033-37
38	FO0033-38
39	FO0033-39
40	FO0033-40
41	FO0033-41
42	FO0033-42
43	FO0033-43
44	FO0033-44
45	FO0033-45
46	FO0033-46
47	FO0033-47
48	FO0033-48
49	FO0033-49
50	FO0033-50

#### Insole Set Wide (Fits Erfurt, Flensburg, Dresden, Leipzig, Kassel, Nurnberg)

Size	Product Code
35	FO0034-35
36	FO0034-36
37	FO0034-37
38	FO0034-38
39	FO0034-39
40	FO0034-40
41	FO0034-41
42	FO0034-42
43	FO0034-43
44	FO0034-44
45	FO0034-45
46	FO0034-46
47	FO0034-47
48	FO0034-48
49	FO0034-49
50	FO0034-50



# Measurement Guide

## Standard

(Flensburg, Erfurt  
Heidelberg and Leipzig)

Inner Shoe Dimensions (mm)			
Size	A	B	C
35	95	66	227
36	96	67	233
37	98	68	239
38	100	69	246
39	101	70	251
40	102	71	257
41	105	73	266
42	107	74	272
43	108	75	280
44	109	75	288
45	111	76	292
46	112	77	299
47	113	78	304
48	114	79	310
49	116	80	315
50	118	81	320

## Standard

(Dresden, Kassel  
and Nürnberg)

Inner Shoe Dimensions (mm)			
Size	A	B	C
35	93	63	239
36	94	63	243
37	95	65	248
38	96	66	256
39	99	67	264
40	100	69	268
41	102	70	272
42	104	70	280
43	105	70	290
44	108	72	295
45	108	73	300
46	110	75	306
47	113	75	315
48	113	77	317

## Wide


(Flensburg and Erfurt)

Inner Shoe Dimensions (mm)			
Size	A	B	C
35	105	69	233
36	107	69	240
37	109	70	245
38	110	71	250
39	111	72	255
40	112	73	260
41	112	74	268
42	113	75	273
43	115	76	283
44	120	77	290
45	123	78	295
46	123	78	301
47	125	81	307
48	126	82	315
49	129	84	321
50	130	85	328

## Wide

(Dresden, Kassel  
and Nürnberg)


Inner Shoe Dimensions (mm)			
Size	A	B	C
35	103	67	243
36	104	68	247
37	105	70	254
38	108	70	259
39	109	71	264
40	109	71	270
41	110	72	276
42	111	73	283
43	113	74	288
44	115	75	295
45	117	76	301
46	118	78	308
47	120	80	313
48	122	82	318



## Flensberg Available direct from UK stock

### Features

- Asymmetric, spacious diabetic shoe
- Available in standard and wide fittings
- Dorsal zip allows donning and doffing without using the hook and loop fasteners.
- Integrated toe cap for additional stability and protection
- 30mm thick, rigid and extra strong rocker bottom sole with pitch prevents flexion of the metatarsophalangeal joints, ensuring the patient's foot is stabilised and immobilised
- Each shoe is supplied with three 5mm thick EVA insoles (STANDARD FO0053-XX and WIDE FO0054-XX). The top insole has an antimicrobial coating. The insoles can be used as an interim solution for pressure relief until being replaced with a diabetic foot orthotic.



#### Flensberg Standard

Size	Product Code (Left)	Product Code (Right)
35	TS5311-35/L	TS5311-35/R
36	TS5311-36/L	TS5311-36/R
37	TS5311-37/L	TS5311-37/R
38	TS5311-38/L	TS5311-38/R
39	TS5311-39/L	TS5311-39/R
40	TS5311-40/L	TS5311-40/R
41	TS5311-41/L	TS5311-41/R
42	TS5311-42/L	TS5311-42/R
43	TS5311-43/L	TS5311-43/R
44	TS5311-44/L	TS5311-44/R
45	TS5311-45/L	TS5311-45/R
46	TS5311-46/L	TS5311-46/R
47	TS5311-47/L	TS5311-47/R
48	TS5311-48/L	TS5311-48/R
49	TS5311-49/L	TS5311-49/R
50	TS5311-50/L	TS5311-50/R


#### Flensberg Wide

Size	Product Code (Left)	Product Code (Right)
35	TS5911-35/L	TS5911-35/R
36	TS5911-36/L	TS5911-36/R
37	TS5911-37/L	TS5911-37/R
38	TS5911-38/L	TS5911-38/R
39	TS5911-39/L	TS5911-39/R
40	TS5911-40/L	TS5911-40/R
41	TS5911-41/L	TS5911-41/R
42	TS5911-42/L	TS5911-42/R
43	TS5911-43/L	TS5911-43/R
44	TS5911-44/L	TS5911-44/R
45	TS5911-45/L	TS5911-45/R
46	TS5911-46/L	TS5911-46/R
47	TS5911-47/L	TS5911-47/R
48	TS5911-48/L	TS5911-48/R
49	TS5911-49/L	TS5911-49/R
50	TS5911-50/L	TS5911-50/R

# Dresden Available direct from UK stock

## Features

- Asymmetric, spacious diabetic shoe
- Available in standard and wide fittings
- The firm outer fabric can be separated from the soft inner fabric to eliminate possible pressure points
- 30mm thick, rigid and extra strong rocker bottom sole with pitch prevents flexion of the metatarsophalangeal joints, ensuring the patient's foot is stabilised and immobilised
- Supplied with a cutting template but without insoles. If insoles are required please see insole sets on page 22



### Dresden Standard

Size	Product Code (Left)	Product Code (Right)
35	TS5811-35/L	TS5811-35/R
36	TS5811-36/L	TS5811-36/R
37	TS5811-37/L	TS5811-37/R
38	TS5811-38/L	TS5811-38/R
39	TS5811-39/L	TS5811-39/R
40	TS5811-40/L	TS5811-40/R
41	TS5811-41/L	TS5811-41/R
42	TS5811-42/L	TS5811-42/R
43	TS5811-43/L	TS5811-43/R
44	TS5811-44/L	TS5811-44/R
45	TS5811-45/L	TS5811-45/R
46	TS5811-46/L	TS5811-46/R
47	TS5811-47/L	TS5811-47/R
48	TS5811-48/L	TS5811-48/R


### Dresden Wide

Size	Product Code (Left)	Product Code (Right)
35	TS5411-35/L	TS5411-35/R
36	TS5411-36/L	TS5411-36/R
37	TS5411-37/L	TS5411-37/R
38	TS5411-38/L	TS5411-38/R
39	TS5411-39/L	TS5411-39/R
40	TS5411-40/L	TS5411-40/R
41	TS5411-41/L	TS5411-41/R
42	TS5411-42/L	TS5411-42/R
43	TS5411-43/L	TS5411-43/R
44	TS5411-44/L	TS5411-44/R
45	TS5411-45/L	TS5411-45/R
46	TS5411-46/L	TS5411-46/R
47	TS5411-47/L	TS5411-47/R
48	TS5411-48/L	TS5411-48/R

## Nürnberg Available direct from UK stock

### Features

- Asymmetric, spacious diabetic shoe
- Available in standard and wide fittings
- The firm outer fabric can be separated from the soft inner fabric to avoid possible pressure points
- 30mm thick, rigid and extra strong rocker bottom sole with pitch prevents flexion of the metatarsophalangeal joints, ensuring the patient's foot is stabilised and immobilised
- Supplied with a cutting template but without insoles. If insoles are required please see insole sets on page 22.



### Nurnberg Standard

Size	Product Code (Left)	Product Code (Right)
35	TS4311-35/L	TS4311-35/R
36	TS4311-36/L	TS4311-36/R
37	TS4311-37/L	TS4311-37/R
38	TS4311-38/L	TS4311-38/R
39	TS4311-39/L	TS4311-39/R
40	TS4311-40/L	TS4311-40/R
41	TS4311-41/L	TS4311-41/R
42	TS4311-42/L	TS4311-42/R
43	TS4311-43/L	TS4311-43/R
44	TS4311-44/L	TS4311-44/R
45	TS4311-45/L	TS4311-45/R
46	TS4311-46/L	TS4311-46/R
47	TS4311-47/L	TS4311-47/R
48	TS4311-48/L	TS4311-48/R

### Nurnberg Wide


Size	Product Code (Left)	Product Code (Right)
35	TS4411-35/L	TS4411-35/R
36	TS4411-36/L	TS4411-36/R
37	TS4411-37/L	TS4411-37/R
38	TS4411-38/L	TS4411-38/R
39	TS4411-39/L	TS4411-39/R
40	TS4411-40/L	TS4411-40/R
41	TS4411-41/L	TS4411-41/R
42	TS4411-42/L	TS4411-42/R
43	TS4411-43/L	TS4411-43/R
44	TS4411-44/L	TS4411-44/R
45	TS4411-45/L	TS4411-45/R
46	TS4411-46/L	TS4411-46/R
47	TS4411-47/L	TS4411-47/R
48	TS4411-48/L	TS4411-48/R



## Kassel Available direct from UK stock

### Features

- Asymmetric, spacious diabetic shoe
- Available in standard and wide fittings
- The firm outer fabric can be separated from the soft inner fabric to avoid possible pressure points
- Protective forefoot cap can be attached to the toe area to keep the foot dry and warm
- 30mm thick, rigid and extra strong rocker bottom sole with pitch prevents flexion of the metatarsophalangeal joints, ensuring the patient's foot is stabilised and immobilised
- Supplied with a cutting template but without insoles. If insoles are required please see insole sets on page 22.



### Kassel Standard

Size	Product Code (Left)	Product Code (Right)
35	TS5421-35/L	TS5421-35/R
36	TS5421-36/L	TS5421-36/R
37	TS5421-37/L	TS5421-37/R
38	TS5421-38/L	TS5421-38/R
39	TS5421-39/L	TS5421-39/R
40	TS5421-40/L	TS5421-40/R
41	TS5421-41/L	TS5421-41/R
42	TS5421-42/L	TS5421-42/R
43	TS5421-43/L	TS5421-43/R
44	TS5421-44/L	TS5421-44/R
45	TS5421-45/L	TS5421-45/R
46	TS5421-46/L	TS5421-46/R
47	TS5421-47/L	TS5421-47/R
48	TS5421-48/L	TS5421-48/R

### Kassel Wide

Size	Product Code (Left)	Product Code (Right)
35	TS5321-35/L	TS5321-35/R
36	TS5321-36/L	TS5321-36/R
37	TS5321-37/L	TS5321-37/R
38	TS5321-38/L	TS5321-38/R
39	TS5321-39/L	TS5321-39/R
40	TS5321-40/L	TS5321-40/R
41	TS5321-41/L	TS5321-41/R
42	TS5321-42/L	TS5321-42/R
43	TS5321-43/L	TS5321-43/R
44	TS5321-44/L	TS5321-44/R
45	TS5321-45/L	TS5321-45/R
46	TS5321-46/L	TS5321-46/R
47	TS5321-47/L	TS5321-47/R
48	TS5321-48/L	TS5321-48/R



## Protective Forefoot Cap Available direct from UK stock

### Features

- Easily attached to Kassel Standard and Wide
- Keeps the foot dry and warm

#### Protective Forefoot Cap

Size	Product Code	Size	Product Code
35-36	FO9520-XXS	43-44	FO9520-L
37-38	FO9520-XS	45-46	FO9520-XL
39-40	FO9520-S	47-48	FO9520-XXL
41-42	FO9520-M	49-50	FO9520-XXXL



## Bremen Available direct from UK stock

Particularly suited for patients with wounds at the foot/forefoot area on the sole of the foot.


### Features

- Spacious post-op shoe with a generous instep fastener for easier donning and doffing
- Integrated toe cap for added stability and protection
- Universal, fits left or right
- 30mm thick, rigid and extra strong rocker bottom sole with pitch
- Supplied with an antimicrobial coated EVA insole (FO0064-XX)



#### Bremen

Size	Inner Shoe Dimensions		Product Code
	A	C	
XXS	89	209	TS1151-XXS
XS	92	230	TS1151-XS
S	100	236	TS1151-S
M	102	254	TS1151-M
L	106	268	TS1151-L
XL	109	285	TS1151-XL
XXL	110	301	TS1151-XXL




## Erfurt Available direct from UK stock

Particularly suited for patients with wounds at the foot/forefoot area on the sole of the foot, swellings or joint pain. Also suitable for those who need a closed shoe with a fastener system offering wider, more generous opening and a tall shaft.

### Features

- Asymmetric, spacious post-op shoe with a generous instep fastener
- Available in standard and wide fittings
- Dorsal zip for easier donning and doffing
- Integrated toe cap for added stability and protection
- 30mm thick, rigid and extra strong rocker bottom sole with pitch
- Supplied with three 5mm thick antimicrobial coated EVA insoles (Standard FO0053-XX and Wide FO0054-XX)



### Erfurt Standard

Size	Product Code (Left)	Product Code (Right)
35	TS5310-35/L	TS5310-35/R
36	TS5310-36/L	TS5310-36/R
37	TS5310-37/L	TS5310-37/R
38	TS5310-38/L	TS5310-38/R
39	TS5310-39/L	TS5310-39/R
40	TS5310-40/L	TS5310-40/R
41	TS5310-41/L	TS5310-41/R
42	TS5310-42/L	TS5310-42/R
43	TS5310-43/L	TS5310-43/R
44	TS5310-44/L	TS5310-44/R
45	TS5310-45/L	TS5310-45/R
46	TS5310-46/L	TS5310-46/R
47	TS5310-47/L	TS5310-47/R
48	TS5310-48/L	TS5310-48/R
49	TS5310-49/L	TS5310-49/R
50	TS5310-50/L	TS5310-50/R


### Erfurt Wide

Size	Product Code (Left)	Product Code (Right)
35	TS5410-35/L	TS5410-35/R
36	TS5410-36/L	TS5410-36/R
37	TS5410-37/L	TS5410-37/R
38	TS5410-38/L	TS5410-38/R
39	TS5410-39/L	TS5410-39/R
40	TS5410-40/L	TS5410-40/R
41	TS5410-41/L	TS5410-41/R
42	TS5410-42/L	TS5410-42/R
43	TS5410-43/L	TS5410-43/R
44	TS5410-44/L	TS5410-44/R
45	TS5410-45/L	TS5410-45/R
46	TS5410-46/L	TS5410-46/R
47	TS5410-47/L	TS5410-47/R
48	TS5410-48/L	TS5410-48/R
49	TS5410-49/L	TS5410-49/R
50	TS5410-50/L	TS5410-50/R

# Leipzig Available direct from UK stock

## Features

- Asymmetric, spacious diabetic shoe
- Hook and loop fasteners, and dorsal zip, facilitate generously wide opening for easier donning and doffing
- Integrated toe cap for more stability and protection in the toe area
- 30mm thick, flexible sole with pitch supports a natural heel to toe movement while walking
- Supplied with three 5mm thick antimicrobial coated EVA insoles (FO0053-XX)



## Leipzig

Size	Product Code (Left)	Product Code (Right)
35	TS5512-35/L	TS5512-35/R
36	TS5512-36/L	TS5512-36/R
37	TS5512-37/L	TS5512-37/R
38	TS5512-38/L	TS5512-38/R
39	TS5512-39/L	TS5512-39/R
40	TS5512-40/L	TS5512-40/R
41	TS5512-41/L	TS5512-41/R
42	TS5512-42/L	TS5512-42/R
43	TS5512-43/L	TS5512-43/R
44	TS5512-44/L	TS5512-44/R
45	TS5512-45/L	TS5512-45/R
46	TS5512-46/L	TS5512-46/R
47	TS5512-47/L	TS5512-47/R
48	TS5512-48/L	TS5512-48/R
49	TS5512-49/L	TS5512-49/R
50	TS5512-50/L	TS5512-50/R



## Heidelberg Available direct from UK stock


### Features

- Asymmetric, spacious diabetic shoe
- Provides a variable volume
- Partially flexible upper fabric is elastic and adapts to different demands on space
- 30mm thick, flexible sole with pitch supports a natural heel to toe movement while walking
- Supplied with three 5mm thick antimicrobial coated EVA insoles (FO0053-XX)



### Heidelberg

Size	Product Code (Left)	Product Code (Right)
35	TS4552-35/L	TS4552-35/R
36	TS4552-36/L	TS4552-36/R
37	TS4552-37/L	TS4552-37/R
38	TS4552-38/L	TS4552-38/R
39	TS4552-39/L	TS4552-39/R
40	TS4552-40/L	TS4552-40/R
41	TS4552-41/L	TS4552-41/R
42	TS4552-42/L	TS4552-42/R
43	TS4552-43/L	TS4552-43/R
44	TS4552-44/L	TS4552-44/R
45	TS4552-45/L	TS4552-45/R
46	TS4552-46/L	TS4552-46/R
47	TS4552-47/L	TS4552-47/R
48	TS4552-48/L	TS4552-48/R
49	TS4552-49/L	TS4552-49/R
50	TS4552-50/L	TS4552-50/R





HEIDELBERG  
UNIVERSITÄT



**Camino**  
Modular  
Collection





# Modular Measurement Guide

## Measuring Guide

All measurements should be taken with the patient sitting down.

### 1. Foot length

Use the Camino foot measuring stick to measure the length of the foot.

### 2. Joint width

The joint width is taken obliquely across the foot from the first to the fifth metatarsal joints at the widest point.

### 3-5. Joint, instep and ankle circumferences


The tape should be pulled tightly enough to compress the tissues a little and then released sufficiently to allow the normal elasticity of the tissue to bring the tape back to the correct measurement tension. The leg circumference should be taken at the correct height. See table on pages 74 and 75.

### 6. Heel width

The heel width measurement should be taken at the widest part of the heel.


### 7. Toe depth

Use the measuring stick to measure the toe depth.



1. Foot length
2. Joint width
3. Joint circumference
4. Instep circumference
5. Leg circumference (boots only)
6. Heel width
7. Toe depth

## Stock/Modular Footwear Measurement Process




# Modular Last Additions

All additions must be in multiples of 3mm only, e.g 3mm, 6mm, 9mm,12mm and 15mm, except when an increase in length is requested. 4mm for half a size and 8mm for a full size. A round toe last is used as default. A square toe last has a straightened medial border providing more room for the hallux. Please note that stock lasts can only have additions added to them, they cannot be reduced in dimension.

## Option 1

### Standard Shoe

Any standard last, without modifications to last or upper, with a choice of shoes and boots.



## Further Options

Special pattern changes can be requested to make modifications to the upper to improve the overall fit i.e lower under ankle section or lower opening .


Special pattern change	CA0103
Square toe option	CA0105
Commando sole unit	CA0106

## Option 2


### CA0100

Any one single modification per last as illustrated below.


**Not to exceed 6mm per last.**




A Increased instep depth




B Increased joint circumference




C Increased toe depth




D Increased width on lateral or medial side




E Increased heel width



F Increased forefoot depth



G Increased through depth



H Increased toe length


Footwear ordered that exceeds the specifications detailed in Option 3 will be treated as Made to Measure footwear and charged accordingly.

## Option 3

### CA0102

Several modifications in multiples of 3mm, 6mm, 9mm or 12mm. Allows you to choose any combination of modifications as illustrated in below.

Not to exceed 15mm per last.



Example:	D Increased width on lateral or medial side	= 6mm
	E Increased heel width	= 3mm
	G Increased through depth	= 6mm
	Total	= 15mm

Footwear can be supplied with up to 3 x 3mm shock absorbing insoles maximum. Please state the number required on your order form.



# Materials

Camino footwear is manufactured from the finest quality materials available.

**Calf leather** has a fine smooth grain. It is soft but sturdy and is suitable for most default footwear.

**Diabetic soft leather** is a thinned calf, making it softer and more supple than standard calf. It is not suitable for use with stiffeners and for patients requiring strong sturdy footwear.

**Distressed leather** is a strong, thicker leather than calf and suitable for the manufacture of more durable footwear. The distressed leather is smooth in appearance.

**Siciliana leather** is a thick leather but is soft and supple due to the tanning process producing a grained finish.

**Strong leather** is the strongest and most durable leather in our range and is ideal for moisture resistance.

**Nubuck leather** has a softer appearance due to the buffed surface but is a robust material with good longevity. It cannot be polished but can be cleaned like suede.

**Stretch leather** can be used for aprons and side panels on certain styles. It is made from natural leather. Stretch leather cannot be used to manufacture a full vamp or full shoe.

## Leather thickness.

Ladies and Gents shoes/boots are manufactured in calf leather, which is 1.4 - 1.6mm thickness. The thickness of Strong Leather, Distressed Leather Val, Nubuck and Siciliana leathers is 1.8 - 2mm. Soft leather is 1.2 - 1.4mm thick.

Please note that the thickness is not the only value that determines the strength of leather; leathers of the same thickness can offer very different properties due to the different tanning processes which break the fibers changing the malleability.

## Synthetic upper materials.

Neoprene can be used to manufacture the apron on some styles to provide a soft stretchy upper. Ideal for painful deformed toes.

**Vegan synthetic leathers** meet the standards as set out and certified by Inescop.

**Corduroy Uppers** are only available on house shoes and will have a fabric lining which will be matched to the colour of the corduroy. Corduroy is a durable cloth textile with a raised parallel cord effect.

**Toe Puffs** stiffen and retain the toe box shape of the footwear.

Three types of toe puffs are available for Made to Measure and Modular footwear:

1. Standard toe puffs are the default and are used in most footwear types, stiffening a larger part of the toe area at the sides and over the dorsum.
2. Winged toe puffs do not extend over the top of the second, third and fourth toes and are ideal for sensitive feet and clawed toes.
3. Rimmed toe puffs hold the front of the shoe up but do not extend back over the dorsum of the toes or along the sides of the toe box. They should only be used when feet are very sensitive or deformities are irritated by an alternative toe puff.

**Heel stiffeners** help to maintain the shape and strength at the rear part of the footwear.

A standard heel stiffener covers the medial and lateral sides of the heel and is joined at the back.

Extended heel stiffeners extend along the medial, lateral or both sides of the foot to behind the joint area

High stiffeners provide additional strength on the medial, lateral or both sides of a boot. As default the height will be to below the collar. Other heights can be specified.

Toe puffs and heel stiffeners are made from a synthetic material as default. The material can be doubled for increased strength.

Custom heel stiffeners and toe puffs can be specified in Made to Measure footwear made from different materials.



# Leather Options

The Camino Collection is manufactured using only the finest quality leather and materials. Please note that actual colours may vary slightly from the illustrated colour tiles.

## Calf Leather



Black



Navy Blue




Dark Brown



Taupe




Cream



White




Caribe Blue



Mid Brown



Dark Burgundy



Light Burgundy




Red



Purple



Olive



Polisander




Grey




Pink

## Diabetic Softee Leather



Black



Navy Blue



Dark Brown




Mid Brown




Taupe


## Nubuck




Black



Navy



Dark Brown




Cream




Pink

## Strong Leather



Black



Navy

## Distressed Leather Val




Leather - Val  
Colour - Nut



Leather - Val  
Colour - Saddale



Leather - Val  
Colour - Soan



Leather - Val  
Colour - Linen



Leather - Val  
Colour - Black



Leather - Val  
Colour - Espanol



Leather - Val  
Colour - Gorbi

## Siciliana Leather



Black



Brown



Blue



White



Silver

## Vegan Synthetic Leather



Black



Navy Blue



Pink



White



Taupe



Dark Brown




Mid Brown




Grey

# Panel Options for Specified Styles


## Stretch Leather




Black Pielax




Navy Pielax




Dark Brown Pielax




Black Ante-Tejus




Ante Cobra Black




Halley Black




Anteraiz



Ante Cobra Marron




Polisandro Ante-Tejus



Vanilla Ante -Tejus




Ante Cobra Vanilla



Viventi Black-Grey


## Neoprene




Black



Navy



Dark Brown




Taupe



Burgundy

## Cordura



Black



Navy




Silver



Taupe

## Linings




Plush Black



Plush Brown



Plush Beige




Leather Brown



Leather Grey




Black Sheepskin  
Synthetic



White Sheepskin  
Synthetic



Leather Beige



Leather White

## Corduroy

Only available for house shoes. The lining colour will match the upper.



Black



Navy



Brown



Grey



Gold

# Sole Unit and Punch Pattern Options

## Modular and Made to Measure Sole Unit Options



### OPTION 1

#### Separate Sole and Heel

Pyramid EVA midsole with textured rubber top piece.  
Heel height 22mm as default at mid-heel with a 2mm +/- tolerance, size dependent.  
Fore part height 8mm.



### OPTION 2

#### Semi Wedge

Pyramid EVA midsole with lightweight rubber top piece and semi wedge heel.  
Heel height 22mm as default at mid-heel with a 2mm +/- tolerance, size dependent.  
Fore part height 6mm.



### OPTION 3

#### Through Sole Unit

Pyramid EVA midsole with textured rubber top piece.  
Heel height 22mm as default at mid-heel with a 2mm +/- tolerance, size dependent.  
Fore part height 5mm.



### OPTION 4

#### Commando Sole Unit

EVA/rubber midsole.  
Heel height 25mm as default at mid-heel with a 2mm +/- tolerance, size dependent.  
Fore part height 10mm.  
Available in:  
**A)** Light Weight  
**B)** Heavy Weight  
**C)** Heavy Weight Through Sole

## Punch Pattern Options

Only use punch patterns on 2nd issue footwear onwards, when the fitting is 100% confirmed. This reduces the chance of it being pulled off centre after adjustments are made. Punch pattern cannot be added at fitting.



Ladies Pattern 1



Ladies Pattern 2



Ladies Pattern 3



Ladies Pattern 4



Ladies Pattern 5



Gents Pattern 1



Gents Pattern 2



Gents Pattern 3



Gents Pattern 4



**OPTION 5**

**Semi Wedge Textured**

EVA midsole with textured rubber top piece and semi wedge heel.

Heel height 22mm as default at mid-heel with a 2mm +/- tolerance, size dependent.

Fore part height 5mm.



**OPTION 6**

**Separate Sole and Heel Textured**

EVA midsole with textured rubber top piece.

Heel height 22mm as default at mid-heel with a 2mm +/- tolerance, size dependent.

Fore part height 5mm.



**OPTION 7**


**Trainer Sole Unit With Negative Heel**

Light weight flexible textured construction with good grip.

Available in Black/White.

Heel height 22mm as default at mid-heel with a 2mm +/- tolerance, size dependent.

Fore part height 5mm.



**OPTION 8**

**Non Standard Sole Unit**

This option can be manufactured as a Separate Sole and Heel, Semi Wedge or Through Sole Unit.

This can be manufactured from a variety of materials. Please specify on your order.

## Mock Welt




Gents Welt



Ladies Welt

## Toe Shapes



# Camino **Unisex** **House** Collection



# Camino **Unisex House** Collection

## Tranquil (Product Code CA0079V)

Tranquil is a house shoe available with a low opening with a lay over hook and loop fastening.

- Tips:**
1. The soft corduroy fabric makes this shoe ideal for indoor use.
  2. This low opening shoe is ideal for feet with variable oedema and when easy access is required.
  3. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in lay on hook and loop

**Materials:** Corduroy upper only

**Sole Units:** All



Shown in  
Grey

## Relax (Product Code CA0080V)

Relax is a house shoe available with a standard opening with either a lay over or reverse pull hook and loop fastening.

- Tips:**
1. The soft corduroy fabric makes this shoe ideal for indoor use.
  2. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull or lay on hook and loop

**Materials:** Corduroy upper only

**Sole Units:** All



Shown in  
Navy

## Calm (Product Code CA0081V)

Calm is a house shoe available with a low opening with either a leather lay over or reverse pull hook and loop fastening.

- Tips:**
1. The soft corduroy fabric makes this shoe ideal for indoor use.
  2. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull or lay on hook and loop

**Materials:** Corduroy upper only

**Sole Units:** All



Shown in  
Gold



# Ladies Camino Modular Collection

All styles in this section are manufactured as illustrated. If there are to be any changes to the upper design or quantity of straps required then a special pattern change will need to be requested at the point of ordering.



# Ladies Camino Modular Collection

## Coral (Product Code Lace CA0007L, Hook & Loop CA0007V)

Coral is a 3 eyelet Derby plain front shoe.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. The standard opening means this style is not recommended for oedema or when easy access is required.
  3. No seams at the joints make it suitable for feet with painful or sensitive 1st or 5th joints.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in  
Taupe  
and Dark  
Burgundy

## Nottingham (Product Code Lace CA0035L, Hook & Loop CA0035V)

Nottingham is a 4 eyelet Derby plain front shoe.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. This low opening shoe is ideal for feet with variable oedema and when easy access is required.
  3. Punch patterns cannot be added on the vamp.
  4. No seams at the joint make it suitable for feet with painful or sensitive 1st or 5th joints.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown  
in Light  
Burgundy  
and Navy

## Moonstone (Product Code Lace CA0064L, Hook & Loop CA0064V)

Moonstone is a 4 eyelet Derby plain front shoe.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. This low opening shoe is ideal for feet with variable oedema and when easy access is required. The opening extends 15-20mm lower than Nottingham.
  3. No seams at the joint make it suitable for feet with painful or sensitive 1st or 5th joints.
  4. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in  
Black and  
Pink

## Larimar

(Product Code Lace CA0014L, Hook & Loop CA0014V)

Larimar is a 5 eyelet Derby plain front shoe.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. This low opening shoe is ideal for feet with variable oedema and when easy access is required.
  3. No seams at the joints make it suitable for feet with painful or sensitive 1st or 5th joints.
  4. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in  
Taupe and  
Navy

## Sunstone

(Product Code Lace CA0078L, Hook & Loop CA0078V)

Sunstone is a 5 eyelet Derby plain front shoe with contrasting trim quilted effect quarters.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. This low opening shoe is ideal for feet with variable oedema and when easy access is required.
  3. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in  
Navy Blue  
and Taupe  
with Cream  
detail

## Jasper

(Product Code Lace CA0013L, Hook & Loop CA0013V)

Jasper is a 5 eyelet Oxford type shoe with decorative apron stitch.

- Tips:**
1. The U opening front makes access and adjustability more difficult so this style is not recommended for oedema or when easy access is required.
  2. No seams at the joints make it suitable for feet with painful or sensitive 1st or 5th joints.
  3. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown  
in Dark  
Brown with  
Medium  
Brown  
detail and  
Taupe

# Ladies Camino Modular Collection

## Jade (Product Code Lace CA0012L, Hook & Loop CA0012V)

Jade is a 5 eyelet Oxford type shoe with decorative upper stitch.

- Tips:**
1. The standard opening means this style is not recommended for oedema or when easy access is required.
  2. Punch patterns cannot be added on the vamp.
  3. No seams at the joint make it suitable for feet with painful or sensitive 1st or 5th joints.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in  
Red and  
Taupe  
Nubuck

## Ruby (Product Code Lace CA0016L, Hook & Loop CA0016V)

Ruby is a 3 eyelet Blucher type shoe with stitched apron design.

- Tips:**
1. The standard opening means this style is not recommended for oedema or when easy access is required.
  2. No seams at the joints make it suitable for feet with painful or sensitive 1st or 5th joints.
  3. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown  
in Cream  
and Red

## Pearl (Product Code Lace CA0015L, Hook & Loop CA0015V)

Pearl is a 3 eyelet Derby type shoe with stitched apron design and fold pattern.

- Tips:**
1. The standard opening means this style is not recommended for oedema or when easy access is required.
  2. This style is unsuitable for sensitive or severely clawed toes due to the seams on the apron.
  3. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown  
in Navy  
Nubuck  
and Cream

## Suffolk

(Product Code Lace CA0039L, Hook & Loop CA0039V)

Suffolk is a 5 eyelet Oxford type shoe with decorative apron stitch.

- Tips:**
1. The stretch side panels provide extra ease for painful toes.
  2. The U opening front makes access and adjustability more difficult so this style is not recommended for oedema or when easy access is required.
  3. No seams at the joints make it suitable for feet with painful or sensitive 1st or 5th joints.
  4. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with neoprene  
Side panels are only available in stretch leather

**Sole Units:** All



Shown  
in Navy  
and Dark  
Brown

## Amethyst

(Product Code Hook & Loop CA0003V)

Amethyst is a T-bar shoe.

- Tips:**
1. A single strap secures the foot inside the shoe so it is not suitable for pronated or pes planus feet.
  2. This style is not suitable for oedema that may swell through the openings.
  3. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  4. The opening makes access and adjustability more difficult so this style is not recommended for oedema or when easy access is required.

**Fastenings:** Available in reverse pull hook and loop

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown  
in Cream  
Nubuck  
and Navy

## Shropshire

(Product Code Code Hook & Loop CA0043V)

Shropshire is a Mary Jane style shoe.

- Tips:**
1. A single strap secures the foot inside the shoe so it is not suitable for pronated or pes planus feet.
  2. This style is not suitable for oedema that may swell through the openings.
  3. Minimal seams on vamp makes the style ideal for low risk diabetic patients or those with sensitive feet.
  4. The opening allows easy access.
  5. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown  
in Olive  
and Dark  
Burgundy

# Ladies Camino Modular Collection

## Slate (Product Code Hook & Loop CA0065V)

Slate is a Mary Jane style shoe.

- Tips:**
1. A single strap secures the foot inside the shoe so it is not suitable for pronated or pes planus feet.
  2. This style is not suitable for severe oedema that may swell through the opening.
  3. The opening allows access easy.
  4. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop

**Materials:** Only available in colour options shown. Apron and edging is a stretched fabric.

**Sole Units:** All



Shown in Grey and Black with patterned stretched vamps

## Essex (Product Code Lace CA0038L, Hook & Loop CA0038V)

Essex is a 3 eyelet Derby type shoe with a stretch leather apron.

- Tips:**
1. The soft stretch front is ideal for sensitive feet or those with toe deformities.
  2. This low opening shoe is ideal for feet with variable oedema and when easy access is required.
  3. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Materials stretch leather apron and calf rearfoot and toe bumper

**Sole Units:** All



Shown in Navy and Dark Brown with stretched vamps

## Cheshire (Product Code Hook & Loop CA0041V)

Cheshire is a T-bar shoe with stretch leather side panels.

- Tips:**
1. A single strap secures the foot inside the shoe so it is not suitable for pronated or pes planus feet.
  2. This style is not suitable for oedema that may swell through the openings.
  3. The stretch side panels provide extra ease for painful toes.
  4. The opening makes access and adjustability more difficult so this style is not recommended for oedema or when easy access is required.
  5. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop

**Materials:** Stretched leather vamp and calf rearfoot

**Sole Units:** All



Shown in Black and Dark Brown

## Rutland

(Product Code Lace CA0042L, Hook & Loop CA0042V)

Rutland is a Monk bar shoe with a stretch leather apron.

- Tips:**
1. It offers easy access and the elastic inserts on the straps allow some adjustability for instep oedema.
  2. The soft stretch front is ideal for sensitive feet or those with toe deformities.
  3. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Materials stretch leather apron and calf rearfoot and toe bumper

**Sole Units:** All



Shown in Black and Light Burgundy

## Carnelian

(Product Code Lace CA0005L, Hook & Loop CA0005V)

Carnelian is a 3 eyelet Blucher type plain front shoe.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. The standard opening means this style is not recommended for oedema or when easy access is required.
  3. No seams at the joints make it suitable for feet with painful or sensitive 1st or 5th joints.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in Dark Brown and Caribe Blue

## Stafford

(Product Code Lace CA0036L, Hook & Loop CA0036V)

Stafford is a 3 eyelet Derby type shoe with a neoprene vamp.

- Tips:**
1. The soft stretch front is ideal for sensitive feet or those with toe deformities.
  2. This shoe is ideal for feet with variable oedema and when easy access is required.
  3. No seams at the joints make it suitable for feet with painful or sensitive 1st or 5th joints.
  4. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Neoprene vamp and calf rearfoot

**Sole Units:** All



Shown in Black and Dark Brown

# Ladies Camino Modular Collection

## Hertford (Product Code Lace CA0037L, Hook & Loop CA0037V)

Hertford is a 3 eyelet Derby type shoe with a stretch leather apron.

- Tips:** 1. *The soft stretch front is ideal for sensitive feet or those with toe deformities, especially bunions.*  
2. *Punch patterns cannot be added on the vamp.*

**Fastenings:** Hook and loop, and lace  
**Materials:** Materials stretch leather apron and calf rearfoot and toe bumper  
**Sole Units:** All



Shown in  
Taupe with  
Ante Cobra  
Vanilla  
stretched  
vamp and  
Black with  
Pielax  
stretched  
leather vamp

## Agate (Product Code Lace CA0002L, Hook & Loop CA0002V)

Agate is a 3 eyelet Blucher type shoe with decorative apron stitch which creates a slimming appearance.

- Tips:** 1. *No seams at the joints make it suitable for feet with painful or sensitive 1st or 5th joints.*  
2. *Punch patterns cannot be added on the vamp.*

**Fastenings:** Available in reverse pull hook and loop, and lace  
**Materials:** Not available with stretch leather or neoprene  
**Sole Units:** All



Shown in  
Red and  
Olive

## Titanite (Product Code Lace CA0066L, Hook & Loop CA0066V)

Titanite is a 5 eyelet Derby trainer with contrasting colour heel insert and punch pattern detail on the sides.

- Tips:** 1. *This low opening trainer is ideal for feet with variable oedema and when easy access is required.*
2. *Punch patterns cannot be added on the vamp.*

**Fastenings:** Hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** Trainer



Hook and loop shown in White with Rose Nubuck detail. Lace shown in White with Black detail

## Tigereye (Product Code Lace CA0017L, Hook & Loop CA0017V)

Tigereye is a 5 eyelet Oxford type trainer with contrasting colour stripes.

- Tips:** 1. *The U opening front makes access and adjustability more difficult so this style is not recommended for oedema or when easy access is required.*
2. *Punch patterns cannot be added on the vamp.*

**Fastenings:** Hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** Trainer



Shown in Black with White stripe and White with Black stripe

## Silver (Product Code Lace CA0067L, Hook & Loop CA0067V)

Silver is a 5 eyelet trainer with contrasting colour inserts

- Tips:** 1. *The U opening front makes access and adjustability more difficult so this style is not recommended for oedema or when easy access is required.*
2. *Punch patterns cannot be added on the vamp.*

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** Trainer



Shown in Black and Grey, and White Calf with Silver Siciliana insert

# Ladies Camino Modular Collection

## Malachite (Product Code Lace CA0082L, Hook & Loop CA0082V)

Malachite is a 6 eyelet Derby trainer with contrasting colour stripes.

- Tips:** 1. This low opening trainer is ideal for feet with variable oedema and when easy access is required.
2. Punch patterns cannot be added on the vamp.



- Fastenings:** Available in reverse pull hook and loop, and lace
- Materials:** Not available with stretch leather or neoprene
- Sole Units:** Trainer

Shown in Black and Grey, and Navy and Red

## Northumberland (Product Code Lace CA0047L, Hook & Loop CA0047V)

Northumberland is a 6 eyelet trainer boot with contrasting colour stripes.

- Tips:** 1. The U opening front makes access and adjustability more difficult so this style is not recommended for oedema or when easy access is required.
2. Punch patterns cannot be added on the vamp.



- Fastenings:** Available in reverse pull hook and loop, and lace
- Materials:** Not available with stretch leather or neoprene
- Sole Units:** Trainer

Shown in Black with White stripe and White with Navy stripe

## Berkshire (Product Code Lace CA0048L, Hook & Loop CA0048V)

Berkshire is a 7 eyelet Derby trainer boot with contrasting colour ankle/heel insert.

- Tips:** 1. This low opening shoe is ideal for feet with variable oedema and when easy access is required.
2. Punch patterns cannot be added on the vamp.



- Fastenings:** Available in reverse pull hook and loop, and lace
- Materials:** Not available with stretch leather or neoprene
- Sole Units:** Trainer

Shown in White with Black stripe and Navy with White stripe. Also available without stripes upon request

## Bronze

(Product Code Lace CA0068L, Hook & Loop CA0068V)

Bronze is a 6 eyelet boot with contrasting colour ankle/heel insert. Ideal for diabetics.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. This low opening boot is ideal for feet with variable oedema and when easy access is required.
  3. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in  
Taupe and  
Polisander

## Gold

(Product Code Lace CA0069L, Hook & Loop CA0069V)

Gold is a low opening Derby boot with quilted effect quarters.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. This low opening boot is ideal for feet with variable oedema and when easy access is required.
  3. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in  
Polisander  
and Olive  
Green

## Zircon

(Product Code Lace CA0020L, Hook & Loop CA0020V)

Zircon is a 5 eyelet Derby plain front boot.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. No seams at the joints make it suitable for feet with painful or sensitive 1st or 5th joints.
  3. This mid opening boot is suitable for feet with mild oedema.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in  
Black and  
Taupe

# Ladies Camino Modular Collection

## Sussex (Product Code Lace CA0050L, Hook & Loop CA0050V)

Sussex is a 5 eyelet Derby boot.

- Tips:**
1. The soft stretch front is ideal for sensitive feet or those with toe deformities.
  2. This low opening shoe is ideal for feet with variable oedema and when easy access is required.
  3. No seams at the joints make it suitable for feet with painful or sensitive 1st or 5th joints.
  4. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Neoprene vamp and calf rearfoot

**Sole Units:** All



Shown in Black and Dark Brown

## Topaz (Product Code Lace CA0070L, Hook & Loop CA0070V)

Topaz is a walking boot style with a combination of eyelets and lace loops and contrasting colour inserts.

- Tips:**
1. The standard opening means this style is not recommended for oedema or when easy access is required.
  2. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in Black with Siciliana Black Inserts and Navy with Taupe Inserts

## Opal (Product Code Lace CA0018L, Hook & Loop CA0018V)

Opal is a 5 eyelet Derby plain front boot.\*

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. The standard opening means this style is not recommended for oedema or when easy access is required.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All

\*Also available as a bootie, please state when ordering.



Shown in Black and Dark Burgundy





A man with a beard and a knit hat is smiling next to a beagle dog in a field. The man is wearing a dark jacket with a white fur collar. The dog is wearing a red collar with a black tag. The background is a bright, open field.

# **Gents** Camino Modular Collection

**All styles in this section are manufactured as illustrated.  
If there are to be any changes to the upper design or  
quantity of straps required then a special pattern change  
will need to be requested at the point of ordering.**

# Gents Camino Modular Collection

## Dolomite (Product Code Lace CA0021L, Hook & Loop CA0021V)

Dolomite is a 4 eyelet Derby plain front shoe. Finished with laid on heel counter and brogue punch pattern vamps.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. The standard opening means this style is not recommended for oedema or when easy access is required.
  3. No seams at the joints make it suitable for feet with painful or sensitive 1st or 5th joints.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in Dark Brown and Black

## Onyx (Product Code Lace CA0022L, Hook & Loop CA0022V)

Onyx is a 3 eyelet Derby plain front shoe.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. The standard opening means this style is not recommended for oedema or when easy access is required.
  3. No seams at the joints make it suitable for feet with painful or sensitive 1st or 5th joints.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in Black and Dark Burgundy

## Huntingdon (Product Code Lace CA0052L, Hook & Loop CA0052V)

Huntingdon is a 3 eyelet Derby type shoe with a neoprene vamp.

- Tips:**
1. The soft stretch front is ideal for sensitive feet or those with toe deformities.
  2. This low opening shoe is ideal for feet with variable oedema and when easy access is required.
  3. No seams at the joints make it suitable for feet with painful or sensitive 1st or 5th joints.
  4. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Neoprene vamp and calf rearfoot

**Sole Units:** All



Shown in Dark Brown and Black

## Durham

(Product Code Lace CA0053L, Hook & Loop CA0053V)

Durham is a 4 eyelet Derby plain front shoe.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. This low opening shoe is ideal for feet with variable oedema and when easy access is required.
  3. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in Dark Brown and Black

## Wiltshire

(Product Code Lace CA0030L, Hook & Loop CA0030V)

Wiltshire is a 4 eyelet Blucher style, seamless vamp with stitch apron design.

- Tips:**
1. The standard opening means this style is not recommended for oedema or when easy access is required.
  2. The stitched apron design means this style is not suitable for toe deformities or sensitive feet.
  3. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in Light Burgundy and Dark Brown

## Herefordshire

(Product Code Lace CA0055L, Hook & Loop CA0055V)

Herefordshire is a 4 eyelet Derby type shoe with a stretch leather apron.

- Tips:**
1. The soft stretch front is ideal for sensitive feet or those with toe deformities, especially bunions.
  2. This low opening shoe is ideal for feet with variable oedema and when easy access is required.
  3. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Stretch leather apron and calf rearfoot and toe bumper

**Sole Units:** All



Shown in Black and Dark Brown

# Gents Camino Modular Collection

## Buckingham (Product Code Lace CA0054L, Hook & Loop CA0054V)

Buckingham is a 3 eyelet Derby plain front shoe.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. The standard opening means this style is not recommended for oedema or when easy access is required.
  3. No seams at the joints make it suitable for feet with painful or sensitive 1st or 5th joints.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in Navy and Dark Brown

## Merseyside (Product Code Lace CA0056L, Hook & Loop CA0056V)

Merseyside is a 4 eyelet Blucher type plain front shoe.

- Tips:**
1. This mid opening shoe is suitable for feet with mild oedema.
  2. A stitched apron is not suitable for toe deformities
  3. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in Grey and Dark Brown

## Dorset (Product Code Lace CA0029L, Hook & Loop CA0029V)

Dorset is a 4 eyelet Derby style shoe with toe cap detail.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. The standard opening means this style is not recommended for oedema or when easy access is required.
  3. No seams at the joints make it suitable for feet with painful or sensitive 1st or 5th joints.
  4. Punch patterns cannot be added on the vamp.
  5. A stitched toe cap means this style is not suitable for toe deformities.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in Black and Mid Brown

## Rhyolite (Product Code Lace CA0023L, Hook & Loop CA0023V)

Rhyolite is a 5 eyelet Oxford type trainer with contrasting colour stripes.

**Tips:** 1. *The U opening front makes access and adjustability more difficult so this style is not recommended for oedema or when easy access is required.*

2. *Punch patterns cannot be added on the vamp.*



Shown in Navy with White stripe and White with Black stripe

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** Trainer

## Cambridge (Product Code Lace CA0057L, Hook & Loop CA0057V)

Cambridge is a 6 eyelet trainer boot with contrasting colour stripes.

**Tips:** 1. *The U opening front makes access and adjustability more difficult so this style is not recommended for oedema or when easy access is required.*

2. *Punch patterns cannot be added on the vamp.*



Shown in Black with White stripe and White with Navy stripe

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** Trainer

## Bedford (Product Code Lace CA0058L, Hook & Loop CA0058V)

Bedford is a 6 eyelet Derby trainer with contrasting colour stripes.

**Tips:** 1. *This low opening trainer is ideal for feet with variable oedema and when easy access is required.*

2. *Punch patterns cannot be added on the vamp.*



Shown in Black with Black Cordura and White leather stripes and Navy Blue with Navy Cordura and white leather stripes

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or Neoprene. All styles come with cordura panel at the quarters and the toe piece and are only available in Black, Navy, Silver and Taupe.

**Sole Units:** Trainer

# Gents Camino Modular Collection

## Pyrite (Product Code Lace CA0071L, Hook & Loop CA0071V)

Pyrite is a 5 eyelet trainer with contrasting colour inserts and separate toe bumper.

- Tips:** 1. *This mid opening trainer is suitable for feet with mild oedema.*  
2. *Punch patterns cannot be added on the vamp.*

**Fastenings:** Available in reverse pull hook and loop, and lace  
**Materials:** Not available with stretch leather or neoprene  
**Sole Units:** Trainer



Shown in White with Navy insert and Black with White insert

## Steel (Product Code Lace CA0072L, Hook & Loop CA0072V)

Steel is a 5 eyelet trainer contrasting colour inserts.

- Tips:** 1. *This mid opening trainer is suitable for feet with mild oedema.*

**Fastenings:** Available in reverse pull hook and loop, and lace  
**Materials:** Not available with stretch leather or neoprene  
**Sole Units:** Trainer



Shown in Blue with White Insert and Black with Light Burgundy insert

## Graphite (Product Code Lace CA0073L, Hook & Loop CA0073V)

Graphite is a Derby boot style with a combination of eyelets and metal hooks and contrasting colour heel/ankle inserts.

- Tips:** 1. *Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.*  
2. *The standard opening means this style is not recommended for oedema or when easy access is required.*

**Fastenings:** Available in reverse pull hook and loop, and lace  
**Materials:** Not available with stretch leather or neoprene  
**Sole Units:** All



Shown in Val Nut with a Black Cuff and Val Saddle With a Black Cuff

## Lava

(Product Code Lace CA0074L, Hook & Loop CA0074V)

Lava is a walking boot style with a combination of eyelets and lace loops and contrasting colour inserts.

- Tips:** 1. The standard opening means this style is not recommended for oedema or when easy access is required.
2. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in Black and Val Nut with a Val Saddale Cuff

## Sandstone

(Product Code, Lace CA0075L, Hook & Loop CA0075V)

Sandstone is a walking boot style with metal loop eyelets and contrasting colour inserts.

- Tips:** 1. The U opening front makes access and adjustability more difficult so this style is not recommended for oedema or when easy access is required.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in Val Nut with a Val Soan Cuff and Val Soan With a Val Saddale Cuff

## Granite

(Product Code Lace CA0076L, Hook & Loop CA0076V)

Granite is a Derby type boot.

- Tips:** 1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
2. The mid opening makes access and adjustability more difficult so this style is not recommended for oedema or when easy access is required.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in Mid Brown and Black

# Gents Camino Modular Collection

## Gloucester (Product Code Lace CA0032L, Hook & Loop CA0032V)

Gloucester is a 6 eyelet Derby boot.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. This low opening boot is ideal for feet with variable oedema and when easy access is required.
  3. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in  
Black and  
Saddale

## Fife (Product Code Lace CA0063L, Hook & Loop CA0063V)

Fife is a 6 eyelet Derby boot with contrasting colour ankle/heel insert.

- Tips:**
1. This mid opening shoe is suitable for feet with mild oedema.
  2. Punch patterns cannot be added on the vamp.
  3. Great walking boot when accompanied with a commando sole unit.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in  
Black and Nut/  
Cream Nubuck  
combination  
Shown with  
optional  
commando  
sole unit

## Hampshire (Product Code Lace CA0031L, Hook & Loop CA0031V)

Hampshire is a 7 eyelet Blucher type boot.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. This mid opening boot is suitable for feet with mild oedema.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in  
Black and  
Mid Brown

## Humber

(Product Code Lace CA0059L, Hook & Loop CA0059V)

Humber is a 6 eyelet Derby trainer boot with contrasting colour stripes.

- Tips:**
1. This mid opening boot is suitable for feet with mild oedema.
  2. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or Neoprene. All styles come with cordura panel at the quarters and the toe piece and are only available in Black, Navy, Silver and Taupe.

**Sole Units:** Trainer



Shown in White with Silver Cordura and Black leather stripes and Val Soan and cream with Taupe Cordura

## Tyneside

(Product Code Lace CA0060L, Hook & Loop CA0060V)

Tyneside is a 5 eyelet Derby boot.

- Tips:**
1. The soft stretch front is ideal for sensitive feet or those with toe deformities.
  2. This low opening boot is ideal for feet with variable oedema and when easy access is required.
  3. No seams at the joints make it suitable for feet with painful or sensitive 1st or 5th joints.
  4. Punch patterns cannot be added on the vamp.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Neoprene vamp and calf rearfoot

**Sole Units:** All



Shown in Black and Dark Brown

## Leicester

(Product Code Lace CA0061L, Hook & Loop CA0061V)

Leicester is a 6 eyelet Derby boot with brogue punch patterns.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. This mid opening boot is suitable for feet with mild oedema.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in Black and Dark Brown.

# Gents Camino Modular Collection

## Quartz (Product Code Lace CA0024L, Hook & Loop CA0024V)

Quartz is a 5 eyelet Derby plain front boot.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. The opening makes access and adjustability more difficult so this style is not recommended for oedema or when easy access is required.
  3. No seams at the joints make it suitable for feet with painful or sensitive 1st or 5th joints.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in  
Black and  
Dark Brown

## Tourmaline (Product Code Lace CA0077L, Hook & Loop CA0077V)

Tourmaline is a 4 metal loop eyelet Derby boot.

- Tips:**
1. Minimal seams on vamp makes the style ideal for diabetic patients or those with sensitive feet.
  2. The standard opening makes access and adjustability more difficult so this style is not recommended for oedema or when easy access is required.

**Fastenings:** Available in reverse pull hook and loop, and lace

**Materials:** Not available with stretch leather or neoprene

**Sole Units:** All



Shown in  
Grey and  
Val Nut



# Boot Last Measurement Guide

## Ladies

Style/Size		2	3	4	5	6	7	8	9
Opal Bootee	Height	107	110	113	116	119	122	125	128
	Top circumference	230	237	243	250	256	263	269	276
Opal	Height	142	146	149	153	156	160	163	167
	Top circumference	234	241	248	255	262	269	276	282
Zircon	Height	112	115	118	121	124	127	130	133
	Top circumference	226	233	240	247	254	261	268	275
Northumberland	Height	99	102	105	108	111	114	117	120
	Top circumference	234	241	247	254	260	267	274	280
Berkshire	Height	109	112	115	118	121	124	127	130
	Top circumference	222	229	235	242	248	255	261	268
Sussex	Height	112	115	118	121	124	127	130	133
	Top circumference	226	233	240	247	254	261	268	275
Gold	Height	99	102	104	106	109	112	115	118
	Top circumference	220	230	236	244	250	258	266	274
Bronze	Height	105	108	111	114	116	119	121	124
	Top circumference	232	236	244	250	258	266	272	278
Topaz	Height	103	106	109	112	114	117	120	123
	Top circumference	224	232	240	248	256	264	272	278

## Gents

Style/Size		5	6	7	8	9	10	11	12
Quartz	Height	122	125	128	131	134	137	140	143
	Top circumference	234	241	274	254	260	267	273	280
Hampshire	Height	117	120	123	126	129	132	135	138
	Top circumference	250	257	264	271	278	285	292	299
Gloucester	Height	107	110	112	115	117	120	122	125
	Top circumference	230	236	242	248	254	260	266	272
Leicester	Height	123	126	129	132	135	138	141	144
	Top circumference	228	234	240	246	252	258	264	270
Fife	Height	109	112	115	118	121	124	127	130
	Top circumference	230	237	243	250	256	263	269	276
Cambridge	Height	118	121	124	127	130	133	136	139
	Top circumference	240	246	252	258	264	270	276	282
Humber	Height	118	121	124	127	130	133	136	139
	Top circumference	240	246	252	258	264	270	276	282
Tyneside	Height	122	125	128	131	134	137	140	143
	Top circumference	234	241	247	254	260	267	273	280
Tourmaline	Height	119	123	125	128	131	134	137	140
	Top circumference	236	244	250	254	260	268	274	280
Graphite	Height	120	123	127	130	133	136	139	142
	Top circumference	240	246	254	260	266	272	278	284
Granite	Height	121	124	127	129	132	135	138	141
	Top circumference	236	244	252	256	262	268	274	282
Sandstone	Height	119	122	125	128	130	133	136	139
	Top circumference	236	242	248	256	262	270	272	282
Lava	Height	119	120	123	127	129	132	135	138
	Top circumference	236	244	250	256	262	268	276	282

Measurements for all sections on these pages are shown in mm.

# Last Measurement Guide

Lasts are available in both square-toe and round-toe options. Please note that if unspecified, round-toe will be supplied by default.

## Ladies

Width/Size		2	3	4	5	6	7	8	9
Extra Narrow	Joint width	79	81	83	85	87	89	91	93
	Joint circumference	201-219	207-225	214-232	220-238	226-244	233-251	239-257	245-263
	Instep circumference	197-215	203-221	210-228	216-234	222-240	229-247	235-253	241-259
	Toe depth	31	32	33	34	35	36	37	38
	Heel Width	47	49	49.5	50	52	53	54	56
Narrow	Joint width	83	85	87	89	91	93	95	97
	Joint circumference	210-228	216-234	223-241	229-247	235-253	242-260	248-266	254-272
	Instep circumference	206-224	212-230	219-237	225-243	231-249	238-256	244-262	250-268
	Toe depth	31	32	33	34	35	36	37	38
	Heel Width	48	49	50	52	53	54	55	57
Standard	Joint width	89	91	93	95	97	99	101	103
	Joint circumference	226-244	232-250	239-257	245-263	251-269	258-276	264-282	270-288
	Instep circumference	211-229	217-235	224-242	230-242	236-254	243-261	249-267	255-273
	Toe depth	31	32	33	34	35	36	37	38
	Heel Width	49	50	51	53	54	55	56	57
Wide	Joint width	94	96	98	100	102	104	106	108
	Joint circumference	236-254	242-260	249-267	255-273	261-279	268-286	274-292	280-298
	Instep circumference	217-235	223-241	230-248	236-254	242-260	249-267	255-273	261-279
	Toe depth	31	32	33	34	35	36	37	38
	Heel Width	52	53	54	55	56	57	59	60
Deep	Joint width	89	91	93	95	97	99	101	103
	Joint circumference	236-254	242-260	249-267	255-273	261-279	268-286	274-292	280-298
	Instep circumference	217-235	223-241	230-248	236-254	242-260	249-267	255-273	261-279
	Toe depth	37	38	39	40	41	42	43	44
	Heel Width	52	53	54	55	56	57	59	60

### Finished Back and Malleoli Standard Measurements

Width/Size		2	3	4	5	6	7	8	9
Extra Narrow, Narrow, Standard and Wide	Back Height	72	72	74	76	78	80	82	84
	Side Height	57	57	59	61	63	65	67	69
Deep	Back Height	78	78	80	82	84	86	88	90
	Side Height	60	60	62	64	66	68	70	72

Tolerance on back height (+/-5mm).

## Gents

Width/Size		5	6	7	8	9	10	11	12
Narrow	Joint width	89	92	93	95	97	99	100	102
	Joint circumference	225-243	231-249	238-256	244-262	250-268	257-275	263-281	269-287
	Instep circumference	215-233	221-239	228-246	234-252	240-259	247-265	253-271	259-277
	Toe depth	37	38	39	40	41	42	43	44
	Heel Width	55	57	58	60	61	63	64	65
Standard	Joint width	95	97	99	101	103	105	107	109
	Joint circumference	237-255	243-261	250-268	256-274	262-280	269-287	275-293	281-299
	Instep circumference	227-245	233-251	240-258	246-264	252-270	259-277	265-283	271-289
	Toe depth	37	38	39	40	41	42	43	44
	Heel Width	59	60	61	63	64	65	66	67
Wide	Joint width	101	103	105	107	109	111	113	115
	Joint circumference	249-267	255-273	262-280	268-286	274-292	281-299	287-305	293-311
	Instep circumference	235-253	241-259	248-266	254-272	260-278	267-285	273-291	279-297
	Toe depth	37	38	39	40	41	42	43	44
	Heel Width	62	64	65	66	68	69	70	71
Deep	Joint width	95	97	99	101	103	105	107	109
	Joint circumference	249-267	255-273	262-280	268-286	274-292	281-299	287-305	293-311
	Instep circumference	235-253	241-259	248-266	254-272	260-278	267-285	273-291	279-297
	Toe depth	43	44	45	46	47	48	49	50
	Heel Width	62	64	65	66	68	69	70	71

### Finished Back and Malleoli Standard Measurements

Width/Size		5	6	7	8	9	10	11	12
Narrow, Standard and Wide	Back Height	76	78	80	82	84	86	88	90
	Side Height	61	63	65	67	69	71	73	75
Deep	Back Height	82	84	86	88	90	92	94	96
	Side Height	64	66	68	70	72	74	76	78

Tolerance on back height (+/-5mm).

# Fitting Solutions

Area of Shoe	Difficulty Identified	Visual Signs	Possible Causes	Solutions
Heel Counter	Too High	Rubbing under either medial or lateral malleoli.	Height of malleoli not measured and stated on the draft. Ankle oedema. Abnormal malleoli shapes and heights.	Draw a line on the shoe counter with a leather pen, to identify where it needs lowering. Provide a measurement in mm, and request that the counter is lowered. This can be altered at fitting stage. On a finished item, add 3 mm cork heel lift under inlay to lift heel higher in the counter to relieve any pressure. Ensure malleolar heights are recorded on future orders.
		Catching the Achilles tendon.	Presence of 'heel bumps' or bursae over Achilles tendon	Soften the counter with a leather softener and use a ring and ball stretcher to create a 'bump' in counter. Encourage gradual 'wearing in' over a few days. Made to Measure items can be ordered with a cavity in the counter stiffener to accommodate the bump or a straighter heel seam. Supply a cast.
	Too Low	Does not grip the foot properly and foot comes out of the shoe too easily. Feels unsafe in wearing.	Counter depth sits too low on the foot.	On unfinished footwear, request counter depth to be increased and give measurements in mm.
			Inlays may be too thick under heel, and lifting foot out.	Reduce inlay thickness.
		Inside heel raise is too high for the counter.	Reduce inside heel raise height, and add balance to external heel height, or change to external shoe raise.	
		Excessive oedema affects how the heel sits down in the counter.	Straighten the back curve to accommodate oedema in Made to Measure footwear or and increase the width of heel seat in modular.	
Heel to Ball Length	Too Long	Footwear causes pain/ discomfort; foot gets 'tired'. Creases too far back. Check for wear marks on sole behind flexion point.	Foot is flexing too far back in shoe, incorrect fitting.	Use a last which has a shorter heel-to-ball measurement. If Made to Measure, re-last with a shorter heel-to-ball measurement after fitting.
	Too Short	Heavily creased vamp below facings. Toe puff may be creased as well. Toes often sore as shoes too short. 'Turned up' toe end. Heavy wear marks placed forward of flexion point.	Foot is flexing too far forward in shoe, incorrect fitting.	Use a last with a longer heel to ball measurement. If Made to Measure, re-last with a longer heel to ball fitting. (This problem is rarely found with Made to Measure footwear if drafts and measurements have been taken correctly).
1st MTP Joint	Limited Range of Movement (ROM)	Twisted or angled creases that seem unusual. Wear marks heavy to lateral side.	Check first metatarsal-phalangeal joint: ROM may be fixed or severely limited.	Add a rocker sole to compensate for the reduced movement.
Instep	Too Narrow	Facings pulled wide apart and possibly at wrong angle. Very tight; may not get shoe onto foot.	Oedema changes or incorrect instep circumference and width measurements supplied.	Increase instep circumference after fitting. Reduce insole thickness to increase volume at instep.
		Quarters bulge on medial side.	Excess pronation may cause poor foot position in shoe.	Use a corrective orthosis. If correction is not possible, make Made to Measure footwear and supply a cast so footwear can accommodate pronation.
		Incorrect style: does not suit foot shape.	Wrong style?	Use a low opening style to improve access.
	Too Wide	Facings feel loose; excess leather, particularly along medial side of quarter.	Fitting too wide. Entry probably too wide as well.	Use a low opening style to improve access
		Facings butt or overlap. Foot not retained adequately into counter.	Incorrect long heel measurement, or reduction in measurement due to oedema or change of dressings, etc. at instep	Change to narrower fitting, check foot on inlay for heel seat and instep fitting. Measure heel width: if it is narrow, it could also affect the instep.
			Inlays too thin.	Check measurements. Reduce instep and heel circumference after fitting. Add tongue padding on completed shoes.
		Inlays too thin.	Increase inlays but beware as this may compromise the fit elsewhere.	

Area of Shoe	Difficulty Identified	Visual Signs	Possible Causes	Solutions
Vamp and Toe Box	Too Wide	Feels big.	Change in width or circumference measurements, or incorrect measurements supplied.	Check draft: and measurements. Remeasure if necessary. Can be re-lasted after fitting.
				Soften the counter with a leather softener and use ring and ball stretcher to create a 'bump' in counter. Encourage gradual 'wearing in' over a few days. Made to Measure items can be ordered with a cavity in the counter stiffener to accommodate the bump.
		Deep fold creases on flexion.	Excess leather.	Feel position of the forefoot in vamp. Draw with a leather pen where the excess is situated. Request reduction of insole board in mm. Check depth as this may also need to be reduced.
	Too Deep	Feels 'baggy'. Excess leather across dorsum. Deep creases on flexion.	Change in measurement; reduction of oedema.	On unfinished footwear, request counter depth to be increased and give measurements. Stock items: Consider a different design. Reduce inside heel raise height, and add balance to external heel height, or change to external shoe raise only.
		Base of facings bulge up from vamp.	Instep/entry is also too deep.	Request modular or Made to Measure only. Straighten the back curve to accommodate oedema and increase width of heel seat.
	Too Narrow	Feels tight.	Toe box shape of shoe not compatible with foot shape.	Compare foot shape with shoe; request change in toe box shape, if necessary, to rounded or square. Use inlay to check for width.
		Forefoot squeezed.	Swelling may have occurred.	Refit once swelling has reduced, or monitor swelling until stabilized.
		Stretching and creases visible.	Increase in toe deformities.	Remeasure and order a different size Stretch footwear if possible and thin the insole to give more room.
		Toe puff may press into toes. Toes bulge into upper on lateral border.	Too narrow in toe box, or toe box wrong shape for foot.	Request change of toe shape and re-last the upper. (Note: increasing the width will also increase the circumference.)
	Too Shallow	May not get shoe on.	Incorrect measurement.	Check draft, remeasure and request increase in depth giving measurement in mm.
		Very tight.	Change in measurements.	Request increase in depth, giving measurement in mm.
		Toes rub on vamp leather.	No toe box depth measurement given on draft.	Request increase in depth, giving measurement in mm.
		One high toe may be pushing up on the upper.	Position of high toe not given on draft.	Mark tight area with leather pen and add addition to last if footwear is at fitting stage. If not at a fitting stage provide new measurements for a new pair.

Distributed by

**Blatchford**

