



PRO-FIT LINER

PFGFR | PFGDT

Ensuring the perfect fit.

DESCRIPTION

ALPS® Pro-Fit is a prosthetic liner made with OptiGel that ensures excellent conformability to the residual limb. At the same time, its new and improved distal end encapsulation allows the user to experience a highly secure fit, resulting in enhanced control, comfort, and safety.

PART #	SUSPENSION
PFGDT	LOCKING
PFGFR	CUSHION

Recommended for K2-K4 activity levels.



TOP FEATURES & BENEFITS


PREMIUM GREY FABRIC

OPTIGEL

CUSHION & LOCKING


- ✓ New encapsulation allows the liner to fit better on the residual limb.
- ✓ The superior conformability of the OptiGel allows for safe suspension of the prosthetic leg.
- ✓ Improves the control of the prosthetic leg and reduces pistoning, while still providing superior protection of the residual limb.

CROSS SECTION DIAGRAM



HOW TO SIZE

Measurements taken 6cm from the distal end



SIZING CHART

MEASURED SIZE (CM)	LOCKING 3MM OR 6MM UNIFORM	CUSHION 3MM OR 6MM UNIFORM
20-23 cm	PFGDT20-(3 or 6)	PFGFR20-(3 or 6)
24-25 cm	PFGDT24-(3 or 6)	PFGFR24-(3 or 6)
26-27 cm	PFGDT26-(3 or 6)	PFGFR26-(3 or 6)
28-31 cm	PFGDT28-(3 or 6)	PFGFR28-(3 or 6)
32-37 cm	PFGDT32-(3 or 6)	PFGFR32-(3 or 6)
38-43 cm	PFGDT38-(3 or 6)	PFGFR38-(3 or 6)
44-53 cm	PFGDT44-(3 or 6)	PFGFR44-(3 or 6)

Blatchford distribute these products on behalf of ALPS® who have supplied all technical content for this publication. Please refer to the product IFU for additional information.